COLESTEROL HDL ENZIMÁTICO

REF **K015**

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of HDL Cholesterol. Colorimetric enzymatic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic - Trinder

The VLDL Lipoproteins (Very Low Density Lipoprotein), and LDL (Low Density Lipoprotein) and the Kilomicrons are precipitate with the mixture of Phosphotungsten Acid and Magnesium Chloride. After centrifugation, the Cholesterol binded to the High Density Lipoprotein (HDL) is determined by enzymatic colorimetric method.

REAGENTS

Number 1 - Standard - Store between 2 and 8°C. Contains: Cholesterol 40 mg/dL and n-Butanol P.A.

Number 3 - Precipitate Reagent - Store between 2 and 8°C. Contains: Phosphotungstic Acid 1,5mmol/L and Magnesium Chloride 100 mmol/L.

PRESENTATION

Reagent	Volume
Nº 1	2 mL
Nº 2	25 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water bath at 37°C, watch or stopwatch, pipettes, test tubes, Moreagent Cholesterol kit, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Water level in water bath should be above the level of reagents in test tubes.

6- The use of Enzymatic Reagent of **Bioclin** Cholesterol Monoreagent kit is of fundamental importance to a good performance of the system.

7- Determine the factor periodically in every batch of the product.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolyses, collected after fasting period of 12 to 14 hours. HDL Cholesterol is stable in serum for up to 7 days between 15°C and 30°C and 14 days between 2 and 8°C.

PROCESS DESCRIPTION

TECHNIQUE

Bioclin recommends, use the Monoreagent Cholesterol (Enzymatic Reagent*) kit and as control serum, Biocontrol N and P Bioclin Kits.

Precipitation of VLDL, LDL and Kilomicrons

In a mixer tube pipette:

Serum 250 µL

Precipitant Reagent Nº 2 250 µL

Agitate in Vortex or manually during 1 minute. Centrifuge in 3500 rpm during 15 minutes. This step is of fundamental importance to obtain the correct result. The limpid supernatant contains HDL Cholesterol and shall be pipette immediately after centrifugation, to avoid wrongly elevated results.

Determination of HDL Cholesterol

Mark 3 test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Supernatant	---	---	50 µL
Reagent Nº 1	---	50 µL	---
Enzymatic Reagent*	1,0 mL	1,0 mL	1,0 mL

Homogenize well and put it in water bath 37 °C for 5 minutes. Read the absorbance of Sample and Standard at 500 nm (490 - 540 nm), matching the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

Since it is require for the technique the diluted serum at 1:2 with the Precipitate Reagent, it is necessary to correct this dilution to obtain the correct results. For that reason, during the calculations, the standard value is corrected to the double of its real concentration, which means, 80 mg/dL.

$$\text{HDL (mg/dL)} = \frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times 80$$

Reaction follows the Lambert-Beer Law, the calibration factor can be used.

$$\text{Calibration Factor} = \frac{(80 \text{ mg/dL})}{\text{Standard Absorbance}}$$

$$\text{mg/dL} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

The VLDL and LDL Cholesterol concentration can be calculated according to the Friedewald equation. This calculation could be used for samples which Triglycerides concentration does not pass 400 mg/dL.

$$\text{VLDL Cholesterol} = \text{Triglycerides} / 5$$

$$\text{LDL Cholesterol} = \text{Cholesterol total} - (\text{HDL} + \text{VLDL})$$

For lipemic sample or cloudy supernatant, dilute the sample 1:2 with Sodium Chloride 0,85% and repeat the precipitation. Multiply the result obtained by 2. If the supernatant remains cloudy, the sample can not be used.

PROCEDURE LIMITATIONS

Keep the ratio of 1:1 between Sample and Precipitating. Some substances, such as Ascorbic Acid, Hemoglobin above 150 mg/dL, Bilirubin above 20 mg/dL and Triglycerides above 600 mg/dL, interfere with the reaction.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 1951.

REFERENCE VALUES

The reference values in mg/dL, for this method, were obtained through the determination of HDL Cholesterol in healthy populations of male and female.

Adults (≥ 20 years)	Desirable	> 60 mg/dL
	Limit	40 - 60 mg/dL
	Low	< 40 mg/dL

Childrens and Adolescents (2 to 19 years)	Desirable	> 45 mg/dL
--	-----------	------------

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE**QUALITY CONTROL****Accuracy****COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY**

The HDL Cholesterol Enzymatic kit was compared with another commercially available kit of measurement of HDL Cholesterol. 7 analyses were conducted and the results were evaluated, obtaining the linear equation of $Y = 1,029X - 2,114$ with a linear correlation coefficient of 0,999. With these results we can conclude that the kit shows good methodological specificity.

PRECISION**REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	35,15	50,70	70,75
Standard Deviation (mg/dL)	0,67	0,73	0,72
Coefficient of Variation (%)	1,91	1,45	1,01

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	34,90	50,72	70,77
Standard Deviation (mg/dL)	0,22	0,03	0,03
Coefficient of Variation (%)	0,62	0,06	0,04

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of HDL Cholesterol. The average found was 1,3307 mg/dL, with standard deviation of 0,0005 mg/dL. The sensitivity, that indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 1,3321 mg/dL.

DIAGNOSTIC SIGNIFICANCE

The determination of High Density Lipoproteins (HDL) and Low Density Lipoprotein (LDL) play an important role in prevention and therapy of coronary heart disease.

The HDL are capable of captating the Cholesterol of the ppheriferal tissue and reconducting to the liver in order to degradate and excrement it. This mechanism denominated reverse transport of Cholesterol, constitute the theory of HDL anti-atherogenicity, where low concentration of HDL are directly related with the incidence of cardiovascular diseases.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP **OBELIS S.A.**

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1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for HDL Cholesterol Enzymatic kit: 10269360082

Review: August/2014

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED