# Bioclin

## BILIRRUBINA

REF K005

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#### FINALIDADE

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for determination of Direct and Total Bilirubin. Colorimetric test, for *in vitro* diagnostic only.

#### PRINCIPLE OF ACTION

**Methodology:** Sims - Horn

The Bilirubin through the reaction of attachment with the Sulphanilic Diazotized Acid forms a complex of red color (Azobilirubin). The Total Bilirubin (Direct and Indirect) is dosed in the presence of an accelerator (Caffeine and Benzoate) that allows the solubilization of the Indirect Bilirubin (not conjugated) insoluble in water. The Direct Bilirubin (conjugated) is dosed in a wet environment.

#### REAGENTS

**Number 1 - Sodium Nitrate** - Store between 15 and 30°C. Contains: Sodium Nitrate 70 mmol/L.

**Number 2 - Sulphanilic Reagent** - Store between 15 and 30°C. Contains: Sulphanilic Acid 6 mmol/L and Chloridric Acid 130 mmol/L.

**Number 3 - Accelerator** - Store between 15 and 30°C. Contains: Caffeine 0,2 mmol/L, Sodium Benzoate 0,4 mol/L and Sodium Acetate 0,7 mol/L.

#### PRESENTATION

Reagent	Volume
Nº 1	5 mL
Nº 2	125 mL
Nº 3	250 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometers or colorimeter, watch and stopwatch, pipettes, test tubes, Bilirubin Calibration Bioclin kit, they can be found at markets specialized in laboratories of clinical analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 15 and 30°C. Protect from light and avoid moisture.

#### SPECIAL CARE

**1- For professional *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain exact results.

**3-** Water used in material cleaning must be recent and free of contaminants.

**4-** Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

**5-** The calibration shall be repeated periodically to verify any change in the colorimeter response or spectrophotometer.

**6-** Hemolysis, even discrete, interferes in the dosage.

**7-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**8-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**9-** Do not use the product in case of damaged packaging.

**10-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum obtained free of hemolysis. The analyte is stable for 3 days when the sample is stored between 2 and 8°C protected from light. The sample for therapeutical control must be collected always in the same time.

#### PROCESS DESCRIPTION

##### DIAZO REAGENT PREPARATION

Add 1 drop of Reagent Nº 1 at 1,5 mL of Reagent Nº 2. Homogenized well and use only on the day of preparation.

#### TECHNIQUE

##### MACRO TECHNIQUE

Label 3 test tubes: "B" (Blank), "BD" (Direct Bilirubin) and "BT" (Total Bilirubin) and proceed as follows:

	Blank	BD	BT
<b>Reagent Nº2</b>	500 $\mu$ L	---	---
<b>Diazo Reagent</b>	---	500 $\mu$ L	500 $\mu$ L
<b>Dist. Water</b>	4,5 mL	4,5 mL	---
<b>Reagent Nº3</b>	---	---	4,5 mL
<b>Sample</b>	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L

Homogenize by mixing, wait 5 minutes and read the Total and Direct Bilirubin absorbancies at 525 nm (500 - 550 nm), adjusting the zero with the Blank. The color is stable for 30 minutes.

The values in mg/dL for Direct and Total Bilirubin are obtained by the Calibration Factor MACRO, acquired with the Bilirubin Calibration Bioclin kit.

##### MICRO TECHNIQUE

Used in bilirubin dosages with concentration superior to 10 mg/dL (jaundiced recently born serum). Use a 0,05 mL micropipette to pipette the serum.

Label 3 test tubes and proceed as follows:

	Blank	BD	BT
<b>Reagent Nº 2</b>	150 $\mu$ L	---	---
<b>Diazo Reagent</b>	---	150 $\mu$ L	150 $\mu$ L
<b>Distilled Water</b>	1,8 mL	1,8 mL	---
<b>Reagent Nº 3</b>	---	---	1,8 mL
<b>Sample</b>	50 $\mu$ L	50 $\mu$ L	50 $\mu$ L

Homogenize by mixing, wait 5 minutes and read the Direct and Total Bilirubin absorbancies at 525 nm (500 - 550 nm), adjusting the zero with the Blank. The color is stable for 30 minutes.

The values in mg/dL for Direct and Total Bilirubin are obtained by the Calibration Factor MACRO, acquired with the Bilirubin Calibration Bioclin kit.

#### CALCULATIONS

Using the Bilirubin Calibration Bioclin\* kit, effectuate the calculation factor. Carry out distinct calculations for Macro technique and Micro technique.

Calibration Factor =  $\frac{\text{Standard Concentration}^* (10 \text{ mg/dL})}{\text{Standard Absorbance}^*}$

mg/dL = Sample Absorbance x Calibration Factor

Indirect Bilirubin = Total Bilirubin – Direct Bilirubin

#### Example:

Satndard Average Absorbance = 0,480

Total Bilirubin Sample Absorbance = 0,045

Direct Bilirubin Sample Absorbance = 0,012

Calibration Factor =  $\frac{10 \text{ mg/dL}}{0,480} = 20,83$

Total Bilirubin =  $0,045 \times 20,83 = 0,93$

Direct Bilirubin =  $0,012 \times 20,83 = 0,25$

Indirect Bilirubin =  $0,93 - 0,25 = 0,68$

#### PROCEDURE LIMITATIONS

Do not use the Macro technique for Bilirubin values below 10 mg/dL. Some compounds as Anabolic Steroids, Ascorbic Acid, Salicylates and Vitamin A can elevate the Bilirubin values. The Caffeine and the Penicillin can interfere in the dosages reducing the Bilirubin.

#### INTERFERENTES

Samples with Hemolysis, even discrete, interferes in the dosage of Bilirubin.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

Bioclin recommends the usage of the Bioclin Bilirubin Calibration kit, which is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 916A.

#### REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of Bilirubin in healthy populations of male and female.

Adults	
<b>Total</b>	0,1 - 1,2 mg/dL
<b>Direct</b>	up to 0,4 mg/dL
<b>Indirect</b>	0,1 - 0,8 mg/dL

Newborns and Childrens		
<b>Bilirubin Total</b>	<b>24 h</b>	< 8,8 mg/dL
	<b>2º day</b>	1,3 - 11,3 mg/dL
	<b>3º day</b>	0,7 - 12,7 mg/dL
	<b>4º to 6º day</b>	10,1 - 12,6 mg/dL
	<b>Childrens &gt; 1 month</b>	0,2 - 1,0 mg/dL
<b>Bilirubin Direct</b>	<b>Newborns and Childrens</b>	up to 0,4 mg/dL

To convert the values of mg/dL to mmol/L (SI) multiply by 0,0171.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

**PRODUCT PERFORMANCE****QUALITY CONTROL****Accuracy****RECOVERY**

The recovery analysis was made with 5 sample determinations. The accuracies were calculated, and if they are found in a good concordance with the reference values, obtaining a recovery between 95 and 103%.

**COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY**

The Bilirubin kit was compared with other commercially available methods for measurement of Total and Direct Bilirubin. 07 analyzes were performed and the results were evaluated.

For Total Bilirubin the linear equation obtained was  $Y = 0,974X + 0,025$  and correlation coefficient 0,998. For Direct Bilirubin the linear equation obtained was  $Y = 1,029X + 0,004$  and correlation coefficient 0,997.

**Precision****REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Direct Bilirubin			
	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0,22	0,32	0,84
Standard Deviation (mg/dL)	0,01	0,01	0,01
Coefficient of Variation (%)	3,61	2,47	0,94

Total Bilirubin			
	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0,43	1,04	3,23
Standard Deviation (mg/dL)	0,01	0,01	0,01
Coefficient of Variation (%)	2,23	1,02	0,37

**REPRODUCIBILITY**

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Direct Bilirubin			
	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0,22	0,32	0,84
Standard Deviation (mg/dL)	0,01	0,01	0,01
Coefficient of Variation (%)	2,59	1,79	0,69

Total Bilirubin			
	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0,43	1,03	3,23
Standard Deviation (mg/dL)	0,01	0,01	0,01
Coefficient of Variation (%)	1,33	0,56	0,18

**Sensitivity**

The sensitivity was calculated from twenty (20) determinations of a sample free of Total and Direct Bilirubin.

The average of the Total Bilirubin found was 0,03600 mg/dL with standard deviation of 0,02542 mg/dL. The sensitivity, which indicates the method detection limit, corresponds to the average plus 3 times the standard deviation, and is equal to 0,11227mg/dL.

The average of the Direct Bilirubin found was 0,0330 mg/dL with standard deviation of 0,0047 mg/dL. The sensitivity, which indicates the method detection limit, corresponds to the average plus 3 times the standard deviation, and is equal to 0,0471 mg/dL.

**Linearity**

The reaction is linear up to 15,0 mg/dL for the Macro Technique and 45 mg/dL for Micro Technique. For samples with higher values dilute the sample with NaCl 0,85% and repeat the dosage. Multiply the result by the dilution factor.

**DIAGNOSTIC SIGNIFICANCE**

The Bilirubin lift on plasma in the presence of hepatic injuries, bile obstruction or when the speed of destruction of red globule is increased. The increase of the Indirect Bilirubin is observed on the hemolytic syndrome, neonatal jaundice, Cligler-Najjar syndrome, Gilbert disease.

The Direct Bilirubin is increased on acute and chronicle hepatitis, the toxic reactions to many drugs (Chlorpromazine, Arsenicals Organics, Methyltestosterone), in bile obstruction treatment.

**NUMBER OF TESTS**

110 Tests / 300  $\mu$ L of Sample

276 Tests / 50  $\mu$ L of Sample

**BIBLIOGRAPHIC REFERENCES**

- 1 - JENDRASSIK, H. T. e Grof, P. : Biochem. Zentr. 297:81, 1938.
- 2 - MALLOY, H. T. e EVELYN, K. A. - J. Biol. Chem. 119:481, 1937.
- 3 - MARTINEK, R. G., Clin. Chem. Acta 13, 161, 1966.
- 4 - POWELL, W. N., Am. J. Clin. Path. 8, 55, 1944.
- 5 - SIMS, F. H. e Horn, C., Am. J. Clin. Path. 29, 412, 1958.
- 6 - TONKS, D. B., Clin. Chem. 9, 217, 1983.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

**EC REP OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium


**CUSTOMER SERVICE**

Customer Advisory Service  
Phone: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for Bilirubin kit: 10269360083

Review: May/2013

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED