URÉIA UV

REF K056

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of the Urea. Kinetic test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Fixed time Kinetic

Urea is hydrolyzed to NH_3 and CO_2 by Urease.

The Glutamate Dehydrogenase (GLDH) in the presence of NH_3 and α -Ketoglutarate, oxidizes NADH to NAD^+ .

The oxidation of NADH to NAD^+ , as measured by the decrease in absorbance is proportional to the concentration of urea in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Tris Buffer (pH 7,6) 100 mmol/L, ADP 0,7 mmol/L, α -Ketoglutarate 9 mmol/L, Sodium Azide 15,38 mmol/L, Urease \geq 6500 U/L, Glutamate Dehydrogenase \geq 1100 U/L.

Number 2 - Coenzyme - Store between 2 and 8°C. Contains: NADH 0,32 mmol/L.

Number 3 - Standard - Store between 2 and 8°C. Contains: Urea 70,0 mg/dL.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2	Reagent N° 3
1	2 x 80 mL	1 x 40 mL	1 x 3 mL
2	1 x 40 mL	1 x 10 mL	1 x 3 mL
3	2 x 40 mL	2 x 10 mL	1 x 3 mL
4	4 x 40 mL	4 x 10 mL	1 x 3 mL
5	2 x 40 mL	1 x 20 mL	1 x 3 mL
6	4 x 40 mL	2 x 20 mL	1 x 3 mL
7	4 x 60 mL	4 x 15 mL	1 x 3 mL
8	5 x 60 mL	5 x 15 mL	1 x 3 mL
9	6 x 60 mL	6 x 15 mL	1 x 3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatches, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Handle Reagent N° 1 with care, for it contains Sodium Azide which can be irritating to skin and mucous membranes.

6- Determinate the factor periodically and each batch of product.

7- Heavy metals and fluoride (above 2 mg/dL) are inhibitors of Urease.

8- No smoking near the place where they carry out the measurement, because the smoke contains ammonia fumes that contaminate the sample and glassware, leading to falsely increase.

9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

10- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

11- Do not use the product in case of damaged packaging.

12- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum, plasma (EDTA or heparin) and urine. Do not use anticoagulants containing ammonia.

Urea is stable in plasma or serum for 7 days between 2 and 8°C and for 90 days at 10°C negative.

Moderate hemolysis and Bilirubin up to 20 mg/dL do not produce significant changes in the results. For determination of Urea in the urine sample collected 24 hours in a vial containing 2,0 mL of HCl 50% (v/v). Centrifuge the sample before starting the technique and to conduct analysis within a few hours, because the urea excreted in urine is easily decomposed by bacterial action.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENT

Mix 4 parts of Reagent N° 1 (Buffer) with 1 part Reagent N° 2 (Coenzyme). The Working Reagent is stable for 05 days between 20 and 30°C and 30 days between 2 and 8°C.

REACTION CONDITIONS

Is indispensable conditions the use a thermostated cuvette at 37°C (25 to 30°C), 1 cm light path and read at 340 nm (334 - 365 nm).

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

Place in a test tube 1,0 mL of Working Reagent and add 10 μL of Sample or Standard, mix and transfer to a thermostated cuvette at 37°C. Trigger the stopwatch simultaneously and measure the absorbance at 30 and 90 seconds at 340 nm (334 - 365 nm). Differences in absorbance (ΔA) between two times, the Standard and the Sample, will be used to calculate the results.

For determination of Urea in urine, following the above technique use sample diluted 1:50 (100 μL urine + 4,9 mL distilled or deionized water). Multiply the result by 50.

CALCULATIONS

$$\text{Urea (mg/dL)} = \frac{\Delta A \text{ of Sample}}{\Delta A \text{ of Standard}} \times 70$$

As the reaction follows the Beer-Lambert Law, the calibration factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (70 mg/dL)}}{\Delta A \text{ of Standard}}$$

$$\text{Urea (mg/dL)} = \Delta A \text{ of sample} \times \text{Calibration Factor}$$

Urine

$$\text{Urea (g/24h)} = \frac{\Delta A \text{ of Sample} \times \text{Factor} \times 50 \times \text{Volume (L)}}{100}$$

The results for serum or plasma are expressed in mg/dL and urine g/24h.

PROCEDURE LIMITATIONS

Contamination of water, glass-making environment and/or the sample with ammonia, may lead to falsely elevated results. Some drugs such as Aminoglycosides, Cephalosporins, Purinol, Methyldopa, Furosemide, Propranolol can produce interference in the results, increasing the values of Urea.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

the kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 912.

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of urea in healthy populations of male and female.

Serum or plasma	
1 day - 12 month	2 - 34 mg/dL
1 - 13 years	8 - 36 mg/dL
Adults	15 - 40 mg/dL
Urine	
All ages	26 - 43 g/24h

To convert the values to mg/dL to mmol/L (SI), multiply by 0,166. These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Urea UV kit was compared with another method for dosage of Urea, commercially available. 07 analyzes were conducted and the results were evaluated. The linear equation obtained was $Y = 1,02X - 1,0842$ with correlation coefficient equal to 09918. With these results, we can conclude that the kit shows good specificity methodology.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	21,85	33,80	17,50
Standard deviation (mg/dL)	0,88	1,24	0,83
Coefficient of variation (%)	4,01	3,67	4,73

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	21,93	33,17	17,45
Standard deviation (mg/dL)	0,08	0,55	0,09
Coefficient of variation (%)	0,35	1,66	0,50

Sensibility

The sensibility was calculated from 20 determinations of a sample free of Urea. The average found was 2,140 mg/dL, with standard deviation of 0,278 mg/dL. The sensibility, which indicates the detection limit of the method, corresponds the average plus 3 times the standard deviation, and is equal to 2,974 mg/dL.

Linearity

The reaction is linear up to 300 mg/dL. For samples with values above 300 mg/dL or optical density above 0,8, dilute it with distilled or deionized water, repeat the dosage and multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Urea, the main product of the catabolism of proteins and amino acids, has its serum concentration affected by diet and hydration status, providing a rough indication of the state of renal function.

Elevations values of plasma Urea are classified as:

Pre-Cause kidney: excretion resulting from defects observed in cardiac decompensation, hemorrhagic shock, acute dehydration, high protein catabolism (burns, fever).

Cause renal: as a consequence of acute or chronic kidney disease (glomerulonephritis, pyelonephritis, tubular necrosis) with levels of Urea plasma of 300 mg/dL or more.

Post-kidney causes: usually resulting from a urinary tract obstruction, can occur in kidney lithiasis tumors by compression of the bladder.

The decrease in serum Urea occurs only in a few situations such as acute liver failure, starvation, in the last trimester of pregnancy.

BIBLIOGRAPHIC REFERENCES

- 1 - BALLETER, W. D. E.; BUSITMAN, C. S.; TIDWELL, P.W., Anal. Chim., 33-59.
- 2 - WINDMANN, F. K.; TURNER, K., Clin. Chem, 1987, 21:1754-1770.
- 3 - BERGMAYER, HU., Methods of Enzymatic Analysis, vol. 9, VCH Publishers, 1985, 449-453.
- 4 - TALKE, H.; SCHUBERT, G. E., Klin, 1965, 43-174.
- 5 - CARL, A. B. and EDWARD R. A., Tietz Textbook of Clinical Chem. 2a Ed., 1994, 1528-1531.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Urea UV kit: 10269360068

Review: June/2014

UNIVERSAL SYMBOLOGY