Bioclin

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REF K093

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for quantitative determination of Lipase in human serum or plasma. Colorimetric test, for *in vitro* diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic Colorimetric

The Lipase chromogenic substrate 1,2 - O - Dilauril - Rac - Glycerol - 3 - Acid - Glutaric - (6 - Methylresorufin) - Ester is cleaved by the catalytic action of alkaline lipase forming 1,2 - O - Dilauril - Rac - Glycerol and an intermediate unstable, Glutaric Acid - (6 - Methylresorufin) Ester. This spontaneously decompose to form Glutaric Acid in alkaline solution and Methylresorufin. The intensity of red color is formed directly proportional to the Lipase activity and can be determined photometrically.

1,2 - O - Dilauril - Rac - Glycerol - 3 - Acid Glutaric - (6 - Methylresorufin) - Ester - Lipase to 1,2 - O - Dilauril - Rac - Glycerol + Glutaric - (6 Methylresorufin) Acid.

Glutaric Acid - (6 - Methylresorufin) - spontaneous decomposition → Glutaric Acid and Methylresorufin.

REAGENTS

Number 1 - Buffer Reagent - Store between 2 and 8°C. Contains: BICIN tampon (N, N - Bis (2 - Hydroxyethyl) - Glycine) 40 mmol/L (pH 8,0), Colipase 1,7 mg/L, Na - Deoxycholate 16,4 mmol/L, Taurodesoxicolato 3,4 mmol/L, Calcium Chloride 12,0 mmol/L, detergent and preservative.

Number 2 - Substrate - Store between 2 and 8°C. Contains: Tartrate tampon 1,5 mmol/L (pH 4,0), 1,2 - O - Dilauril - Rac - Glycerol - 3 - Glutaric Acid - (6 - Methylresorufin) - Ester 0,13 mmol/L, Taurodesoxicolato 3,4 mmol/L, detergent and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 20 mL	1 x 5 mL
2	2 x 20 mL	2 x 5 mL
3	3 x 20 mL	3 x 5 mL
4	4 x 20 mL	4 x 5 mL
5	5 x 20 mL	5 x 5 mL
6	1 x 40 mL	1 x 10 mL
7	2 x 40 mL	2 x 10 mL
8	4 x 40 mL	4 x 10 mL
9	2 x 40 mL	1 x 20 mL
10	4 x 40 mL	2 x 20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes and tips, test tubes, watch or stopwatch, Biocal and Biocontrol Bioclin. Can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transportation temperature should be between 2 and 8°C. The transport temperature between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in the cleansing of the material should be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Reagent N° 2 must be protected from light.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or heparinized plasma can be used. Serum and plasma are stable for up to 1 week between 2 and 30°C and 1 year when frozen at -20°C.

TECHNIQUE

For the use of Lipase, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol N and P Bioclin kits.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

PROCESS LIMITATIONS

Studies have been performed to determine the level of interference of different compounds: Conjugate Bilirubin and non conjugate do not produce interference up to 60 mg/dL; Hemoglobin does not produce significant interference up to 500 mg/dL; Lipids do not produce significant interference up to 1000 mg/dL; there is a small correlation between the Lipemia index and the concentration of Triglycerides.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitored by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The calibration kit must be carried on using the BIOCAL calibrator, which is traceable to the reference material BCR-694 (Community Bureau of Reference).

REFERENCE VALUES

Adults: < 60 U/L

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS, METHODOLOGICAL SPECIFICITY AND SENSITIVITY

Lipase Automated kit was compared with other method commercially available for dosage of Lipase. 15 analyzer were performed and the results were evaluated. The linear equation obtained was $Y = 1,0119X - 1,6264$ and correlation coefficient 0,9584. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	42,25	48,90	30,85
Standard Deviation (U/L)	1,33	1,68	1,35
Coefficient of Variation (%)	3,15	3,44	4,37

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	41,05	49,42	30,83
Standard Deviation (U/L)	1,31	0,53	0,08
Coefficient of Variation (%)	3,19	1,06	0,25

Sensitivity

Limit of detection is 3 U/L.

Linearity

Test is linear up to 300 U/L.

DIAGNOSTIC SIGNIFICANCE

Lipase is an enzyme produced almost exclusively by pancreas. In the presence of pancreatitis, the lipase is released at high concentrations in the blood, persisting for much longer than the Amylase, allowing confirmation of late diagnosis disease. The values of Lipase in pancreatitis can reach from 5 to 10 times the upper limit of reference, while the other diseases Lipase levels are always below 3 times the values of reference. Lipase increases in acute pancreatitis, pancreatitis disease (cancer of pancreas), obstruction of the pancreatic ducts by calculation or malignancies. Since the excretion of Lipase is via the kidneys, kidney failure in general, Lipase levels are high. In intestinal infarction, in acute cholangitis and in small bowel obstruction increased values of Lipase can also be observed.

The lipase may be reduced physiologically during pregnancy, in the course of tuberculosis and other infectious diseases.

BIBLIOGRAPHIC REFERENCES

- 1 - Steinberg WM, Goldstein SS, Davies ND et al. Diagnostics assays in acute pancreatitis [Review]. Ann Intern Med 1985; 102:576-580.
- 2 - Panteghini M et al. Diagnostic value of four assays for lipase determination in serum: A comparative reevaluation. Clin Biochem 1991; 24:497-503.
- 3 - TIETZ, N. W., Chlinical Guide to Laboratory Tests., 3ª edição. Philadelphia, PA: WB Saunders, 1995: 865.
- 4 - Neumann U et al. New substrates for the optical determination of lipase. EP 207252 (1987).
- 5 - Gargouri Y, Julien R, Bois A et al. Studies on the detergent inhibition of pancreatic lipase activity. J of Lipase Research 1983; 24:1336-1342.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium


CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA resgistration for Lipase Automated kit: 10269360158

Review: August/2013

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED