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REF K131

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for the quantitative determination of Potassium in human serum. Enzymatic test, for *In vitro* diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic.

The Potassium is spectrophotometrically determinate through a kinetic assay using the Pyruvate Kinase dependent of potassium. The Pyruvate generated is converted into Lactate accompanying a conversion of NADH to NAD. The correspondent reduction of the optical density of 380nm is proportional to the concentration of the potassium in the serum.

REAGENTS

Reagent N° 1: Buffer - Store between 2 and 8°C. Contains: LDH, NADH, substrate and preservative.

Reagent N° 2: Enzymatic Reagent - Store between 2 and 8°C. Contains: Pyruvate Kinase and preservative.

Reagent N° 3: Standard Level 1 - Store between 2 and 8°C. Contains: Potassium Chloride (concentration varies to each lot)

Reagent N° 4: Standard Level 2 - Store between 2 and 8°C. Contains: Potassium Chloride (concentration varies to each lot)

PRESENTATION

Presenta-tion	Reagent N° 1	Reagent N° 2	Reagent N° 3	Reagent N° 4
5	40 mL	10 mL	3 mL	3 mL
6	2 x 40 mL	2 x 10 mL	3 mL	3 mL
7	3 x 40 mL	3 x 10 mL	3 mL	3 mL
8	4 x 40 mL	4 x 10 mL	3 mL	3 mL
9	30 mL	10 mL	3 mL	3 mL
10	2 x 30 mL	2 x 10 mL	3 mL	3 mL
11	3 x 30 mL	3 x 10 mL	3 mL	3 mL
12	4 x 30 mL	4 x 10 mL	3 mL	3 mL
13	60 mL	20 mL	3 mL	3 mL
14	2 x 60 mL	2 x 20 mL	3 mL	3 mL
15	3 x 60 mL	3 x 20 mL	3 mL	3 mL
16	4 x 60 mL	4 x 20 mL	3 mL	3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Although the substrate contains preservative, must avoid bacterial contamination.

5- Do not ingest. Avoid contact with eyes and skin.

6- The samples must be treated as potentially infectious.

7- The potassium controls are human serum. They were tested and presented negative results for HIV-1, HIV-2, HBV and HCV. However, the controls must be treated as a potentially infectious material.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.
11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum.

Do not use hemolyzed serum. The analyte is stable for 4 days between 2 and 8°C and 14 days at -20°C.

PROCESS DESCRIPTION

REAGENT PREPARATION

The reagents are ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

The kit is indicated for biochemistry automated analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

RESULTS

The results are expressed in mmol/L.

PROCEDURE LIMITATIONS

Periodically test the calibration to detect possible changes in the response of the instrument.

INTERFERENCE

The test does not suffer interference from the substances with the top concentrations of:

Na⁺ 150 mmol/L

NH₄⁺ 0,5 mmol/L

Ca²⁺ 7,5 mmol/L

P_i 2,0 mmol/L

Ascorbic Acid 10,0 mmol/L

Zn²⁺ 0,5 mmol/L

Fe³⁺ 0,5mmol/L

Cu²⁺ 0,5 mmol/L

Triglycerides 1000 mg/dL

Hemoglobin 500 mg/dL

Bilirubin conjugated 20 mg/dL

Bilirubin non conjugated 15 mg/dL

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standards are traceable to the reference material NIST (National Institute of Standards and Technology) SRM 918B.

REFERENCE VALUES

The reference values in mmol/L, for this method, were obtained through the determination of potassium in healthy populations of male and female.

Serum..... 3,5 a 5,1 mmol/L
13,7 - 19,9 mg/dL

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE**QUALITY CONTROL****Accuracy****COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY**

The Potassium Ions Automated was compared with another method commercially available to measure Potassium. 07 analysis were conducted and the results were evaluated. The linear equation was obtained $Y = 1,001X - 0,008$, with a correlation coefficient equal to 0,999. With these results we can conclude that the kit has good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mmol/L)	3,019	4,150	5,719
Standard Deviation (mmol/L)	0,009	0,008	0,009
Coefficient of Variation (%)	0,282	0,183	0,153

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mmol/L)	3,020	4,151	5,721
Standard Deviation (mmol/L)	0,002	0,001	0,002
Coefficient of Variation (%)	0,051	0,024	0,035

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Potassium. The average found was of 0,863 mmol/L with standard deviation of 0,004 mmol/L. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,876 mmol/L.

Linearity

The Potassium Ions Automated kit is linear between 2 and 8 mmol/L.

DIAGNOSTIC SIGNIFICANCE

In health subjects, the level of the liquid extra cellular potassium is regulated to be maintained at 3,5 - 5,1 mmol/L. Small deviation in the normal levels may cause severe consequences in one's health. The monitoring of the potassium serum concentration is important to both routine and emergency routines.

BIBLIOGRAPHIC REFERENCES

1. Wu, A.H.B., ed. Tietz clinical guide to laboratory tests, 4th edition, p. 880. W.B. Saunders Company, St. Louis (2006).
2. Bergmeyer, H.U., Gawehn, K., and Grassl, M. (1974) in *Methods of Enzymatic Analysis*. Second Edition, Volume I, 509-510, Academic Press, Inc., New York.
3. M.N. Berry, R. D. Mazzachi, M. Pejakovic, and M. J. Peake Enzymatic Determination of Potassium in Serum. *CLIN. CHEM.* 35/5, 817-820 (1989).

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Potassium Ions Automated kit: 10269360209

Review: September/2015

UNIVERSAL SYMBOLOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED