GLICOSE MONOREAGENTE

REF K082

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of glucose. Colorimetric enzymatic test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Colorimetric Enzymatic – GOD - PAP

The Glucose is enzymatically oxidized by Oxidase Glucose (GOD) according to the following reaction:

The Peroxide of Hydrogen in the presence of Peroxidase (POD) reacts with 4 - Aminoantipyrine and Phenol, forming a cherry red chromogen which color intensity is proportional to the glucose concentration.

REAGENTS

Reagent N° 1 - Enzymatic Reagent - Store between 2 and 8°C. Contains: Tampon < 36 mmol/L, Phenol < 20 mmol/L, 4 - Aminoantipyrin < 5 mmol/L, Oxidase Glucose > 10.000 U/L, Peroxidase > 700 U/L, stabilizer, surfactant and preservative.
Reagent N° 2 - Standard - Store between 2 and 8°C, well sealed. Contains: Glucose 100,0 mg/dL (5,56 mmol/L) and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	250 mL	3 mL
2	2 x 250 mL	3 mL
3	4 x 250 mL	3 mL
4	5 x 20 mL	3 mL
5	5 x 40 mL	3 mL
6	10 x 40 mL	3 mL
7	2 x 60 mL	3 mL
8	4 x 60 mL	3 mL
9	6 x 60 mL	3 mL
10	8 x 60 mL	3 mL
11	10 x 60 mL	3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watches and stopwatches, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found in the specialized market of clinical analysis laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8°C. The transport temperature between 15 and 30°C, must not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Water level in water bath should be above the level of reagents in test tubes.
- Handle the reagents with care. Reagent N° 1, containing sodium azide, irritating to skin and mucous membranes.
- The development of a pinkish color in reagent N° 1 do not interfere in reagents quality and stability, once the correspondent Blank be used and dosages in regular bases.
- Determinate the factor periodically and in every batch.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Plasma (fluoride), serum, liquor, liquid (ascilic, pleural and synovial).
 The use of Bioclin's anticoagulant Fluoride is recommended because it is a glucoses inhibitor. Use 1 drop for each 3 mL of blood. The serum can only be used in case it is centrifuged and dosed immediately after collection. In other biological liquids add an glucose inhibitor in the same proportion described for the blood, centrifuging the sample before the dosage begins.

PROCESS DESCRIPTION

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 tests tubes with the letters B (Blank), A (Sample), P (Standard) and proceed as follows:

	Blank	Standard	Sample
Sample	--	--	10 µL
Reagent N° 2	--	10 µL	--
Reagent N° 1	1,0 mL	1,0 mL	1,0 mL

Homogenize well and place in water-bath at 37°C for 10 minutes.

Read the sample and standard absorbance at 505 nm (490 - 550 nm), matching the zero with the blank. The color is stable for 30 minutes.

CALCULATIONS

$$\text{Glucose (mg/dL)} = \frac{\text{Sample Absorbance} \times 100}{\text{Standard Absorbance}}$$

Example:

$$\begin{aligned} \text{Sample Absorbance} &= 0,347 \\ \text{Standard Absorbance} &= 0,350 \end{aligned}$$

$$\text{Glucose (mg/dL)} = \frac{0,347 \times 100}{0,350} = 99$$

The reaction follows the Lambert-Beer Law, the calibration factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (100 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{mg/dL} = \text{Absorbance Sample} \times \text{Calibration Factor}$$

Example:

$$\text{Calibration Factor} = \frac{100}{0,350} = 286$$

$$\text{Glucose (mg/dL)} = 0,347 \times 286 = 99$$

The results are express in mg/dL.

PROCEDURE LIMITATIONS

The proposed method is not indicated for glucose dosage in urine.

INTERFERENCES

Samples with concentration up to 20 mg/dL of Bilirubin, 750 mg/dL of triglycerides and 160 mg/dL of Hemoglobin do not produce significant interference.

In case of interferences produced by the samples, also carry out the sample of Blank, in order to minimize the action of interferences. Proceed as follows:

Mark 1 tube with Sample of Blank and place 1,0 mL of Sodium Chloride 0,85% with 10 µL of sample. Determine the absorbance at 490 - 550 nm, matching the zero with distilled or deionized water. Subtract the obtained absorbance, of the Sample tube absorbance. Calculate the concentration multiplying the result of Calibration Factor.

The use of highly reductive medicine as such as Ascorbic Acid (Vitamin C) interfere in the reaction, because they compete with the consumption of H₂O₂, supplying values falsely lowered. For this reason, must suspend its use per less than 12 hours before the sample collection.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 917.

REFERENCE VALUES

The reference values in mg/dL for this method were obtained through the determination of glucose in healthy populations of male and female.

Plasma	Premature	20 - 60 mg/dL
	0 to 1 day	40 - 60 mg/dL
	> 1 day	50 - 80 mg/dL
	Children and adults	65 - 99 mg/dL
Liquor		50 - 70 mg/dL

To convert the values of mg/dL into mmol/L (SI), multiply the obtained results by 0,0556.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy RECOVERY

The recovery analysis was made with 5 sample determinations. The accurate was calculated and was found in a good parity with the reference values, obtaining a recovery between 94% and 104%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

Glucose Monoreagent kit was compared with other commercially available methods for measurement of glucose. For this test, 42 analyses are performed and the results were evaluated. The linear equation obtained was Y = 1,005X + 0,465 and the correlation coefficient 0,998. With these results, we can conclude that the kit shows good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	240,75	115,01	73,59
Standard Deviation (mg/dL)	3,72	4,15	1,92
Coefficient of variation (%)	1,55	3,61	2,61

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	242,29	114,58	74,06
Standard Deviation (mg/dL)	2,56	0,60	0,56
Coefficient of variation (%)	1,06	0,52	0,76

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of glucose. An average value of 0,55 mg/dL, with standard deviation of 0,25 mg/dL. The sensitivity that indicates the method detection limit corresponds the average plus three times the standard deviation and is equal to 1,31 mg/dL.

Linearity

The reaction is linear up to 500 mg/dL. For samples with values above 500 mg/dL or optical density above 0,8, dilute the sample with Sodium Chloride 0,85%, repeat the dosage and multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The glycemc homeostasis is controlled by the action of several hormones, especially insulin, who keeps the equilibrium of the glucose concentration. Hormonal alterations and other factors generate a variation in the homeostasis, leading to a hyper and hypoglycemic.

The Hyperglycemia occurs in several types of Mellitus Diabetes, where frequently are retinopathy, kidney injuries, neuropathies and arteriosclerosis. The Mellitus Diabetes is classified in: Mellitus diabetes insulin dependent (Type I), Mellitus diabetes non insulin dependent (Type II), Mellitus diabetes associated at certain conditions and syndromes (classified previously as secondary diabetes) and Pregnancy diabetes.

The Hypoglycemic (HG) levels of glycemia that leads to its manifestations are extremely variable. The manifestations can occur during fasting or postprandial. Occurs fasting HG in insulinome, non pancreatic tumors, hepatic diseases, hypoadrenalin (Addison's disease), hypopituitarism (Simmond's disease), storage of infirmity of glycogenic stored (Von Gierke's diseases). The postprandial

hypoglycemia occurs due to the reactive cause (symptoms of HG 1 to 3 hours after the meals), it can also be of alimental origin or in consequence of mellitus diabetes Type II and the intolerance Glucose.

BIBLIOGRAPHIC REFERENCES

- 1 - TRINDER, P., Ann. Clin. Biochem., 6-24, 1969.
- 2 - TONKS, D. B., Quality Control in Clinical Laboratories, 1983.
- 3 - LOTT, J. A.; TURNER, K., Clin. Chem., 1987, 21:1754, 1970.
- 4 - CARL, A. B. and EDWARD R. A., Tietz Textbook of Clinical Chem. 2nd ed., 928-997,
- 5 - Quibasa: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium


CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Glucose Monoreagent kit:
10269360136

Review: August/2014

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED