## COLINESTERASE

REF K094

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#### FINALIDADE

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- 1 - Recommendations of the German Society for Clinical Chemistry. Standardization of methods for the estimation of enzyme activities in biological fluids: Standard method for the determination of Cholinesterase activity. J Clin Chem Clin Biochem 1992; 30: 163-170.
- 2 - Thomas L, Clinical laboratory diagnostics. 1st ed Frankfurt: THBooks Verlagsgesellschaft; 1998. p-65-71.

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for determination of activity of Cholinesterase. Kinetic enzymatic test, only for *in vitro* diagnostic use.

#### PRINCIPLE OF ACTION

**Methodology:** Kinetic Enzymatic (DGKC).

On the catalytic action of Cholinesterase, the Butyrylthiocholine is hydrolyzed to Thiocholine and Butyrate:

The Thiocholine reduces the Hexacyanoferrate (III), yellow, Hexacyanoferrate (II), colorless. The decrease in absorbance is measured at 405 nm.

#### REAGENTS

**Number 1 - Buffer Reagent** - Store between 2 and 8°C. Contains: Pyrophosphate (pH 7,60), Potassium Hexacyanoferrate (III) and preservative.

**Number 2 - Substrate** - Store between 2 and 8°C. Contains: Butyrylthiocholine and preservative.

#### PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 20 mL	1 x 5 mL
2	2 x 20 mL	2 x 5 mL
3	1 x 40 mL	1 x 10 mL
4	2 x 40 mL	2 x 10 mL
5	2 x 40 mL	1 x 20 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated Spectrophotometer, pipettes, tips, test tubes, watch or stopwatch, can be found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

6- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

7- Do not use the product in case of damaged packaging.

8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum or heparinized plasma or EDTA can be used. Serum and plasma are for up to 1 week if kept at room temperature, 2 weeks if kept between 2 and 8°C and 6 months if frozen.

#### PROCESS DESCRIPTION

#### REACTION CONDITIONS

It is indispensable condition the use of thermostated cuvette at 37°C with 1cm optical path and 405 nm reading.

#### TECHNIQUE

	Blank	Sample
Sample	-	20 µL
Distilled Water	20 µL	-
Reagent 1	1000 µL	1000 µL
Homogenize and incubate for approximately 3 minutes at 37°C.		
Reagent 2	250 µL	250 µL

Homogenize again, incubate for exactly 2 minutes at 37°C and read the absorbance (405nm). Measure the change in absorbance per minute ( $\Delta A/\text{min}$ ) after 1, 2 and 3 exact minutes.

$$\Delta A/\text{min} = [\Delta A/\text{min Sample}] - [\Delta A/\text{min Blank}]$$

#### CALCULATIONS

Considering the optical path of 1 cm and reaction temperature 37°C the calculation will be:

$$\text{Activity (U/L)} = \Delta A/\text{min} \times 68500$$

#### PROCEDURE LIMITATIONS

Does not produce interference significant: Bilirubin up to 40 mg/dL, Hemoglobin up to 1000 mg/dL and lipemia up to 1400 mg/dL Triglyceride equivalents.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### REFERENCE VALUES

Men	4620 - 11500 U/L
Women	3930 - 10800 U/L

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

#### PRODUCT PERFORMANCE

##### QUALITY CONTROL

##### Accuracy

##### COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Cholinesterase kit was compared with other method commercially available to measure Cholinesterase. 07 analyzes were conducted and the results were evaluated. The linear equation obtained was  $Y = 1,0093X - 54,406$  and correlation coefficient 0,9998. With these results we can conclude the kit shows good methodological specificity.

##### Precision

##### REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	3643,00	4262,00	5323,50
Standard Deviation (U/L)	3,32	3,66	2,82
Coefficient of Variation (%)	0,09	0,09	0,05

**REPRODUCIBILITY**

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (U/L)</b>	1635,95	4483,80	7722,35
<b>Standard Deviation (U/L)</b>	3,58	3,11	3,87
<b>Coefficient of Variation (%)</b>	0,22	0,07	0,05

**Sensitivity**

Minimum limit of detection is 50 U/L.

**Linearity**

Linear test at 25000 U/L.

**DIAGNOSTIC SIGNIFICANCE**

(ChE) Serum Cholinesterase (pseudocholinesterase II or cholinesterase, EC3.1.1.8) is found primarily in the liver but also in white matter of the brain, pancreas, heart and serum. ChE serum levels may be determined as an indicator of hepatic poisoning by inhalation or skin contact by some organophosphorus compounds (including some insecticides neurological or gases) in liver diseases, or before of anesthesia with succinylcholine thus eliminating a deficiency congenita these enzymes that could lead to an apnea prolonged due to slow degradation of miorelaxant. A decrease of 15 - 25% is observed in mild poisoning, a decrease of 25 - 35% is observed in moderate poisoning, a decrease of 30 - 50% is observed in severe poisoning, but is also observed in patients with severe or chronic hepatitis. A decrease of 50 - 70% occurs in patients with an advanced cirrhosis and carcinoma liver metastasis.

**BIBLIOGRAPHIC REFERENCES**

- 1 - Recommendations of the German Society for Clinical Chemistry. Standardization of methods for the estimation of enzyme activities in biological fluids: Standart method for the determination of Cholinesterase activity. J Clin Chem Clin Biochem 1992; 30: 163-170.
- 2 - Thomas L, Clinical laboratory diagnostics. 1st ed Frankfurt: THBooks Verlagsgesellschaft; 1998. p-65-71.
- 3 - Hallbach J, Klinische Chemie fur den Einstieg. 1st ed Stuttgart: Thieme ; 2001. p. 143-144.
- 4 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

**OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium

**CUSTOMER SERVICE**

Customer Advisory Service  
Phone: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for Cholinesterase kit: 10269360159

Review: August/2015

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED