FINALIDADE

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of the MB fraction Creatine Kinase (CK). Kinetic test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Kinetic UV

The process involves the measurement of the activity of CK in the presence of an antibody against fraction M. This antibody completely inhibits the activity of CK-MM and M fraction of CK-MB, without, however, affect the activity of the B subunit of CK-MB and CK-BB. Based on the premise that the BB dimer is practically nonexistent in the peripheral blood, the residual enzyme activity found will be corresponding only to the fraction B of CK-MB. As the M fraction and fraction B have similar activities, we have the activity of CK-MB corresponds to the value found multiplied by two.

CK - Creatine Kinase

G-6-PDH - Glucose-6-Phosphate-Dehydrogenase

HK - Hexokinase

The speed of reduction of NADP⁺ to NADPH is proportional to the activity of CK in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Imidazole Acetate (pH 6,7) 100 mmol/L, Glucose 20 mmol/L, EDTA 2 mmol/L, NADP⁺, 2 mmol/L, Hexokinase 3500 U/L, Magnesium Acetate 10 mmol/L, N-Acetylcysteine 20 mmol/L, Polyclonal Anti CK-M sufficient to inhibit up to 2,000 U/L CK-MM.

Number 2 - Enzyme - Substrate - Store between 2 and 8°C. Contains: Glucose-6-Phosphate-Dehydrogenase 2000 U/L, Creatine Phosphate 30 mmol/L, ADP 2 mmol/L, AMP 5 mmol/L, Diadenosine Pentaphosphate 10 mmol/L.

PRESENTATION

Reagent	Presentation				
	1	2	3	4	5
N°1	1 x 40 mL	2 x 40 mL	4 x 40 mL	2 x 40 mL	4 x 40 mL
N°2	1 x 10 mL	2 x 10 mL	4 x 10 mL	1 x 20 mL	2 x 20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatch, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers, that may alter the results significantly.

5- It is important for the proper performance of the test, a strict control of time, temperature and pH.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma collected with heparin or EDTA. The sample is stable for 24 hours between 15 and 30°C and a week between 2 and 8°C. Do not use hemolyzed sample.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENT

Mix 4 parts of reagent N°1 with 1 part of reagent N°2. Stable for 1 day at temperatures between 15 to 30°C and 15 days between 2 to 8°C.

REACTION CONDITIONS

It is indispensable the use of a thermostated cuvette at 37°C, optical path of 1cm and reading at 340 nm.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

Working Reagent	1,0 mL
Sample	40 µL

Homogenize and transfer to a thermostated cuvette at 37°C and wait for 5 minutes. Make the initial reading, triggering the stopwatch. Repeat the readings at 1, 2 and 3 minutes.

CALCULATIONS

Calculate the mean difference in absorbance per minute ($\Delta A/\text{min}$.) and use this value to calculate the result.

CK-MB (U/L) 340 nm = $\Delta A/\text{min} \times 8254$

TECHNICAL FOR CK-MB GREATER THAN 20% OF TOTAL CK VALUE

If the measurement of activity of CK-MB exceeds 20% of total CK activity, the presence of macro CK should be suspected and the following technique to determine the value of CK-MB real should be performed.

Example:

CK Total = 150 U/L

20% CK Total = 30 U/L

If CK-MB > 30 U/L, performing technique below.

- Make the CK-NAC (total) and CK-MB before heating the serum.
- Heating the serum at 45°C for 20 minutes.
- Test the CK-MB again.

Value of CK-MB (U/L) = initial CK-MB - CK-MB*

CK-MB* = value found after heating the serum

Note: See Procedure Limitations

PROCEDURE LIMITATIONS

The method quantifies all CK-BB isoenzyme in serum. Normally the activity of this isoenzyme is small. However, if a significant increase of this isoenzyme occurs, the activity of CK-MB is overrated. Moreover, macromolecular two isoenzymes have been considered as false positive for CK-MB: a macro CK type 1 and type 2 macro CK. The type 1 is a complex of CK-BB or CK-MM bound to IgA or IgG. The macro CK type 2 is an oligomeric complex of mitochondrial origin and is associated with neoplasias. If the measurement of activity of CK-MB exceeds 20% of total CK activity, one should suspect the presence of these isoenzymes macromolecules.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The calibration kit can be made using the theoretical calibration factor based on the molar absorptivity of the NADP, or through the calibrator BIOCAL. The Bioclin recommends using BIOCAL calibrator which is traceable to the reference material ERM-AD455 / IFCC and the IFCC reference method.

REFERENCE VALUES

The reference values for this method were obtained by determining CK-MB in healthy populations of male and female.

37°C	
Men/Women	0 - 25 U/L

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE**QUALITY CONTROL****Accuracy****METHOD COMPARISON AND METHODOLOGICAL SPECIFICITY**

The CK MB UV Kit was compared with other methods for CK-MB dosage available in the market. 42 analysis were performed and the results were evaluated. The linear equation obtained was $Y = 0,994X - 0,485$, with linear correlation coefficient equal to 0,996. With these results we can conclude that the kit shows good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	51,50	253,45	155,68
Standard Deviation (U/L)	1,13	3,13	1,79
Coefficient of Variation (%)	2,2	1,2	1,1

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	52,21	247,12	155,25
Standard Deviation (U/L)	1,40	5,76	2,49
Coefficient of Variation (%)	2,69	2,33	1,61

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of CK-MB. The average found was 0,78 U/L with standard deviation of 0,52 U/L. The sensitivity, that indicates the method detection limit, correspond to the average plus 3 times the standard deviation, and is equal to 2,33 U/L.

Linearity

The method is linear up to 600 U/L. If the variation of molar extinction per minute is larger than 0,150, repeat the test using a dilution of 1:10 with saline. The result obtained is multiplied by 10.

DIAGNOSTIC SIGNIFICANCE

The CK-MB elevations are suffering myocardial diseases, such as myocardial infarction, myocarditis. Then the infarct (3 to 5 hours), the serum concentration of this enzyme rises, reaches a peak after 12 - 24 hours and remains elevated for a short period (2 to 3 days). As this is a specific fraction to the heart muscle, this enzyme has great significance for monitoring of myocardial infarction.

BIBLIOGRAPHIC REFERENCES

- 1 - SCANDINAVIAN Socit For Clinical Chemistry Scand. J. Clin. Lab. INVEST., 1974, 33, 291.
- 2 - Ann. Biol. Clin., 1982, 40,99.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for CK MB UV kit: 10269360073

Review: October/2014

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED