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USAGE INSTRUCTIONS


FUNCTION

Method for quantitative determination of Latent Binding Capacity of the Iron in serum or plasma samples. Colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Colorimetric test

A known concentration of Iron ions is incubated with serum and will bind specifically with the Transferrin in insaturated Iron-binding sites. The exceeding Iron ions are measured through the reaction with Ferene. The difference between the exceeding Iron and the total quantity of iron initially added to the serum is equivalent to the quantity bounded to the Transferrin, which will be the Latent Binding Capacity of the Iron – CLLF.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Buffer (pH 8,7), Iron Sulphate II and Tyurea.

Number 2 – Color Reagent - Store between 2 and 8°C. Contains: Ascorbic Acid, Ferene and Tyurea.

Reagent N° 3 - Calibrator - Store between 2 and 8°C. Contains: Serum proteins and 0.9% Sodium Azide.

Potentially Infectious.

Note: The concentration of serum proteins varies with the batch.

Reagent N° 4 - Control - Store between 2 and 8°C. Contains: Serum proteins and 0.9% Sodium Azide.

Potentially Infectious.

Note: The concentration of serum proteins varies with the batch.

PRESENTATION

Presentation	Reag. 1	Reag. 2	Reag. 3	Reag. 4
1	1 x 20 mL	1 x 5 mL	1 mL	1 mL
2	2 x 20 mL	2 x 5 mL	1 mL	1 mL
3	3 x 20 mL	3 x 5 mL	1 mL	1 mL
4	4 x 20 mL	4 x 5 mL	1 mL	1 mL
5	5 x 20 mL	5 x 5 mL	1 mL	1 mL
6	1 x 40 mL	1 x 10 mL	1 mL	1 mL
7	2 x 40 mL	2 x 10 mL	1 mL	1 mL
8	4 x 40 mL	4 x 10 mL	1 mL	1 mL
9	2 x 40 mL	1 x 20 mL	1 mL	1 mL
10	4 x 40 mL	2 x 20 mL	1 mL	1 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry auto-chemistry analyzers. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 2 and 8°C. The transport temperature between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- For professional *in vitro* diagnostic use only.
- Reagent n° 2 must be kept protected from light.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Calibration should be repeated periodically to verify any changes in the response of the colorimeter or spectrophotometer.
- Hemolysis, even discrete, interferes in the dosage.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Heparinized serum or plasma free of hemolysis. The sample is stable for 8 months between 2 - 8°C, protected from light, and 4 months between 20 - 25°C. In dosages for therapeutic control the samples must be collected always in the same time.

PROCESS DESCRIPTION

REAGENT PREPARATION

The reagents are ready to use.

TECHNIQUE

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

CALCULATIONS

The obtained value corresponds to the CLLF.

CLLF - Latent Iron-Binding Capacity

CTLF - Total Iron-Binding Capacity

IST - Transferrin Saturation Index

CTLF = CLLF + Iron Serum Concentration

$$\text{IST}(\%) = \frac{\text{Serum Iron}}{\text{CTLF}} \times 100$$

INTERFERENCES

There was no sign of interference of Ascorbic Acid up to 30 mg/dL, conjugated and free Bilirubin up to 60 mg/dL, lipemia up to 2000 mg/dL of Triglycerides, RF up to 350 UI/mL, Copper up to 15 mg/dL and Zinc up to 15 mg/dL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

The reference values for this method were obtained through the determination of the Iron-Binding Capacity in healthy populations of male and female.

Latent Capacity (CLLF)	140 - 280 $\mu\text{g}/\text{dL}$
Total Binding Capacity (CTLF)	250 - 410 $\mu\text{g}/\text{dL}$
Transferrin Saturation (IST)	20 - 50%

To convert the values of $\mu\text{g}/\text{dL}$ to $\mu\text{mol}/\text{L}$ (SI), multiply the obtained results by 0,179.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS

Iron-Binding Capacity Automated kit was compared with other methods for determination of Iron-Binding Capacity dosage commercially available measurement. 07 analyzes were performed and the results were evaluated. The equation linear obtained was $Y = 1,0128X - 2,969$ and the correlation coefficient 0,9987. With these results, we can conclude that the kit shows good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration ($\mu\text{g/dL}$)	179,05	384,60	487,70
Standard Deviation ($\mu\text{g/dL}$)	1,50	4,03	3,60
Coefficient of Variation (%)	0,84	1,05	0,74

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration ($\mu\text{g/dL}$)	155,65	334,03	422,75
Standard Deviation ($\mu\text{g/dL}$)	1,45	3,07	3,80
Coefficient of Variation (%)	0,93	0,92	0,90

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Iron-Binding Capacity. The average found was 5,81 $\mu\text{g/dL}$, with deviation standard of 0,07 $\mu\text{g/dL}$. The sensitivity, which indicates the method limit detection, corresponds the average plus 3 times the standard deviation, and is equal to 6,025 $\mu\text{g/dL}$.

Linearity

The reaction is linear until the concentration of 750 $\mu\text{g/dL}$ or 135 $\mu\text{mol/L}$. For higher values, dilute the sample with NaCl solution at 0,85% and repeat the dosage. Multiply the results by the dilution factor applied.

DIAGNOSTIC SIGNIFICANCE

Several diseases involving Iron in the organism may be evaluated from the determination of the Iron Binding Capacity combined with the determination of Iron Serum. The sum of those two values represents the maximum iron concentration that the proteins are capable of binding. In a iron deficiency, the iron binding capacity is increased, and in severe or malignant inflammatory diseases this binding capacity is reduced.

BIBLIOGRAPHIC REFERENCES

1. FAIRBANKS, V.F.; KLEE, G.G. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
2. WICH, M.; PINGERRA, W.; LEHMANN, P. Clinical

aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.

3. GUDER, W.G.; ZAWTA, B. et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46.

4. DATI, F.; SHUMANN, G.; THOMAS, L.; AGUZZI, F.; BAUDNER, S.; BIENVENU, J. et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34:517-20.

5. THOMAS L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:273-5.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Iron-Binding Capacity
Automation kit: 10269360180

Review: October/2015

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED