## CÁLCIO ARSENAZO III

REF K051

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#### FINALIDADE

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### USAGE INSTRUCTIONS

#### FUNCTION

Quantitative determination of calcium in serum, plasma and urine sample by color reaction. Manual and Automatic application.

#### PRINCIPLE OF ACTION

**Methodology:** Colorimetric Endpoint - Arsenazo III

The calcium reacts with arsenazo III in acid environment forming the blue-colored complex, whose intensity is proportional the calcium concentration in the sample. The absorbance of the reaction product is measured in wavelengths between 600 and 680 nm.

#### REAGENTS

**Reagent Nº 1: Arsenazo** - Store between 2 and 8°C. Contains: Buffer < 100 mmol/L, Arsenazo III 0,2 mmol/L, 8-Hydroxyquinoline < 20 mmol/L, surfactante > 1% and preservative.

**Reagent Nº 2: Standard** - Store between 2 and 8°C. Contains: Calcium Carbonate 10 mg/dL, stabilizer and preservative.

#### PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2
1	50 mL	3 mL
2	2 x 50 mL	3 mL
3	5 x 20 mL	3 mL
4	5 x 40 mL	3 mL
5	10 x 40 mL	3 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Automated instruments and semi-automatic biochemistry, clock or stopwatch, pipettes, test tubes, a water bath at 37°C, Biocontrol N and Biocontrol P Bioclin. It lies in the specialized Articles of Clinical Laboratory Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For professional *in vitro* diagnostic use only.**

**2-** The validity of the reactants is printed on the labels of same. However, once opened, the reagent can suffer from contamination during handling and have stability decreased. Important to note changes (Color change, precipitate formation) that may indicate contamination.

**3-** Nº 2 reagent contains Sodium Azide which can be irritating to skin and mucous membranes. Handle with care.

**4-** To follow strictly the methodology proposed for obtaining accurate results.

**5-** The water used in material cleaning must be recent and free from contaminants.

**6-** Columns deionizer saturated releases alkaline water, several ions and oxidizing and reducing agents, which can significantly alter the results.

**7-** It is extremely important the thorough cleaning of all material (pipettes, tips, test tubes, cuvettes, etc.) to obtain accurate results. Using solution sulphochromic or aqueous solutions of HCl 50% or detergent nonionic. Rinse the material in running water and rinse with distilled/deionized water. Dry the material under heater protected from direct contact with the metal part of (put the material in glass vats).

**8-** Do not allow contact of the sample anticoagulant such as EDTA that is a chelat calcium, citrate that complex Ca<sup>2+</sup>, oxalate and fluoride that precipitate Calcium.

**9-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**10-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**11-** Do not use the product in case of damaged packaging.

**12-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

**Serum or Plasma (Heparin):** Free of hemolysis. Do not use oxalate or EDTA as an anticoagulant by interfering the determination of calcium.

**Urine:** Make the urine collection within 24 hours in a recipient free of calcium. Prior to the collection add into container 10 mL of Nitric acid 50% v/v. When the sample of 24 hours is not acidified during collection, add 20 mL of HCl 6M, homogenize, wait 60 minutes and take the aliquot for the test.

#### PROCESS DESCRIPTION

##### TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Select the wavelength from 600 nm to 680 nm.

Mark 3 test tubes: B (Blank), P (Standard) and A (Sample) and proceed as follows:

	Blank	Standard	Sample
Reagent Nº 1	1,0 mL	1,0 mL	1,0 mL
Reagent Nº 2	---	10 µL	---
Sample	---	---	10 µL

Homogenize and incubate for 2 minutes at 37°C.

Read the absorbance of standard and sample against the blank reactive. The color is stable for 1 hour.

#### CALCULATIONS

**Serum Calcium:**

$$\text{Calcium (mg/dL)} = \frac{\text{Test Absorbance}}{\text{Standard Absorbance}} \times 10$$

As the reaction follows the Beer-Lambert law, the Calibration Factor can be used.

$$\text{Calibration Factor} = \frac{10}{\text{Standard Absorbance}}$$

$$\text{Calcium (mg/dL)} = \text{Absorbance Test} \times \text{Factor}$$

#### Urinary Calcium:

$$\text{Urinary Calcium (mg/dL)} = \frac{\text{Urinary Calcium (mg/dL)} \times \text{Volume (mL)}}{\text{(mg/24h)} \quad 100}$$

#### Ionized Calcium (Cal):

$$\text{Cal (mg/dL)} = \frac{6 \times \text{Ca} - [(0,19 \times \text{P}) + \text{A}]}{(0,19 \times \text{P}) + \text{A} + 6}$$

Ca = Calcium Serum (mg/dL)

P = Total Proteins (g/dL)

A = Albumin (g/dL)

#### PROCEDURE LIMITATIONS

Do not use hemolyzed samples.

Icteric and lipemic samples may interfere with measurements.

Proceed to the dosage as follows:

1 - Dosing the sample normally following the described technique. Record the readings.

2 - Add 0,02 mL of the anticoagulant EDTA Bioclin to the blank and test it. Proceed with the new reading hitting the zero with the blank.

3 - Express the difference multiplying the difference between the two readings by the Calibration Factor.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 915B.

#### REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of calcium in healthy populations of male and female.

Calcium (Serum or Plasma)	
Umbilical Cord	8,2 - 11,2 mg/dL
Premature	6,2 - 11,0 mg/dL
0 - 10 days	7,6 - 10,4 mg/dL
Infants	9,0 - 11,0 mg/dL
2 - 12 years	8,8 - 10,8 mg/dL
Adults	8,8 - 11,0 mg/dL
Calcium Ion	
1 - 18 years	4,80 - 5,52 mg/dL
Adults	4,6 - 5,4 mg/dL
Calcium Urine	
Adults	Up to 180 mg/24h (With calcium-restricted diet: 500 mg/24h)
	Up to 280 mg/24h (Without calcium-restricted diet)

### Critical values

Calcium (serum or plasma): less than 6.0 mg/dL and higher to 14,0 mg/dL.

To convert the values to mg/dL to mmol/L (SI), multiply by 0,25. These values should be used as reference, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

### PRODUCT PERFORMANCE QUALITY CONTROL

#### Accuracy COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY.

The Calcium Arsenazo III kit was compared to other method for calcium dosage available in the market. 42 analysis were performed and the results were analyzed. The obtained linear equation was  $Y = 1,012X - 0,569$ , with linear correlation coefficient equal to 0,993. With these results we can conclude that the kit shows good methodological specificity.

#### Precision REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	6,28	11,89	9,33
Standard Deviation (mg/dL)	0,19	0,10	0,07
Variation Coefficient (%)	3,04	0,88	0,76

#### REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	6,38	11,85	9,21
Standard Deviation (mg/dL)	0,15	0,05	0,10
Variation Coefficient (%)	2,34	0,40	1,11

#### Sensitivity

The sensitivity was calculated from 40 determinations a sample free of Calcium. The average found was 0,022 mg/dL with a standard deviation of 0,018 mg/dL. The sensitivity, which indicates the detection limit of the method, corresponds the average plus 3 times the standard deviation and is equal to 0,074 mg/dL.

#### Linearity

The reaction is linear up to a concentration of 20 mg/dL. For higher values dilute the sample 1:2 or 1:4 with deionized or distilled water and repeat the determination. Multiply the result by the dilution factor.

#### DIAGNOSTIC SIGNIFICANCE

Calcium plays an important role in the composition of the skeleton and reactions as essential in the clotting blood, in maintaining the integrity and permeability cell membranes, the stimulation of the skeletal and heart muscles, in neuromuscular conduction. The hypercalcemia can be seen in hyperparathyroidism, malignant tumor, Paget's disease (osteitis deformans), hypervitaminosis D, in cases of immobilization for fractures. Lower serum calcium can be observed in hypoparathyroidism, renal failure, deficiency in vitamin D (rickets), acute pancreatitis, massive blood transfusion.

#### BIBLIOGRAPHIC REFERENCES

- CONNERTY, H. V.; BRIGGS, A. R., Amer J. Clin. Pathol., 1.966, 45 - 290.
- GITELLMAN, H. J., Anal. Biochem, 1.967, 18 - 521.
- HENRY, R. J., Clinical Chemistry Principles and Technics.
- JANSEN JW, Helbing AR Eur J Clin Chem Clin Biochem 1991;29:197- 201.
- LEARY NO, Pembroke A, Duggan PF. Clin Chem 1992;38: 904 - 8.
- MARTIN, E. W., Harzards of Medication, 1971.
- MORGAN BR, Artiss JD, Zak B. Clin Chem 1993;39:1608-12.
- POTTGEN P, Davis ER. Clin Chem 1976;22: 1752.
- STEM, J., LEWIS V. H. P., Clin. Chim. Acta , 1957, 2 - 576.
- TIETZ,N.W. Tietz Textbook of Clinical Chemistry, Burtis CA, Ashwood ER eds, 2a. edição, Philadelphia: W.B. Saunders Co., 1994.
- TODD; SANFORD; DAVIDSOHN, Diagnósticos Clínicos, 1983.
- TONKS, D. B., Clin. Chem., 1.963, 9 - 217.

#### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under the appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP **OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium


#### CUSTOMER SERVICE

Customer Advisory Service  
Phone: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for Calcium Arsenazo III kit: 10269360126

Review: November/2014

## UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED