BIOPROT U/LCR

REF K108

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Total Protein in urine and liquor samples. Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Colorimetric Test - Pyrogallol Red

The proteins present in urine react to Pyrogallol Red and the Molybdate forming a red color complex, with maximum absorption in 600 nm (580 - 620 nm). The intensity of the color is directly proportional to the protein concentration in the sample.

REAGENTS

Reagent N° 1 - Color Reagent - Store between 2 and 8°C. Contains: Pyrogallol Red < 60 µmol/L, Sodium Molybdate > 40 µmol/L, surfactant, solubilizer, stabilizer and preservative.

Reagent N° 2 - Standard - Store between 2 and 8°C. Contains: Albumin 65 mg/dL and preservative.

Reagent N° 3 - Control - Store between 2 and 8°C. Contains: Proteins* and preservative.

* Value varies with each batch.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2	Reagent N° 3
1	1 x 25 mL	1 x 3,0 mL	1 x 3,0 mL
2	2 x 25 mL	1 x 3,0 mL	1 x 3,0 mL
3	3 x 25 mL	1 x 3,0 mL	1 x 3,0 mL
4	4 x 25 mL	1 x 3,0 mL	1 x 3,0 mL
5	5 x 20 mL	1 x 3,0 mL	1 x 3,0 mL
6	5 x 40 mL	1 x 3,0 mL	1 x 3,0 mL
7	10 x 40 mL	1 x 3,0 mL	1 x 3,0 mL

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the spectrophotometer or colorimeter, watch and stopwatch, pipettes and test tubes. Can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 2 and 8°C. The transport temperature between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- The presence of detergent in the used material can alter the results significantly.

6- The calibration should be repeated periodically, to verify any alteration in the colorimeter or spectrophotometer response.

7- Reagents N° 1, 2 and 3 contain Sodium Azide, irritating to skin and mucous membrane. Handle with caution.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Data Sheet Safety of Chemicals) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Urine (aleatory sample)

24 Hours urine

Brainspinal fluid (liquor)

Preparation of samples:

Urine - There is no necessity of preservatives. Take an aliquot of 20 mL of urine, centrifuge and proceed the technique using the supernatant.

Liquor: Centrifuge the sample and use the supernatant.

In samples of urine stored between 2 and 8°C, the protein is stable for 2 days, 1 day between 20 and 25°C and 1 year at -20°C. In liquor samples, it is stable for 1 day between 20 and 25°C, 6 days between 4 and 8°C and 1 year at -20°C. In dosages for therapeutic control the samples must be collected always in the same time.

PROCESS DESCRIPTION

REAGENT PREPARATION

The reagents are ready to use.

TECHNIQUE

Mark 3 tests tubes with the letters B (Blank), S (Standard), A (Sample) and proceed as follows:

	B	S	A
Sample	--	--	50 µL
Standard	--	50 µL	--
Distilled Water	50 µL	--	--
Reagent N° 1	1000 µL	1000 µL	1000 µL

Homogenize and incubate for 5 minutes at 37°C and read the absorbance in a 600 nm filter (578 to 630 nm), against a Blank reagent in no more than 30 minutes.

CALCULATIONS

The determination of total protein concentration in urine and liquor is calculated by the reaction:

$$\text{Protein (mg/dL)} = \frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times 65$$

Example:

In a total protein dosage in liquor, the following absorbance was found:

$$\begin{aligned} \text{Sample Absorbance} &= 0,201 \\ \text{Standard Absorbance} &= 0,360 \\ \text{Standard Concentration} &= 65 \end{aligned}$$

$$\text{Protein (mg/dL)} = \frac{0,201}{0,360} \times 65 = 36,29 \text{ mg/dL}$$

According to the Lambert-Beer law, you can use the Calibration Factor:

$$\text{Calibration Factor} = \frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$$

Example:

$$\text{Cf} = \frac{65}{0,360} = 180,55$$

$$\begin{aligned} \text{Protein (mg/dL)} &= \text{Sample Absorbance} \times \text{Cf} \\ \text{Protein (mg/dL)} &= 0,201 \times 180,55 = 36,29 \end{aligned}$$

The results are express the mg/dL.

24 hours Urine

Proceed the 24 hours dosage using a aliquot of the total volume. This aliquot must be centrifuge. The concentration in mg/24h is obtained using the dosage result (mg/dL) multiplying this value by the total urine volume (mL) and divide the result by 100.

$$\begin{aligned} \text{24 Hour Urine Protein (mg/24h)} &= \\ \frac{\text{Aliquot dosage (mg/dL)} \times \text{24 Hours volume (mL)}}{100} \end{aligned}$$

Example:

In a dosage of total protein in 24 hour urine, the following result was obtained:

$$\begin{aligned} \text{Sample dosage (mg/dL)} &= 5,0 \\ \text{24 Hours Urine volume (mL)} &= 2000 \end{aligned}$$

$$\text{24 Hours Urine Protein (mg/24h)} = \frac{5,0 \times 2000}{100} = 100 \text{ mg/24 hours}$$

INTERFERENCES

Blood presence in urine or liquor produce wrongly elevated results. The presence of detergent in the used material can altered significantly the results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D.

REFERENCE VALUES

The reference values for this method were obtained through the determination of the Total Protein in healthy populations of male and female.

Random Urine	1 to 15 mg/dL
24 Hours Urine	24 to 141 mg/24h
Liquor	15 to 45 mg/dL

To convert the values of mg/dL to g/L (SI), multiply the obtained results by 0,01. These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE**QUALITY CONTROL****Accuracy****COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY**

The Bioprot U/CRL was compared with other commercially available method for dosage of Total Protein in urine. 41 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1,000X - 0,080$, and correlation coefficient 0,99999. With these results, we can conclude that the kit shows good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/24h)	32,84	15,25	27,97
Standard Deviation (mg/24h)	0,46	0,17	0,17
Coefficient of Variation (%)	1,40	1,13	0,59

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/24h)	31,96	15,44	27,89
Standard Deviation (mg/24h)	0,77	0,16	0,10
Coefficient of Variation (%)	2,42	1,06	0,37

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of Total Protein in urine. The average found was 1,484 mg/dL, with standard deviation of 0,093 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 1,764 mg/dL.

Linearity

The reaction is linear until the concentration of 100 mg/dL. For higher values, dilute the sample with NaCl solution at 0,85% and repeat the dosage. Multiply the results by the applied dilution factor.

DIAGNOSTIC SIGNIFICANCE

The proteinuria in high concentration is found in the majority of the kidney problems, such as primary and secondary nephropathies, that can cause the increase of the glomerular permeability or reduction of the kidney tubules absorption. It could also have post kidney causes, such as infections, bleedings or malignant urinary tract diseases. In other acute situations such as fever or physical or psychological stress, high levels in protein can also be found. In liquor, high levels of protein can be found in cases of increase intracranial pressure due to cerebral tumors, intracerebral hemorrhage or traumatism, in inflammations (mainly due to bacterial meningitis) and in case of multiple sclerosis. Increase the permeability of the barrier of cephalorachidian liquid-blood is reflected as the increase proportion of liquor/serum of total protein.

BIBLIOGRAPHIC REFERENCES

- Johnson AM, Rohlfis EM, Silverman LM . Proteinias. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3 ed. Philadelphia: W.B Saunders Company; 1999. p. 477-540.
- Felgenhauer K. Laboratory diagnosis of neurological diseases. In: Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1308-26.
- Boege F. Urinary proteins. In: Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 382-400
- Orsonneau JL, Douet P, Massoubre C, Lustenberger P, Bernard S. An improved pyrogallol red-molybdate method for determining total urinary protein. Clin Chem 1989;35:2233-6.
- Watanabe N, Kamei S, Ohkubo A, Yamanaka M, Ohsawa S, Makino K et al. Urinary protein as measured with a pyrogallol red-molybdate complex. Manually and in a Hitachi 726 automated analyzer. Clin Chem 1986; 32:1551-4.
- Bioclin – Dados de arquivo.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil


 **OBELIS S.A.**
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Bioprot U/CRL kit: 10269360178

Review: August/2014

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED