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REF. K106

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Total Bilirubin in samples serum or plasma. Colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Colorimetric Test

Total Bilirubin, through the coupling reaction with 2,4 Dichloroaniline Diazotised, forms an azocompound, a red colored complex with maximum absorption at 546 nm. The intensity of the color formed is directly proportional to the concentration of Total Bilirubin in the sample.

REAGENTS

Reagent N°1 Buffer - Store between 2 and 8°C. Contains: Phosphate Buffer 40 mmol/L, NaCl 9 g/L, detergent and stabilizer.

Reagent N°2 Color Reagent - Store between 2 and 8°C. Contains: 2,4 Dichlophenyl -Diazonium Salt 1 mmol/L, HCl 30 mmol/L and detergent.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 40 mL	1 x 10 mL
2	2 x 40 mL	2 x 10 mL
3	3 x 40 mL	3 x 10 mL
4	4 x 40 mL	4 x 10 mL
5	5 x 40 mL	5 x 10 mL
6	2 x 40 mL	1 x 20 mL
7	4 x 40 mL	2 x 20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the Biochemistry Automated Equipment, Biocal kit and Biocontrol Bioclin.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Reagent N° 2 must be kept protected from light.
- 3- Water used in cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Hemolyses, even if discreet, may interfere in dosage.

6- Reagent N° 1 contains Sodium Azide, irritant to skin and mucous. Handle reagent with care.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma free of hemolysis.

Sample must be protected from light.

In samples stored between 2 to 8°C and protected from light, bilirubin is stable for three (03) days.

Sample for therapeutic control should always be collected at the same hour.

PROCESS DESCRIPTION

PREPARATION OF REAGENTS

Reagents are ready for use.

TECHNIQUE

For the use of the kit, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol Kit Bioclin.

This kit is for use in automated biochemical analyzers.

Check the application sheet for this equipment on www.bioclin.com.br or through SAC.

INTERFERING

Hemolysis, even if slight, interfere with measurement of Bilirubin. No interference was observed for Ascorbic Acid up to 30 mg/dL and Lipemia up to 1500 mg/L from Triglycerides. According to bibliographic data, some drugs, as Anabolic Steroids, Ascorbic Acid, Salicylates and Vitamin A, can raise Bilirubin values. Caffeine and Penicillins may lead to falsely decreased.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory.

Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 916A.

REFERENCE VALUES

Reference values were obtained through determinations of Total Bilirubin in healthy populations of male and female.

Newborns	24 hs	< 8,8 mg/dL
	2nd day	1,3 - 11,3 mg/dL
	3rd day	0,7 - 12,7 mg/dL
	4th to 6th day	0,1 - 12,6 mg/dL
Children	> 1 month	0,2 - 1,0 mg/dL
Adults		0,1 - 1,2 mg/dL

To convert the values from mg/dL into mmol/L (SI), multiply the results obtained for 0,0171.

These values should be used only as a guide. Each laboratory should establish its own reference range, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

Bilirubin Total Automated kit was compared with other commercially available methods for measurement of Total Bilirubin. 42 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 0,999X + 0,018$ and the correlation coefficient 0,99998. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	1,44	11,90	6,96
Standard Deviation (mg/dL)	0,01	0,03	0,04
Coefficient of Variation (%)	0,91	0,24	0,64

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	1,44	11,96	6,91
Standard Deviation (mg/dL)	0,01	0,06	0,05
Coefficient of Variation (%)	0,86	0,49	0,66

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of Bilirubin Total. The average found was 0,023 mg/dL with a standard deviation of 0,005 mg/dL. The sensitivity, which indicates the method detection limit, correspond the average plus 3 times the standard deviation, is equal to 0,037 mg/dL.

Linearity

The reaction is linear up to the concentration of 30 mg/dL. For higher values dilute the sample with NaCl 0.85% and repeat the dose. Multiply the resulting number by the dilution factor used.

DIAGNOSTIC SIGNIFICANCE

Bilirubin is a breakdown product of hemoglobin in reticuloendothelial system. It is conjugated in the liver to be then excreted in bile.

The level of bilirubin is useful for differential diagnosis of hepatobiliary disease and other causes of jaundice which becomes clinically evident when the concentration of Total Bilirubin greater than 2,5 mg/dL.

The Direct Bilirubin rises in plasma in the presence of hereditary liver diseases such as disease Dubin-Johnson and Rotor, hepatocyte injury (viral, toxic or alcoholic), biliary obstruction (stones or neoplasm), hepatitis acute or chronic and toxic reactions to various drugs (such as Chlorpromazine, organic arsenic compounds and Methyltestosterone, among others). Levels of Direct Bilirubin over 50% of the total bilirubin concentration indicate the existence of post-liver.

Already the indirect Bilirubin elevates with the existence of hemolytic anemia, self-immune hemolysis, blood transfusion, bruising re-absorption, inefficient erythropoiesis, hereditary diseases, such as Gilbert and Crigler-Najar diseases, and in new born jaundice.

The use of drugs that activate the hepatic microsomal system can reduce the concentration of Total and Direct Bilirubin in plasma.

BIBLIOGRAPHIC REFERENCES

- 1 - Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:192-202.
- 2 - Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors, Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
- 3 - Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962; 6:570-8.
- 4 - Bioclin - Dados de arquivo.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium


CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Bilirubin Total Automated kit:
10269360175

Review: September/2014

UNIVERSAL SYMBOLLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED