ALBUMINA MONOREAGENTE

REF K040

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Albumin. Colorimetric test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Bromocresol Green (BCG)

The measurement uses what is called "protein error of indicators". In the presence of Albumin, Bromocresol Green forms a colored complex, which exhibits an absorption spectrum different from the dye in its free state, thus allowing the Albumin dosage.

REAGENTS

Number 1 - Color Reagent - Store between 2 and 8°C. Contains: Bromocresol Green 0,1 mmol/L, Citrate Buffer Solution (pH 3,6) 20 mmol/L, preservative and surfactant.

Number 2 - Standard - Store between 2 and 8°C. Agitate before use. Contains: Albumin 3,8 g/dL, Sodium Azide 15,38 mmol/L.

Presentation	Reagent N° 1	Reagent N° 2
1	250 mL	2 mL
2	500 mL	2 mL
3	5 x 40 mL	2 mL
4	10 x 40 mL	2 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, watch or stopwatch, pipettes, test tubes, Bioclin's Biocontrol N and Biocontrol P. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Reagent N° 2 contains Sodium Azide: handle with care.

6- Do not use plasma.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, The analyte is stable for 03 days between 2 and 8°C.

WORKING REAGENT PREPARATION

The reagent is ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 03 test tubes: B (Blank), P (Standard), A (Sample), and proceed as follows:

	Blank	Standard	Sample
Sample	--	--	10 µL
Standard	--	10 µL	--
Color Reagent	2,5 mL	2,5 mL	2,5 mL

Homogenize well and leave it at rest for 5 minutes. Read the absorbance of Sample and Standard at 630 nm (600 - 640 nm), hitting the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

$$\text{Albumin (g/dL)} = \frac{\text{Sample Absorbance} \times 3,8}{\text{Standard Absorbance}}$$

As the reaction follows the Beer-Lambert Law, the calibration factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (3,8 g/dL)}}{\text{Standard Absorbance}}$$

$$\text{g/dL} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

Results are expressed in g/dL.

PROCEDURE LIMITATIONS

Use an automated pipette and good quality tips for quality pipetting from Standard and Sample, to minimize problems with volume imprecision.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D.

REFERENCE VALUES

The reference values in g/dL, for this method, were obtained through the determination of Albumin in healthy populations.

Children and Teenagers	1 - 30 days	2,6 - 4,3 g/dL
	31 - 182 days	2,8 - 4,6 g/dL
	183 - 365 days	2,8 - 4,8 g/dL
	1 - 18 years	2,9 - 4,7 g/dL
Adults		3,5 - 5,5 g/dL

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Albumin Monoreagent kit was compared with another method commercially available to measure Albumin. 07 analysis were performed and the results were evaluated. The linear equation obtained was $Y = 0,998X + 0,077$ and the correlation coefficient 0,994. With these results we can conclude that the kit shows good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (g/L)	3,41	4,69	5,29
Standard Deviation (g/L)	0,07	0,07	0,08
Coefficient of Variation (%)	2,11	1,53	1,54

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (g/L)	3,41	4,69	5,31
Standard Deviation (g/L)	0,01	0,01	0,02
Coefficient of Variation (%)	0,29	0,21	0,29

Sensitivity

Sensitivity was calculated from 20 determinations of a sample free of Albumin. The average found was of 0,0326 g/dL with standard deviation of 0,0001 g/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,0327 g/dL.

Linearity

The reaction is linear up to 6,0 g/dL. For samples with values above 6,0 g/dL or optical density above 0,8, it is recommended to dilute it with Saline Solution at 0,85%, repeat the dosage and multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Albumin constitutes the main serum protein. Synthesized almost entirely by the liver, has a half-life of approximately two weeks. An increase of albumin may be observed in dehydration, shock and hemoconcentration. Lower values occur in malnutrition, nephrotic syndrome, hepatic failure, glomerulonephritis, multiple myeloma, anemia serious, pregnancy, severe infections and protracted.

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6 - TONKS, D. B., Clin. Chem., 1983, 9,217.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service


Phone.: 0800 0315454

Email: sac@bioclin.com.br

ANVISA registration for Albumin Monoreagent kit:
10269360114

Review: December/2013

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED