## MICROALBUMINÚRIA

REF K078

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#### FINALIDADE

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for quantitative determination of Microalbumin in human urine by turbidimetry analysis, for *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Immunoturbidimetry

Turbidimetric assay for quantification of Microalbumin in urine sample. The Microalbumin present in human urine sample reacts with a specific antibody. Turbidity is induced by the formation of immune complexes, which are directly proportional to its concentration in the sample.

#### REAGENTS

**Number 1: Buffer** - Store between 2 and 8°C. Contains: Saline-Phosphate buffer (pH 7,4), Polyethylene Glycol 40 g/L and Sodium Azide 0,95 g/L.

**Number 2: Antiserum** - Store between 2 and 8°C. Contains: Saline-Phosphate buffer (pH 7,4), Polyclonal Anti-serum Albumin and Sodium Azide 0,95 g/L.

**Number 3: Calibrator** - Store between 2 and 8°C. Contains: Albumin and Sodium Azide 0,9%.

**Warning: The concentration of the Calibrator varies with the lot - See the bottle label.**

**Number 4: Control** - Store between 2 and 8°C. Contains: Albumin and Sodium Azide 0,9%.

**Warning: The concentration of Control varies according to the lot - See the bottle label.**

#### PRESENTATION

Reagent	Presentations			
	1	2	3	4
N°1	25 mL	2 x 25 mL	3 x 25 mL	4 x 25 mL
N° 2	2 mL	2 x 2 mL	3 x 2 mL	4 x 2 mL
N° 3	1 mL	1 mL	1 mL	1 mL
N° 4	1 mL	1 mL	1 mL	1 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

This kit is specific to automatic equipment.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transportation, at temperatures between 15 and 30°C, should not exceed 72 (seventy two) hours. Keep away from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For professional *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain accurate results.

**3-** The water used in the material rinsing must be recent and deionized.

**4-** It is important for the proper performance of the test, a rigorous controlling time and temperature.

**5-** Do not freeze the reagents, the freezing of antibody or diluent can affect their performance.

**6-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**7-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**8-** Do not use the product in case of damaged packaging.

**9-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Turbid urine or with visible sediment must be centrifuged before use. Centrifuge for 10 minutes at 3000 rpm. Use the supernatant for testing.

The Microalbumin in urine is stable for 48 hours if refrigerated at 2 and 8°C. For longer periods, freeze the sample at -20°C for up to 3 months, and only once. All samples are considered potentially infectious, and therefore must be handled following the rules of biosafety.

#### WORKING REAGENT PREPARATION

Before preparing the reagent check the application sheet for the equipment.

In case of application sheet with only one reagent perform the preparation as described below. For bireagent application sheet is not necessary preparation of reagents.

If in doubt contact the SAC (0800 031 5454).

Dilute 1 portion of Antiserum (Reagent N° 2) with 12,5 portions of Buffer (Reagent N° 1).

Example: 1 mL of Antiserum + 12,5 mL of Buffer

The diluted Antiserum, stored in an ambar flask is stable during 7 days in temperatures between 2 and 8°C and 1 day in temperatures between 15 and 30°C.

#### TECHNIQUE

For calibration and control of the reaction, using calibrator and internal controls in the kit.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC.**

#### CALCULATIONS

##### 24 Hour Urine

For determination of Microalbumin in 24-hours urine should be performed dosage Microalbumin in the sample following the original technique of kit. The result found by the technique is given in mg/L. For expression of the result in mg/24 hours using the formula below.

$$\text{Microalbumin (mg/24horas)} = \frac{\text{Microalbumin (mg/L)} \times \text{Urinary Volume (L)}}{24 \text{ hours (L)*}}$$

\*To convert the urine volume of mL to L divide the value given in mL by 1000.

Example:

Results of Microalbumin found in urine samples of 24 hours = 56,0 mg/L

Urinary volume of 24 hours = 1,9 L

Microalbumin (mg/24 hours) = 56,0 x 1,9 = 106,4 mg/24 hours

#### Relationship Albumin/Creatinine (mg/g)

To determine the relationship between Microalbumin and urinary Creatinine due, firstly, to perform dosage Microalbumin and Creatinine in the same urine sample. A urine sample can be random or first morning.

The urine Creatinine values are given in mg/dL. To convert the obtained concentration in mg/dL to g/L, simply divide the result by 100. The result of Microalbumin and Creatinine ratio is given as mg/g.

Example:

Creatinine value found in sample = 120 mg/dL

Microalbumin value found in the sample (mg/L) = 60,0 mg/L

$$\text{Urinary Creatine (g/L) found} = \frac{\text{Value Creatinine}}{100}$$

$$\text{Urinary Creatine (g/L)} = \frac{120 \text{ mg/dL}}{100} = 1,2 \text{ g/L}$$

$$\text{Microalbumin / Creatinine (mg/g)} = \frac{\text{Microalbumin (mg/L)}}{\text{Creatinine (g/L)}}$$

$$\text{Microalbumin / Creatinine (mg/g)} = \frac{60,0 \text{ mg/L}}{1,2 \text{ g/L}} = 50,0 \text{ mg/g}$$

#### PROCEDURE LIMITATIONS

Clinical diagnosis should be made not only with results of one test should be considered at the same time, clinical data of patients.

#### INTERFERENCES

Triglycerides up to 2500 mg/dL, Bilirubin up to 0,1 mg/dL and Hemoglobin up to 125 mg/L did not significantly interfere with results. Some drugs and substances affect the concentration of Microalbumin in urine.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit's calibrator is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D.

#### REFERENCE VALUES

Reference values below should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

#### Random Urine or First Morning

Microalbumin: 0 - 25 mg/L (IFCC)

Relation Microalbumin and Urinary Creatinine (mg/g)		
Reference Value	Men (mg/g)	Women (mg/g)
Normal	< 17	< 25
High	17 - 250	25 - 355
Very High	> 250	> 355

**24 Hour Urine**

Reference Value	Men and Women (mg/24 hours)
Normal	< 30
High	30 - 300
Very High	> 300

**Urine 4, 6 or 12 Hours**

Urines to 4, 6 or 12 hours perform the same calculations used for 24-hour urine, but with the total volume of urine collected during the time interval.

Collection Time	Reference Value Men e Women	
Urine of 4 Hours	< 5 (mg/4 hours)	Or < 20 µg/min
Urine of 6 Hours	< 7,5 (mg/6 hours)	
Urine of 12 Hours	< 15 (mg/12 hours)	

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

**PRODUCT PERFORMANCE****QUALITY CONTROL****Accuracy****RECOVERY**

The recovery analysis was performed with 05 determinations samples. Accuracies were calculated and met in good agreement with the reference values, obtaining a recovery between 91 and 102%.

**COMPARISON OF METHODS AND SPECIFICITY METHODOLOGY**

Microalbumin kit was compared with another method for measurement of Microalbumin commercially available. 07 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 0,995X + 0,2262$  and coefficient correlation 0,9999. With these results we can conclude that the kit shows good methodological specificity.

**Precision****REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	9,26	20,45	58,50
<b>Standard Deviation (mg/L)</b>	0,01	0,02	0,22
<b>Coefficient of Variation (%)</b>	0,14	0,08	0,38

**REPRODUCIBILITY**

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	9,25	20,45	58,53
<b>Standard Deviation (mg/L)</b>	0,01	0,01	0,03
<b>Coefficient of Variation (%)</b>	0,06	0,03	0,05

**Sensitivity**

The sensitivity was calculated from 20 determinations of a sample free of Microalbumin. The average found was 0,0258 mg/L, with standard deviation of 0,0002 mg/L. The sensitivity, which indicates the method detection limit, corresponds to the average plus 3 times the standard deviation and is equal to 0,0264 mg/L.

**Linearity**

The linearity depends on the concentration value of the calibrator used and can not be greater than 400 mg/L. For higher values, dilute sample 1:5 with distilled or deionized water and repeat the determination. Multiply the result by the factor dilution.

**DIAGNOSTIC SIGNIFICANCE**

Diabetic nephropathy, which is accompanied by irreversible damage kidneys and persistent proteinuria is the leading cause of death people with insulin dependent *diabetes mellitus*. A signal early diabetic nephropathy is the presence of Albumin in urine at low concentrations. Early detection of minor kidney damage is very important, because the rapid onset of treatment allows reversion of the condition.

**BIBLIOGRAPHIC REFERENCES**

1. Mount, JN, J. Clin. Pathology, 22, 12 (1986).
2. Schmidt, A., et al., Diabetic Medicine, 5, 126 (1988).
3. Killingsworth, LM; SAVORY, J: Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem. 1971, 17, 936.
4. PESCE, AJ, KAPLAN, LA: Methods in Clinical Chemistry CV Mosby Company, 1987.
5. Hellsing, K., Profiles in the Biological Fluids, 1973, 23: 579.
6. BLOM, M. and HOJME, H., Clinical Chemistry, 1976, 22: 657.
7. HILLS, LP and TIFFANY, IT, Clinical Chemistry, 1980, 26: 1459.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: [bioclin@bioclin.com.br](mailto:bioclin@bioclin.com.br)  
CNPJ: 19.400.787/0001-07 - Made in Brazil

**EC REP OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium


**CUSTOMER SERVICE**

Customer Advisory Service  
Phone: 0800 0315454  
E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

ANVISA registration for Microalbumin kit: 10269360206

Review: October/2013

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED