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### USAGE INSTRUCTIONS

#### FUNCTION

Method of quantitative determination of Rheumatoid Factor (RF). Immunoturbidimetric test, for *in vitro* diagnostic only.

#### PRINCIPLE OF ACTION

**Methodology:** immunoturbidimetry.

The reaction permits you to quantify, using the turbidimetric method, the concentration of RF present in sample. The suspended latex particles covered with human Gamma globulin are agglutinated by Rheumatoid Factors in serum. The agglutination process causes an increase in optical density. The intensity of scattered light is proportional to the concentration of RF.

#### REAGENTS

**Number 1 - Buffer** - store in temperatures between 2 and 8°C. Contains: Sodium Chloride 0,15 mol/L, Tris 50 mmol/L, sodium azide 15,38 mmol/L, surfactant.

**Number 2 - RF Latex** - store in temperatures between 2 and 8°C. Suspension of latex particles sensitized with human gamma globulin. Contains: Sodium azide 15.38 mmol/L.

**Number 3 - Calibrator** - store in temperatures between 2 and 8°C. Contains: Rheumatoid Factor and Sodium azide 15,38mmol/L.

**Warning: The concentration of Rheumatic Factor varies according to the lot - see the bottle label.**

#### PRESENTATION

Reagent	Presentations		
	K058-1	K058-2	K058-3
Reagent Nº 1	18 mL	45 mL	2 x 45 mL
Reagent Nº 2	2 mL	5 mL	2 x 5 mL
Reagent Nº 3	1,0 mL	1,0 mL	2 x 1,0 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Multical and Multicontrol Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. **Do not freeze.** Protect from light and avoid moisture.

#### SPECIAL CARE

1 - For *in vitro* diagnostic only.

2 - Strictly follow the methodology proposed to obtain exact results.

3 - Water used in material cleaning must to be recent and free of contaminants.

4 - Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Reagents 1, 2 and 3 should be handled cautiously since they are subject to biological contamination.

6 - Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.

7 - The reagent has been tested for HIV and HBsAg using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially Infectious.**

8 - We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9 - To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Advisory Service Customer) of Quibasa.

#### SAMPLES

Serum obtained free of hemolysis and severe lipemia. The analyte is stable between 2 and 8 up to 2 days and 3 months at -20 °C.

#### PROCESS DESCRIPTION

##### TECHNIQUE

Bioclin recommends, as calibrator Multical Bioclin Kit and as control Multicontrol Bioclin Kit.

This kit should only be used in automated equipment. **Check programming through the website [www.bioclin.com.br](http://www.bioclin.com.br) or through the SAC: 0800 031 5454.**

##### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

##### TRACEABILITY

The kit's calibrator is traceable to the reference material NIBSC 64/002 (Rheumatoid Arthritis Serum).

#### REFERENCE VALUES

The reference values for this method were obtained through determination of Rheumatoid Factor in healthy male and females populations. Less than 20 UI/mL. These values should be used as guidance and each laboratory should establish its range of reference values, according to the population served.

#### PRODUCT PERFORMANCE

##### QUALITY CONTROL

##### Accuracy

##### RECOVERY

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 95 and 107%.

#### COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The Bioclin Rheumatoid Factor kit was compared with other commercially available methods for determination of Rheumatic Factor dosage. 07 analyses were performed and the results were evaluated. We obtained the linear equation  $Y = 1,0163 + 0,9987X$  and correlation coefficient = 0,998677. With these results we can conclude that the kit has good methodological specificity.

##### Precision

##### REPEATABILITY

20 successive measurements were performed with three samples with different concentration, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (UI/mL)	24,60	150,5	84,80
Standard Deviation (UI/mL)	0,995	2,704	1,673
Coefficient of Variation (%)	4,044	1,797	1,973

	Sample 1	Sample 2	Sample 3
Concentration (UI/mL)	24,58	153,28	84,33
Standard Deviation (UI/mL)	0,08	2,47	0,43
Coefficient of Variation (%)	0,31	1,61	0,51

### Sensitivity

Sensitivity was calculated from 20 determinations of a sample with concentration of 0 (zero) Rheumatic Factor. The average 7,26 UI/mL with a standard deviation of 0,89 UI/mL. The sensitivity, which indicates the Method Detection Limit, corresponds to the average plus 3 times the value of the standard deviation and is equal to 9,919 UI/mL.

### Linearity

The reaction is linear until the highest point in the calibrator curve. For samples with higher values than it is recommended to dilute the sample with sodium chloride 0.85% repeat the dosage and multiply the result by the dilution factor.

### DIAGNOSTIC SIGNIFICANCE

Rheumatoid factor refers to a group of macroglobulins (antiglobulin) that react with the Fc fragment of immunoglobulins IgG. In rheumatoid arthritis the finding of Rheumatoid Factor represents the most important serological data, present in about 75% of patients. Its plasma level can be significantly increased in old age, connective tissue diseases, liver diseases syphilis, tuberculosis, leprosy, bacterial endocarditis, mononucleosis, sarcoidosis, leishmaniasis, measles, cancer, parasitic infestations, blood transfusions, renal transplant, Sjogren syndrome.

### BIBLIOGRAPHIC REFERENCES

- 1 - PESCE, A. J.; KAPLAN, L. A.: Methods in Clinical Chemistry the C. V. Mosby Company, 1987.
- 2 - HELLSING, K., Profides in the Biological Fluids, 1973, 23, 579.
- 3 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976, 22,657.
- 4 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980;26, 1459.
- 5 - GALVIN, J. P. et al; Clin. Lab Assays 73, 4th, 1983.
- 6 - PASSING, H.; BABLOCK, W., J., Clin. Chem. Clin. Biochem 21, 709; 1983.
- 7 - NOKAMA, Y.; NAKAMURA, R. M., J. Clin. Lab. Anal., 1:15; 1987.
- 8 - OKOMURA, J. M.; MIYAGI, J. M.; TERADA, K. J., Clin. Lab. Anal. 4, 231; 1990.
- 9 - MOORE, T. L., Clin. Biochem, 26; 75-84; 1993.
- 10 - SINGER, J. M. Et al Amer J Med 21, 888; 1952.

### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

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