BIOLÁTEX PCR

REF K044

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of C Reactive Protein (CRP), by agglutination of latex particles, without previous dilution of the sample. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Latex

This method is based on a particle agglutination reaction of latex particles covered with Gamma-Globulin anti-CRP, especially treated to prevent nonspecific agglutination. Agglutination is visible in samples with concentrations of CRP equal to or greater than 6 mg/L, according to the references established by the WHO International Standards.

REAGENTS

Number 1 - CRP Latex - Store between 2 and 8°C. **Do not freeze.** Contains: Sensitized latex particles in suspension. Homogenize before using it.

Number 2 - Positive Control - Store between 2 and 8°C. Contains: Serum with concentration equal or greater than 6 mg/L, Sodium Azide 15,38 mmol/L. **Potentially infectious.**

Number 3 - Negative Control - Store between 2 and 8°C. Contains: Serum with concentration less than 6 mg/L, Sodium Azide 15,38 mmol/L.

PRESENTATION

| Reagent | Presentation | |
|---------|----------------------|--------------------|
| | 1 (Economic Packing) | 2 (Normal Packing) |
| Nº 1 | 2 mL | 2 mL |
| Nº 2 | -- | 1 mL |
| Nº 3 | -- | 1 mL |

The normal package is accompanied by spatulas and test cards.

EQUIPMENTS AND OPERATIONAL INPUTS

Slides or dark bottom plates, spatulas, pipettes and watch or stopwatch. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Do not freeze reagents.

5- Always use reagents from same lot.

6- Do not use lipemic serum. Do not use plasma.

7- Handle with care Reagents Nº 2 and 3, they contain Sodium Azide.

8- Negative and positive controls are composed human fluids and were analyzed to detect antibodies Anti HIV and Anti HBs, with results turning negative. However, for greater safety, consider handling as **potentially infectious.**

9- We recommend applying the local, state and federal rules

for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

10- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

11- Do not use the product in case of damaged packaging.

12- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum, with no previous dilution. The analyte is stable for 02 days if kept between 2 and 8°C.

PROCESS DESCRIPTION

TECHNIQUE

In each circle of the slide, place the following:

| | Circle Nº 1 | Circle Nº 2 | Circle Nº 3 |
|-------------------------|-------------|-------------|-------------|
| Negative Control | 20 µL | -- | -- |
| Positive Control | -- | 20 µL | -- |
| Sample | -- | -- | 20 µL |
| Reagent Nº 1* | 20 µL | 20 µL | 20 µL |

*Previously homogenized

Homogenize with the aid of a spatula using the entire length of each circle of the slide. Then, agitate the slide with circular motions for two minutes. Make the reading with an artificial light, using a dark background to facilitate test interpretation.

A clear agglutination indicates the presence of C Reactive Protein at a concentration equal to or greater than 6 mg/L. In this case, perform the semi-quantitative test.

SEMIQUANTITATIVE TEST

1- Perform sample dilutions with saline, starting from initial sample (1:2, 1:4, 1:8, 1:16, 1:32, etc.);

2- Follow the process as described in qualitative test for each of the dilutions.

It will be considered as title the biggest serum dilution that presents agglutination.

RESULTS

Positive: Clear presence of agglutination.

Negative: Absence of agglutination (homogenized suspension).

CALCULATIONS

| Sample | Concentration (mg/L) |
|-------------|----------------------|
| No dilution | 6 |
| 1/2 | 12 |
| 1/4 | 24 |
| 1/8 | 48 |
| 1/16 | 96 |
| 1/32 | 192 |

Results can be expressed in titles or in mg/L.

mg/L = 6 x title from last dilution (dilution nº)

Negative Tests: express the results as negative or minor than 6 mg/L.

PROCEDURE LIMITATIONS

Do not use plasma, hemolysate serum or limpemic, as they may produce nonspecific agglutination.

By correlating methods for determining C Reactive Protein, to verify the sensitivity of the reagents. The obtained results should only be compared when expressed in mg/L.

INTERNAL QUALITY CONTROL

The kit's sensitivity was determined through the reference material NIBSC 85/506 (Human C-Reactive Protein 1st International Standard - WHO International Standard).

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's sensitivity was determined through the reference material NIBSC 85/506 (Human C-Reactive Protein 1st International Standard - WHO International Standard).

REFERENCE VALUES

Up to 6 mg/L

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Biolatex CRP kit was compared with another method for confirmation of inflammatory process. According to results from 10 clinical samples, all methods showed excellent correlation. With these results we can conclude the kit has good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 5 different samples, obtaining the following results:

| Samples | Nº of Repetitions | Expected Results | Obtained Results |
|---------|-------------------|------------------|------------------|
| 01 | 20 | Positive 1/8 | Positive 1/8 |
| 02 | 20 | Negative | Negative |
| 03 | 20 | Negative | Negative |
| 04 | 20 | Positive 1/2 | Positive 1/2 |
| 05 | 20 | Positive 1/32 | Positive 1/32 |

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

| Samples | N° of Repetitions | Expected Results | Obtained Results |
|---------|-------------------|------------------|------------------|
| 01 | 20 | Negative | Negative |
| 02 | 20 | Positive 1/4 | Positive 1/4 |
| 03 | 20 | Negative | Negative |
| 04 | 20 | Positive 1/32 | Positive 1/32 |
| 05 | 20 | Positive 1/16 | Positive 1/16 |

Sensitivity

Clinical or diagnostic sensitivity and methodological sensitivity were determined through samples of 20 human serum truly positive. Tests show 100% positive results, no false negative results were found.

DIAGNOSTIC SIGNIFICANCE

C Reactive Protein is a useful indicator of inflammatory process in activity, whether it is infectious (pneumonia, tuberculosis) or non-infectious (active rheumatic fever, rheumatoid arthritis, lupus erythematosus).

Is also present in several other pathological conditions such as in acute myocardial infarction, neoplasms, intense trauma, viruses and burns.

The determination of plasma concentrations constitutes an effective test to monitor the therapy used and prognosis of inflammation.

NUMBER OF TESTS

Presentation 1 (Economic Packing) 100 Tests / 20 µL Sample / 20 µL Reagent

Presentation 2 (Normal Packing) 100 Tests / 20 µL Sample / 20 µL Reagent

BIBLIOGRAPHIC REFERENCES

- 1 - WARWORTH, E. Wadsworth, Ch. Clin. Chim. Acta, 138, 1984.
- 2 - PEPYS, M. B.; DASH, A. C.; ASHLEY M. J., Clin. Exp. Immunol, 30, 32-37, 1977.
- 3 - DEYO, R. A.; POPE, R. M., PERSELLIN, R. H.; J. Rheumatol, 279, 1980.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Biolatex CRP kit: 10269360093

Review: September/2014

UNIVERSAL SYMBOLOGY

| | | | |
|---|---------------------------------------|---|--|
|  | CATALOG NUMBER |  | MANUFACTURED BY |
|  | BATCH CODE |  | CONTROL |
|  | DATE OF MANUFACTURE |  | POSITIVE CONTROL |
|  | USED BY (last day of month) |  | NEGATIVE CONTROL |
|  | TEMPERATURE LIMITATION (store at) |  | BIOLOGICAL RISK |
|  | CONTAINS SUFFICIENT FOR <N> TESTS |  | INFLAMMABLE |
|  | CONSULT INSTRUCTIONS FOR USE |  | CORROSIVE |
|  | IN VITRO DIAGNOSTIC DEVICE |  | POISON |
|  | EUROPEAN AUTHORIZED REPRESENTATIVE |  | CE MARK |
|  | KEEP AWAY FROM SUNLIGHT |  | DO NOT USE IF PACKAGE IS DAMAGED |