## TRANSFERRINA

REF K114

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## FINALIDADE

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for determination of transferrin in samples of human serum. Immunoturbidimetric test, for *in vitro* use only.

#### PRINCIPLES OF ACTION

**Methodology:** Immunoturbidimetry

Turbidimetric method in which the specific antibodies present in reagent form insoluble complexes with serum transferrin, generating turbidity, with intensity proportional to the quantity of transferrin in sample.

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: PEG 6.000 < 15%, PBS buffer and preservative.

**Number 2 - Antiserum** - Store between 2 and 8°C. Contains: Goat antiserum for Transferrin < 3,8 g/L, PBS buffer and preservative

#### PRESENTATION

Presentation	Reagent	
	Nº 1	Nº 2
1	25 mL	2 mL
2	2 x 25 mL	2 x 2 mL
7	20 mL	5 mL
8	2 x 20 mL	2 x 5 mL
9	3 x 20 mL	3 x 5 mL
10	4 x 20 mL	4 x 5 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the biochemistry automated equipment, Multical kit and Multicontrol Bioclin.

#### TRANSPORTATION AND STORAGE CONDITIONS

Store temperatures must be between 2 and 8°C.

Transportation in temperatures between 15 and 30°C must not exceed 72 (seventy two) hours. Keep away from light and avoid humidity. **Do not freeze.**

#### SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Water used in material cleansing must be recent and free of contaminants.
- 3- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 4- Reagents Nº 1 and Nº 2 contain Sodium Azide, irritant to skin and mucous. Handle with care.
- 5- Do not freeze the reagents, the freezing of antibody or buffer can affect their performance.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Fresh serum or plasma obtained with heparin or EDTA. The analyte is stable for 2 days at 2 to 8 ° C. If you need to store for a longer period, the sample should be frozen.

#### PROCESS DESCRIPTION

#### REAGENTS PREPARATION

Reagents are ready to be used.

#### TECHNIQUE

This kit is indicated to use in automated biochemical analyzers only. **Verify equipment set up through the website [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC.**

#### INTERFERENCE

No interferences to hemoglobin up to 1000 mg/dL. Sodium citrate up to 1000 mg/dL, heparin up to 50 mg/dL, bilirubin up to 20 mg/dL and triglycerides up to 2500 mg/dL.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### REFERENCE VALUES

Reference values were obtained through the determination of transferrin in healthy populations of male and female.

Reference Values: 170 - 340 mg/dL (IFCC)

These values must be used only as a reference. Each laboratory must create his own range of reference, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

#### PRODUCT PERFORMANCE

##### QUALITY CONTROL

#### Accuracy

##### COMPARISON OF METHODS

The kit Transferrin was compared with other methods commercially available for measurement of transferrin. For this test, 41 clinical samples were analyzed. The linear equation  $Y = 0,976X + 3,228$  with correlation coefficient equal to 0,992. With these results, we can conclude that the kit has good methodological specificity.

#### Precision

##### REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

Repeatability	Samples	
	1	2
Average Concentration (mg/dL)	217,57	129,33
Standard Deviation (mg/dL)	2,03	0,93
Coefficient of variation (%)	0,93	0,72

##### REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

Reproducibility	Samples	
	1	2
Average Concentration (mg/dL)	217,82	129,33
Standard Deviation (mg/dL)	0,34	1,13
Coefficient of variation (%)	0,16	0,87

**Sensitivity**

The sensitivity was calculated from 10 determinations of a sample free of Transferrin. The average found was 1,50 mg/dL, with standard deviation of 0,08 mg/dL. The sensitivity, which indicates the method detection limit, corresponds to the average plus 3 times the standard deviation and is equal to 1,74 mg/dL.

**Linearity**

The reaction is linear up to the level of the highest calibrator concentration.

**DIAGNOSTIC SIGNIFICANCE**

Transferrin is a glycoprotein synthesized in the liver and is the main plasma protein for iron transport. Its concentration is correlated with the total capacity of serum iron binding. Dosage levels of transferrin in plasma are useful for differential diagnosis of anemia and to monitor their treatment. Increased levels are found in iron deficiency in pregnancy, administration of estrogen and hereditary hemochromatosis. Decreased levels can be found in hereditary deficiencies, testosterone administration, infection, acute inflammation, some forms of necrosis, tumors, acute malaria and in cases of malnutrition.

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**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

**OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium


**CUSTOMER SERVICE**

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for Transferrin kit: 10269360186

**Review:** June/2015

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED