# ÁCIDO ÚRICO MONOREAGENTE

REF K139

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### FINALIDADE

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### USAGE INSTRUCTIONS


#### FUNCTION

Method for determination of Uric Acid. Enzymatic colorimetric test, only for *in vitro* diagnostic use.

#### PRINCIPLE OF ACTION

**Methodology:** Enzymatic colorimetric (UOD-PAP)

The enzymatic determination of Uric Acid is made according to the following reactions:

The intensity of the cherry color formed is directly proportional to the concentration of Uric Acid in the sample.

#### REAGENTS

**Number 1 - Enzyme Reagent** - Store between 2 and 8°C. Do not freeze. Contains: Buffer  $\geq 100$  mmol/L, 4-Aminoantipyrine  $< 1$  mmol/L, Peroxidase  $< 18.000$  U/L, Uricase  $< 3.000$  U/L, DHBS  $< 5$  mmol/L, Stabilizer, Surfactant and Preservative.  
**Number 2 - Standard** - Store between 2 and 8°C. Contains: Uric Acid 6,0 mg/dL, Buffer  $< 500$  mmol/L, Stabilizer and Preservative.

#### PRESENTATION

Apresentation	Reagent N° 1	Reagent N° 2
1	100 mL	3 mL
2	2 x 100 mL	3 mL
3	4 x 100 mL	3 mL
4	5 x 20 mL	3 mL
5	2 x 40 mL	3 mL
6	3 x 40 mL	3 mL
7	4 x 40 mL	3 mL
8	5 x 40 mL	3 mL
9	6 x 40 mL	3 mL
10	2 x 60 mL	3 mL
11	3 x 60 mL	3 mL
12	4 x 60 mL	3 mL
13	8 x 60 mL	3 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze**

#### SPECIAL CARE

**1- For professional *in vitro* diagnostic use only.**

- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Water level in water bath should be above the level of reagents in test tubes.
- Handle with care Reagent N° 1, containing Sodium Azide, irritating to skin and mucous.
- Slightly pink color development in the Reagent N° 1 does not damage their quality, as long as the corresponding white and standard dosages are used.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance Blank current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum obtained free of hemolysis, Plasma (collected with heparin or EDTA) and Urine. The analyte is stable in plasma, serum or urine for 3 days if kept between 2 and 8°C, and 7 days if kept at 10°C negative.

#### PROCESS DESCRIPTION

#### TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark three test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Sample	---	---	25 $\mu$ L
Reagent N° 2	---	25 $\mu$ L	---
Reagent N° 1	1,0 mL	1,0 mL	1,0 mL

Homogenize well and put in a water-bath at 37°C for 5 minutes. Read the absorbance of Sample and Standard at 505 nm (490 - 540 nm), hitting the zero with the Blank. The color is stable for 30 minutes.

**24 - Hour urine** - Preparation of sample:

Take note of the volume (V) of 24-hour urine in mL, homogenize, separate a rate of 10 mL and adjust the pH between 7 and 9 with 5% NaOH. If the sample is cloudy, warm it in water bath for 10 minutes at 56°C, to dissolve the crystals. Dilute urine 1:10 (100  $\mu$ L of urine + 900  $\mu$ L of distilled or deionized water). Multiply the result by 10. Proceed the dosage as described for serum.

#### CALCULATIONS

$$\text{Uric Acid (mg/dL)} = \frac{\text{Sample Absorbance} \times 6.0}{\text{Standard Absorbance}}$$

As the reaction follows the Beer-Lambert Law, the Calibration Factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (6.0 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{mg/dL} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

#### Urine

$$\text{Uric Acid (mg/24h)} = \frac{\text{mg/dL} \times \text{Urinary Volume (mL)}}{100}$$

Results are expressed in mg/dL.

#### PROCEDURE LIMITATIONS

Use plasma samples containing only EDTA or Heparin. Other anticoagulants interfere with the reaction, with Fluoride Uricase inhibitor.

For serum lipids or Bilirubin levels ( $> 10$  mg/dL), mix 2,0 mL of 0,85% NaCl and 0,05 mL of serum. Measure absorbance at 505 nm, hitting the zero with distilled water. Decrease the absorbance thus obtained in the absorbance of the sample, and calculate the concentration in mg/dL.

The use of highly reducing drugs, such as Ascorbic Acid (Vitamin C), interferes with the reaction, because they compete with the consumption of  $\text{H}_2\text{O}_2$  by providing falsely decreased values. For this reason, you should discontinue its use at least 12 hours before collection of sample.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 913A.

#### REFERENCE VALUES

The reference values, in mg/dL, for this method were obtained through the determination of Uric Acid in healthy populations of male and female.

Serum	Adults	Male	2,5 to 7,0 mg/dL
		Female	1,5 to 6,0 mg/dL
	Children	Male	1,5 to 6,0 mg/dL
		Female	0,5 to 5,0 mg/dL
Urine	250 to 750 mg/24 hours		

To convert the values from mg/dL to mmol/L (IS), multiply by 0,05948.

These values should be used as guidelines, each laboratory should create their own range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

## PRODUCT PERFORMANCE

### QUALITY CONTROL

#### Accuracy

#### COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The Uric Acid Monoreagent kit was compared with other method for determination of Uric Acid dosage commercially available. 42 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 1,006X + 0,035$  and correlation coefficient 0,999. With these results we can conclude that the kit shows good methodological specificity.

#### Precision

##### REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	4,03	8,11	6,34
<b>Standard Deviation (mg/dL)</b>	0,02	0,05	0,04
<b>Coefficient of Variation (%)</b>	0,51	0,58	0,55

##### REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	4,03	8,07	6,16
<b>Standard Deviation (mg/dL)</b>	0,01	0,07	0,15
<b>Coefficient of Variation (%)</b>	0,16	0,87	2,42

#### Sensitivity

Sensitivity was calculated from 40 determinations of a sample free of Uric Acid. The average found was 0,031 mg/dL with standard deviation of 0,004 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 0,043 mg/dL.

#### Linearity

The reaction is linear up to 20,0 mg/dL. For samples with values above 20 mg/dL or optical density above 0,8, it is recommended diluting the sample with distilled water, repeat the dosage and multiply the result by the dilution factor.

#### DIAGNOSTIC SIGNIFICANCE

Uric Acid is the end product of metabolism of Nucleic Acids and Purines. Its concentration in body fluids depends on the balance between production and delimitation through the kidneys, and generated about 400 mg daily.

The biochemical characteristics and the basic condition for the diagnosis of gout is hyperuricemia. From the epidemiological point of view, Uric Acid levels above 7 mg/dL may indicate high risk of gouty arthritis or nephrolithiasis.

The concentration of plasma Uric Acid is elevated in several other clinical conditions such as renal failure, congestive heart failure, toxemia of pregnancy (eclampsia), leukemias and lymphomas, ketoacidosis, Down's syndrome. Often the observation of hyperuricemia associated with factors such as obesity, hypertriglyceridemia, hypertension and mellitus diabetes.

Reduced levels of Uric Acid are observed in Fanconi syndrome, Wilson's disease and in patients treated with Allopurinol, Corticosteroids, ACTH and Phenylbutazone.

#### BIBLIOGRAPHIC REFERENCES

- 1 - TRIVEDI, R. C.; REBAR, L.; BERKA, E.; STRONG, L., Clin. Chem., 1978, 24, 1908.
- 2 - TONKS, D. B., Quality Control in Clinical Laboratories, 1983.
- 3 - TRINDER, P., Ann. Clin. Biochem., 1969, 6, 24.
- 4 - TIETZ Textbook of Clinical Chemistry, 2nd ed., 1994.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

#### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 ( 31 ) 3439.5454 - Fax: +55 ( 31 ) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

#### OBELIS S.A.

Bd. Général Wahis, 53  
1030 Brussels, Belgium

#### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for Uric Acid Monoreagent kit: 10269360228

**Review:** August/2014

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY  
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION  
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT  
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS  
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED  
REPRESENTATIVE

CE MARK

KEEP AWAY  
FROM SUNLIGHT

DO NOT USE IF  
PACKAGE IS  
DAMAGED