HDL DIRETO

REF K071

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FINALIDADE

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- 1 - WILLIAM, P., ROBINSON, D., BAILY, A.. High Density Lipoprotein and Coronary Risk Factor. Lancet, 1:72 (1979).
- 2 - BADIMON, J.J., BADIMON, L., FUESTER, V.. Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma. Fraction in the Cholesterol-Fed Rabbit. Journal of Clinical Investigation, 1990;85:1234-41.
- 3 - BARR, D.P., RUSS, E.M., EDER, H.A.. Protein-Lipid Relationships in Human Plasma. Am. J. Med., 11:480, 1951.
- 4 - GOTTO, A.M.. Lipoprotein Metabolism and Etiology of Hyperlipidemia. Hospital Practice, 23: Suppl. 1,4 1988.
- 5 - Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of HDL Cholesterol. Colorimetric Enzymatic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic - Colorimetric

The HDL Direct kit allows a direct determination without need for prior precipitation or sample treatment. Reagents are used to carry out selective dosage of cholesterol bound to HDL. The surfaces of low density lipoproteins (LDL), very low density (VLDL) and chylomicrons are stabilized by adsorption of polyanion and do not undergo the action of modified enzymes present in reagent. HDL, however, solubilize under the action of a detergent, allowing the enzymatic action on cholesterol attached to it. Thus, the intensity of color formed is proportional to the concentration of HDL cholesterol in the sample.

REAGENTS

Number 1 - Polyanion - Store between 2 and 8°C. Contains: Polyanion 0,2% and 4-Aminophenazone 2,0 mmol/L.

Number 2 - Enzymatic Reagent - Store between 2 and 8°C. Contains: Buffer (pH 7,0) 100 mmol/L, DSBT 24 mmol/L, Cholesterol Esterase < 300 U/L, Cholesterol Oxidase < 300 U/L, Peroxidase < 1000 U/L and detergent.

Number 3 - Calibrator - Store between 2 and 8°C. Contains: Lyophilized Serum with HDL value set. **Potentially infectious.**

Warning: The concentration varies with the lot - see bottle label.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2	Reagent N° 3
1	60 mL	20 mL	2 mL
2	240 mL	80 mL	2 mL
3	2 x 30 mL	2 x 10 mL	2 mL
4	30 mL	10 mL	2 mL
5	4 x 30 mL	4 x 10 mL	2 mL
6	6 x 30 mL	6 x 10 mL	2 mL
7	2 x 60 mL	2 x 20 mL	2 mL
8	4 x 60 mL	4 x 20 mL	2 mL
9	3 x 60 mL	3 x 20 mL	2 mL
23	60 mL	20 mL	5 mL
24	240 mL	80 mL	5 mL
25	30 mL	10 mL	5 mL
26	2 x 30 mL	2 x 10 mL	5 mL
27	4 x 30 mL	4 x 10 mL	5 mL
28	6 x 30 mL	6 x 10 mL	5 mL
29	2 x 60 mL	2 x 20 mL	5 mL
30	3 x 60 mL	3 x 20 mL	5 mL
31	4 x 60 mL	4 x 20 mL	5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, watch or stopwatches, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

6- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

7- Do not use the product in case of damaged packaging.

8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, plasma obtained with heparin. The cholesterol HDL is stable in serum for 07 days between 2 and 8°C.

PROCESS DESCRIPTION

REAGENTS PREPARATION

Before reconstitute the Calibrator should be noted the volume indicated on the bottle and proceed as described below.

Presentation 1 to 9: Open carefully the bottle and add exactly **2 mL of distilled water**, as volume described in the label of the bottle.

Presentation 23 to 31: Open carefully the bottle and add exactly **5 mL of distilled water**, as volume described in the label of the bottle.

Mix. Leave the bottle resting for 30 minutes. Every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid formation of foam.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes: B (Blank), C (Calibrator), S (Sample) and proceed as follows:

	Blank	Calibrator	Sample
Reagent N° 1	300 µL	300 µL	300 µL
Calibrator	---	3 µL	---
Sample	---	---	3 µL

Homogenize and place it in a water bath at 37°C for 5 minutes. Read the absorbance from Sample (S1) and Calibrator (C1) in 550 nm (540 - 560 nm), adjusting the zero with the Blank.

Reagent N° 2	100 µL	100 µL	100 µL

Homogenize and place it in a water bath at 37°C for 5 minutes. Read the absorbance from Sample (S2) and Calibrator (C2) in 550 nm (540 - 560 nm), adjusting the zero with the Blank. Color is stable for 60 minutes.

CALCULATIONS

HDL Cholesterol = $\frac{(A2 - A1)}{(C2 - C1)} \times \text{Calibrator Concentration}$

Results are expressed in mg/dL.

PROCEDURE LIMITATIONS

Do not use anticoagulants containing citrate. Do not expose the reactive to light.

INTERFERENCES

The Ascorbic Acid up to 50 mg/dL, Hemoglobin up to 500 mg/dL, Bilirubin up to 30 mg/dL and Triglycerides up to 1000 mg/dL did not alter the results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's calibrator is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 1951.

REFERENCE VALUES

The reference values in mg/dL, for this method, were obtained through the determination of HDL cholesterol in healthy populations of male and female.

Adults (≥ 20 years)	Desirable	> 60 mg/dL
	Limit	40 - 60 mg/dL
	Low	< 40 mg/dL

Childrens and Adolescents (2 to 19 years)	Desirable	> 45 mg/dL

To convert the data in International System (IS), multiply the results: $\text{mg/dL} \times 0,02586 = \text{mmol/L}$

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

RECOVERY

The recovery analysis was performed with 05 determinations samples. Accuracies were calculated and found themselves in good agreement with the reference values, obtaining a recovery between 95% and 103%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The HDL Direct kit was compared with other method for HDL Cholesterol dosage commercially available. 42 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 0,989X + 0,361$ and the correlation coefficient 0,998. With these results we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average concentration (mg/dL)	16,98	26,18	21,05
Standard Deviation (mg/dL)	0,16	0,38	0,22
Coefficient of Variation (%)	0,9	1,5	1,0

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average concentration (mg/dL)	16,70	26,16	21,13
Standard Deviation (mg/dL)	0,52	0,03	0,67
Coefficient of Variation (%)	3,12	0,11	3,15

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of HDL Cholesterol. The average found was 0,66 mg/dL with a standard deviation of 0,08 mg/dL. The sensitivity, which indicates the method detection limit, correspond the average plus 3 times the standard deviation, and is equal to 0,91 mg/dL.

Linearity

Reaction is linear from 3 to 200 mg/dL. For values higher than 200 mg/dL, repeat the determination diluting the sample with NaCl 0,85%. Multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The plasmatic lipoproteins are spherical particles consisting of Cholesterol, Triglycerides, and Phospholipids and Proteins. These particles are responsible for transportation of the cholesterol from peripheral cells and liver. The relative proportion between proteins and lipids determine the density of lipoproteins - chylomicrons, VLDL, LDL and HDL. Monitoring HDL Cholesterol in serum is of great importance due to inverse correlation between plasmatic concentrations and the risk coronary disease.

BIBLIOGRAPHIC REFERENCES

- 1 - WILLIAM, P., ROBINSON, D., BAILY, A.. High Density Lipoprotein and Coronary Risk Factor. Lancet, 1:72 (1979).
- 2 - BADIMON, J.J., BADIMON, L., FUESTER, V.. Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma. Fraction in the Cholesterol-Fed Rabbit. Journal of Clinical Investigation, 1990; 85: 1234-41.
- 3 - BARR, D.P., RUSS, E.M., EDER, H.A.. Protein-Lipid Relationships in Human Plasma. Am. J. Med., 11:480, 1951.
- 4 - GOTTO, A.M.. Lipoprotein Metabolism and Etiology of Hyperlipidemia. Hospital Practice, 23: Suppl. 1,4 1988.
- 5 - Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for HDL Direct kit: 10269360122

Review: November/2014

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED