Bioclin

FERRO SÉRICO

REF K017

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Iron Serum. Colorimetric test, only for *in vitro* diagnostic.

PRINCIPLE OF ACTION

Methodology: Modified Goodwin

Iron is released from Transferrin in acidic and low to its ferrous state by the action of Hydroxylamine. Subsequently, reacts with Ferrozine leading to the formation of a violet color complex.

REAGENTS

Number 1 - Reduction Buffer - Store between 15 and 30°C. Contains: Succinic Acid 250 mmol/L, Hydroxylamine 140 mmol/L, surfactant and stabilizer.

Number 2 - Color Reagent - Store between 15 and 30°C. Contains: Ferrozine 12 mmol/L in Glacial Acetic Acid solution.

Number 3 - Standard - Store between 15 and 30°C. Contains: Iron 100 µg/dL.

PRESENTATION

Reagent	Volumen
Nº 1	40 mL
Nº 2	2,5 mL
Nº 3	5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watch and stopwatch, pipettes test tubes, Biocontrol N and Biocontrol P Bioclin. The instruments and glassware can be found in the specialized market of clinical analysis laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature shall be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Water level in water bath should be above the level of reagents in test tubes.
- 6- Hemolysis, even if discrete, can interfere with dosage.
- 7- Do not use plasma.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Sample obtained free of hemolysis, stable for 7 days between 2 and 8°C. Therapeutic control samples should always be collected at the same time.

PROCESS DESCRIPTION

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes with letters B (Blank), P (Standard), A (Sample) and proceed as follows:

	Blank	Standard	Sample
Reagent Nº 1	2,0 mL	2,0 mL	2,0 mL
Distilled Water	500 µL	---	---
Sample	---	---	500 µL
Reagent Nº 3	---	500 µL	---

Smoothly homogenize and determine the absorbance of the sample at 560 nm (540 - 580 nm), hitting the zero with the Blank. This will be the absorbance A_1 .

And then add:

	Blank	Standard	Sample
Reagent Nº 2	1 drop	1 drop	1 drop

Smoothly homogenize and incubate in water-bath at 37°C, per 10 minutes. Determinate the absorbancies of Sample and the Standard in 560 nm (540 - 580 nm), matching the zero with the Blank. This will be the absorbance of the sample A_2 .

CALCULATIONS

Iron Serum (µg/dL) = $\frac{(A_2 - A_1)}{\text{Standard Abs.}} \times 100$

The color reaction follows the Lambert Beer Law, it can use the Calibration Factor:

$$\text{Calibration Factor} = \frac{100}{\text{Standard Absorbance}}$$

$$\mu\text{g/dL} = (A_2 - A_1) \times \text{Calibration Factor}$$

The results are expressed in µg/dL.

PROCEDURE LIMITATIONS

The use of water contaminated by Iron increase Blank color intensity, changing the results. Some drugs can raise levels of Serum Iron such as: Chloramphenicol, Estrogens, Ethanol, Methotrexate and Anovulatory. Other drugs may decrease the levels of Serum Iron such as: Acetylsalicylic Acid (high doses), Allopurinol, Cortisone and Corticotrophin.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 937.

REFERENCE VALUES

The reference values in µg/dL for this method were obtained through the determination of Iron in healthy populations of male and female.

Iron Serum (µg/dL)		
Newborn	100 - 250	
Infant	40 - 100	
Children	50 - 120	
Adult	Male	65 - 170
	Female	50 - 170

To convert the values of µg/dL to mmol/L, multiply by 0,179.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

RECOVERY

The recovery analysis was made with 5 sample determinations. The accurate was calculated and was found in a good parity with the reference values, obtaining a recovery between 94% and 104%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

Iron Serum kit was compared with other commercially available methods for measurement of Iron. 07 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1,011X - 1,534$, and the correlation coefficient 0,999. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration ($\mu\text{g/dL}$)	57,95	85,00	192,95
Standard Deviation ($\mu\text{g/dL}$)	0,83	0,79	0,89
Coefficient of Variation (%)	1,42	0,93	0,46

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration ($\mu\text{g/dL}$)	57,92	85,05	193,03
Standard Deviation ($\mu\text{g/dL}$)	0,25	0,05	0,08
Coefficient of Variation (%)	0,43	0,06	0,04

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Iron. The average found was 1,1506 $\mu\text{g/dL}$ and the standard deviation of 0,0005 $\mu\text{g/dL}$. The sensitivity, that indicates the methods detection limit, corresponds the average plus 3 times the standard deviation and is equal to 1,1521 $\mu\text{g/dL}$.

Linearity

The reaction is linear up to 600 $\mu\text{g/dL}$. For larger values than 600 $\mu\text{g/dL}$, dilute the serum with distilled water, repeat dosage and multiply the result by a factor of dilution.

DIAGNOSTIC SIGNIFICANCE

Iron is the most important transition metal and abundant in the human body. It is essential to the synthesis of hemoglobin and myoglobin, besides being present in several biochemical reactions of metabolism.

The serum iron is elevated in sideroblastic anemia, hemolytic anemias, hereditary anemias (Thalassemia Major), acute hepatitis, chronic liver diseases, hemochromatosis and massive blood transfusions.

Low serum iron levels are found in chronic infection diseases, neoplasms, parasites, menstrual period, large bleeding, pregnancy and dietary deficiencies.

NUMBER OF TESTS

20 Tests / 500 μL of Sample / 2mL of Reagent

40 Tests / 250 μL of Sample / 1mL of Reagent

BIBLIOGRAPHIC REFERENCES

- 1 - PETERS T., et coll.: J. Lab and Clin Med. 48:280, 1956.
- 2 - HENRY, R. J., Cannon. D. C. and Winkleman, J. W., eds.: Clinical Chemistry: Principles and Technics, Haper and Row Publishers, 1974.
- 3 - ZAK, B.: Clin Chim. Acta 3: 328, 1958.
- 4 - BRITTENHAM, G.: Clin Chim. Acta. 91: 203, 1979.
- 5 - TONKS, D. B.: Chim. Acta 2: 393, 1983.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
 CEP 31565-130 - Belo Horizonte - MG - Brasil
 Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
 email bioclin@bioclin.com.br
 CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP OBELIS S.A.

Bd. Général Wahis, 53
 1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
 Phone.: 0800 0315454
 Email: sac@bioclin.com.br

ANVISA registration for the Iron Serum kit: 10269360110

Review: April/2013

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED