CÁLCULO RENAL

REF K008

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of the usual components of Kidney Stone. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Bioclin

The physical-chemical analysis of Kidney Stone is of great clinical importance in preventive orientation from new formations. The Bioclin Kidney Stone kit allows the fast, simple and secure identification of Carbonate, Oxalate, Ammonium, Phosphate, Calcium, Magnesium, Urate and Cystine.

Carbonate: Formation of carbonic gas in acidic medium.

Oxalate: Precipitation in an acid medium (pH 3 - 4).

Phosphate: Formation of Molybdenum Blue in the presence of Ammonium Molybdate.

Calcium: Selective precipitation in alkaline medium.

Magnesium: Formation of violaceous, by complexation with specific ligand.

Ammonium: Formation yellowish with Iodine-Iodide Potassium and Mercury.

Urate: Formation of Blue Tungsten through reduction of Sodium Phosphotungstate.

Cystine: Formation of red color in the presence of Sodium Cianetonitroprussiato.

REAGENTS

All the reagents are ready to be used and are stable for 2 years if kept under room temperature. Included the set 01 bottle of Reagent N° 15 - Standard, containing all usual components from Kidney Stones (Oxalate, Calcium, Carbonate, Magnesium, Ammonium, Cystine and Uric Acid).

PRESENTATION

Kit composed with 14 Reagents and 1 Standard, enough to perform 10 complete analysis.

Reagent	Volume
N° 1 - Acid Reagent	10,0 mL
N° 2 - Oxalate	3,0 mL
N° 3 - Phosphorus I	3,0 mL
N° 4 - Phosphorus II	3,0 mL
N° 5 - Alkaline Reagent	2 x 8,0 mL
N° 6 - Calcium	5,0 mL
N° 7 - Magnesium I	7,0 mL
N° 8 - Magnesium II	10,0 mL
N° 9 - Ammonium	5,0 mL
N° 10 - Uric Acid I	6,0 mL
N° 11 - Uric Acid II	6,0 mL
N° 12 - Cystine I	2,0 mL
N° 13 - Cystine II	2,0 mL
N° 14 - Cystine III	2,0 mL
N° 15 - Standard	580 mg

EQUIPMENTS AND OPERATIONAL INPUTS

Porcelain crucible, gas beaks, centrifuge, water bath, watch or stopwatch, pipettes, test tubes. Can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperatures should be between 15 to 30°C. Protect from light and avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- If an analytical balance is not available, you can omit the stones weight, without loss of quality of information provided.

5- There is no need to weigh the amount of stones pulverized that will be used in the tests. Just a small amount in the bottom of the tube 13 x 100 mm.

6- It is convenient to proceed with the complete analysis of Standard that follows the kit, with the purpose of familiarizing yourself with the test.

7- If the stones sent are extremely low, proceed with the analysis reducing it to half the volume of all reagents utilized in step B.1.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Kidney Stones

PROCESS DESCRIPTION

TECHNIQUE

A) PHYSICAL ANALYSIS

Describe the dimension, weight (if possible), form, color, surface (smooth or irregular) of Kidney Stone. During its pulverization, observe the consistence (petrified, soft).

B) CHEMICAL ANALYSIS

1- Preparing the Sample to be analyzed

Pulverize the stones. Transfer a small amount (40 - 50 mg), of homogeneous powder, to a test tube 13 x 100 mm, add 10 drops of Reagent N° 5 and 10 drops of distilled or deionized water. Heat in water bath at 56°C, for 5 minutes, shaking the tube 2 to 3 times during this period. Centrifuge at 3000 rpm for 3 minutes. Transfer the entire supernatant to another tube 13 x 100 mm, and mark it with the letter S (for use in analysis of the supernatant) and mark the tube containing the precipitate with the letter P (for solubilized use).

2- Solubilization of the Precipitate and analysis of Carbonate

CARBONATE: To the tube marked with the letter P, add 10 drops of Reagent N° 1, observing simultaneously for outgasing. If so, the test is positive for Carbonate. Add then 10 drops of distilled water and mix. Heat to the first sign of boiling (In flame directly, hot plate or any other heating device). Leave to cool. Use this solution for analysis of Oxalate, Calcium and Magnesium.

3- Analysis of the Solution of the Precipitate

OXALATE: Take 0,1 mL from sample resulting from solubilization of the precipitate (described above) and add 3 drops of Reagent N° 2. The formation of intense turbidity or white precipitate indicates the presence of Oxalate.

CALCIUM: Take 0,1 mL sample of the resulting solubilization of the precipitate (described above) and add 5 drops of Reagent N° 6. The formation of white precipitate indicates the presence of Calcium.

MAGNESIUM: Transfer 0,05 mL of the resulting sample of solubilization of the precipitate (described above) to a Erlenmeyer and add 20 mL of distilled or deionized water. Add to this solution 1 drop of Reagent N° 5. Mix. This will be the diluted sample. Use a test tube of 12 x 75 mm and add 7 drops of Reagent N° 7 + 10 drops of Reagent N° 8. Mix. Add 0,05 mL of the diluted sample and mix. The appearance of violet color indicates the presence of Magnesium.

4- Analysis of Supernatant

URATA: Transfer 0,1 mL of the supernatant from step B.1 to a test tube 12 x 75 mm. Add 5 drops of Reagent N° 10 and 5 drops of Reagent N° 11. The appearance of intense blue color indicates the presence of Urate.

CYSTINE: Transfer 0,1 mL of supernatant from step B.1 to a test tube 12 x 75 mm. Add 1 drop of Reagent N° 12 and 1 drop of reagent N° 13. Wait 5 minutes. Add 2 drops of Reagent N° 14. The appearance of intense red color indicates the presence of Cystine.

Note: The color disappears quickly formed. Caution to manipulate Reagent N° 13. Contains Cyanide.

AMMONIUM: Transfer 0,1 mL of supernatant from step B.1 to a test tube 12 x 75 mm. Add 10 drops of distilled or deionized water. Mix. Add 5 drops of Reagent N° 9. The formation of precipitate orange - yellow indicates the presence of Ammonium.

Note: Precipitation of the Reagent N° 9 may occur, but it doesn't interfere on the test results.

The color formed can acquire a dark color if there is a high concentration of Ammonium in the sample.

PHOSPHATE: Transfer 0,1 mL of supernatant from step B.1 to a test tube 12 x 75 mm. Add 1 mL of deionized or distilled water and 01 drop of Reagent N° 1. Mix. Add 2 drops of Reagent N° 3. Mix. Add 2 drops of Reagent N° 4. Leave resting in 2 minutes and add 2 drops of Reagent N° 5. Wait for 5 minutes to read the color verification. The appearance of blue color indicates the presence of Phosphate.

RESULTS

A model result to be provided is described as follows:

Material: Kidney Stone

1. Physical Analysis

- 1.1. Weight: 50 mg
- 1.2. Dimensions: 5 x 3 mm
- 1.3. Form: Oval
- 1.4. Color: Dark gray
- 1.5. Surface: Irregular
- 1.6. Consistency: Stony

2. Chemical Analysis

- 2.1. Carbonate: Negative
- 2.2. Oxalate: Positive
- 2.3. Phosphate: Positive
- 2.4. Calcium: Positive
- 2.5. Magnesium: Negative
- 2.6. Ammonium: Negative
- 2.7. Urata: Negative
- 2.8. Cystine: Negative

3. Conclusion

The examined material contain Oxalate and Calcium Phosphate.

Note: The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

NUMBER OF TESTS

Biocal Kidney Stone kit contains reagents for 10 complete analysis.

BIBLIOGRAPHIC REFERENCES

- 1- OHWEILER, O.A.: Química Analítica, 1975.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Kidney Stone kit: 10269360091

Review: April/2015

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED