TRIGLICÉRIDES MONOREAGENTE

REF K117

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FINALIDADE

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



USAGE INSTRUCTIONS

FUNCTION

Method for determination of Triglycerides. Enzymatic colorimetric test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzymatic Colorimetric

O H₂O₂, 4 Aminoantipyrine and p-Chlorophenol in the presence of Peroxidase result in a cherry colored compound whose color intensity is proportional to the concentration of Triglycerides.

REAGENTS

Reagent N° 1 - Enzymatic Reagent - Store between 2 and 8°C. Contains: Buffer, 4-Chlorine Phenol < 5 mmol/L, Lipoprotein Lipase < 5000 U/L, Glycerol Kinase < 3000 U/L, Peroxidase < 5000 U/L, Glycerol-3-Phosphate Oxidase < 5000 U/L, 4-Amineantipyrine < 1 mmol/L, ATP < 5 mmol/L, activator, stabilizer, surfactant and preservative.

Reagent N° 2 - Standard - Store between 2 and 8°C. Contains: Triglycerides 100,0 mg/dL and diluent.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 100 mL	1 x 3 mL
2	2 x 100 mL	1 x 3 mL
3	4 x 100 mL	1 x 3 mL
4	5 x 20 mL	1 x 3 mL
5	5 x 40 mL	1 x 3 mL
6	10 x 40 mL	1 x 3 mL
7	2 x 60 mL	1 x 3 mL
8	4 x 60 mL	1 x 3 mL
9	6 x 60 mL	1 x 3 mL
10	8 x 60 mL	1 x 3 mL
11	10 x 60 mL	1 x 3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Water level in water bath should be above the level of reagents in test tubes.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis (to avoid falsely elevated results) or plasma collected with EDTA or heparin. The compound is stable for up to 3 days between 2 to 8°C, and 30 days at 10°C negative. Blood should be collected after a fasting period of 12 to 14 hours. Lipemic samples must first be diluted with Sodium Chloride 0,85%, in proportion 1:2.

PROCESS DESCRIPTION

TECHNIQUE

The Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

Mark 3 test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Sample	---	---	10 µL
Reagent N°1	1,0 mL	1,0 mL	1,0 mL
Reagent N°2	---	10 µL	---

Homogenize well and put in a water bath at 37°C for 10 minutes. Read the absorbance of Sample and Standard at 500 nm (490 - 540 nm), hitting the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

Triglycerides (mg/dL) = $\frac{\text{Sample Absorbance} \times 100}{\text{Standard Absorbance}}$

As the reaction follows the Beer-Lambert Law, the Calibration Factor can be used.

Calibration Factor = $\frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$ (100 mg/dL)

mg/dL = Sample Absorbance x Calibration Factor

Results are expressed in mg/dL.

PROCEDURE LIMITATIONS

Reducing substances such as Ascorbic Acid, even at low concentrations, and icteric samples with Bilirubin above 5 mg/dL, interfere in the methodology, leading to falsely decreased results. Some substances such as Alcohol, Oral Contraceptives and Estrogens raise Triglyceride values.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 1951.

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of Triglycerides in healthy populations of male and female.

Adults (≥ 20 years)	Desirable	< 150 mg/dL
	Limit	150 – 200 mg/dL
	Elevated	200 – 499 mg/dL
	Very Elevated	≥ 500 mg/dL

Childrens and Adolescents (2 to 19 years)	Desirable	<100 mg/dL
	Limit	100 - 129 mg/dL
	Elevated	≥ 130 mg/dL

To convert the values from mg/dL to mmol/L (SI), multiply by 0,0113.

These values should be used as a guideline only, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

Triglycerides Monoreagent kit was compared with another method for measurement of Triglycerides commercially available. 42 tests were performed and the results were evaluated. The linear equation obtained was Y = 0,969X + 3,595 and coefficient of correlation 0,997. With these results we can conclude that the kit shows good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average (mg/dL)	102,47	173,31	139,40
SD (mg/dL)	0,63	1,08	0,75
CV (%)	0,61	0,62	0,53

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average (mg/dL)	103,00	176,67	139,87
SD (mg/dL)	0,48	1,36	0,54
CV (%)	0,46	0,78	0,38

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of Triglycerides. The average was 1,34 mg/dL, with standard deviation of 0,41 mg/dL. The sensitivity that indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 2,58 mg/dL.

Linearity

The reaction is linear up to the concentration of 900 mg/dL. For samples with values above 900 mg/dL dilute the sample with Sodium Chloride 0,85%, repeat the dosage and multiply the result by dilution factor.

DIAGNOSTIC SIGNIFICANCE

Triglycerides, constituents of many lipoproteins, are found in different concentrations and are of great importance classification and phenotyping of hyperlipoproteinaemia. In several pathologies that occur in hyperlipidemia, triglycerides only are not elevated in type IIa. Values are checked increased in several diseases such as diabetes, diseases cardiovascular disease, pancreatitis, nephrotic syndrome, uremia, hypothyroidism, chronic alcoholism. In screening procedures, the determination of triglycerides, cholesterol, evidence of cooling serum and other parameters can provide consistent data to classification and phenotyping of hyperlipoproteinaemia.

BIBLIOGRAPHIC REFERENCES

- 1 - BUCOLO, G.; DAVID, H., Clin. Chem., 1973, 476.
- 2 - TONKS, D. B., Quality Control in Clinical Laboratories, 1983.
- 3 - MC GOWAN, M. W.; ARTISS, J. D.; STRANDBERGH, D. R.; ZAK, B., Clin. Chem., 1983, 29-538.
- 4 - TRINDER, P.; Ann. Clin. Biochem, 1969, 6-24.
- 5- Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
- 6 - Bioclin - Dados de arquivo.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Triglycerides Monoreagent kit:
10269360193

Review: August/2014

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED