Bioclin

PROTEÍNAS TOTAIS MONOREAGENTE

REF K031

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of total protein. Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Biuret

The peptidic bindings of the protein (-CONH-) reacts with the cupric ions, in an alkaline environment, forming a violet colored complex proportional to the protein content in the environment. The Sodium Tartrate and Potassium stabilizes the reagent and the appropriate concentration of Potassium Iodate prevents its self-reduction.

REAGENTS

Number 1 - Biuret - Store between 2 and 8°C. Contains: Sodium Hydroxide 0,2 mol/L, Sodium and Potassium Tartrate 32 mmol/L, Potassium Iodate 6 mmol/L, Cooper Sulphate 12 mmol/L.

Number 2 - Standard - Store between 2 and 8°C. Contains: Albumin 4g/dL, stabilizer and Sodium Azide 15,38 mmol/L.

PRESENTATION

| Presentation | Reagent N° 1 | Reagent N° 2 |
|--------------|--------------|--------------|
| 1 | 250 mL | 2 mL |
| 2 | 500 mL | 2 mL |
| 3 | 5 x 40 mL | 2 mL |
| 4 | 10 x 40 mL | 2 mL |

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometers or colorimeter, watch and stopwatch, pipettes and test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized in clinical analysis laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature shall be between 2 and 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not Freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Reagent N° 2 is caustic, handle carefully. Reagent N° 1 contains sodium azide, irritating to skin and mucous membrane.
- 6- Do not use plasma.
- 7- Determine the factor periodically and in each product batch.
- 8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis and liquids (pleural, synovial and ascitic). The analytical is stable for 07 days between 2 and 8°C and 60 days at -10°C.

PROCESS DESCRIPTION

REAGENT PREPARATION

The reagents are ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes: B (Blank), P (standard), A (Sample) and proceed as follows:

| | Blank | Standard | Sample |
|----------|--------|----------|--------|
| Sample | --- | --- | 50 µL |
| Standard | --- | 50 µL | --- |
| Biuret | 2,5 mL | 2,5 mL | 2,5 mL |

Homogenize well and let it rest for 10 minutes. Read the sample and standard absorbances at 545 nm (510 - 550 nm), matching the zero with the blank. The color is stable for 30 minutes.

CALCULATIONS

$$\text{Total Protein (g/dL)} = \frac{\text{Sample absorbance}}{\text{Standard Absorbance}} \times 4$$

As the reaction follows the Lambert-Beer law, the Calibration Factor can be used.

$$\text{Calibration factor} = \frac{\text{Standard Concentration (4 g/dL)}}{\text{Standard Absorbance}}$$

$$\text{g/dL} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

The results are express in g/dL

PROCEDURE LIMITATIONS

The described method is not applied to the protein dosage in urine and liquor. For this samples use the Bioclin Bioprot U/LCR kit.

INTERFERENCES

The hemolyzed sample, jaundice, lipemic or medicine such as: digitalytics, corticosteroids, anabolic steroids interfere in the dosages. The estrogen, laxatives and oral contraceptives interfere, producing reduced values of the proteins.

For samples with bilirubin above 40 mg/dL and Hemoglobin above 180 mg/dL make the Blank of the sample as follows: mixture 2,5 mL of NaCl 0,85% with 0,05 mL of sample. Determine the absorbance at 545 nm matching the zero with distilled or deionized water. Reduce the absorbance obtained from the absorbance test and calculate the concentration.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D.

REFERENCE VALUES

The reference values in g/dL, for this method were obtained through the determination of the total protein in healthy populations of male and female.

| Children and Teenagers | |
|------------------------|----------------|
| Premature | 3,6 - 6,0 g/dL |
| Newborn | 4,6 - 7,0 g/dL |
| 7 days - 1 year | 4,4 - 7,6 g/dL |
| 1 - 2 years | 5,6 - 7,5 g/dL |
| Over 3 years | 6,0 - 8,0 g/dL |
| Adults | |
| | 6,0 - 8,0 g/dL |

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

METHOD COMPARISON AND METHODOLOGICAL SPECIFICITY

The Total Protein Monoreagent kit was compared with other method for total protein dosage, commercially available. 07 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 0,970X + 0,246$ and the correlation coefficient 0,995. With these results we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

| | Sample 1 | Sample 2 | Sample 3 |
|-------------------------------------|----------|----------|----------|
| Average Concentration (g/dL) | 5,49 | 6,69 | 8,41 |
| Standard Deviation (g/dL) | 0,07 | 0,09 | 0,08 |
| Coefficient of Variation (%) | 1,36 | 1,27 | 0,90 |

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

| | Sample 1 | Sample 2 | Sample 3 |
|-------------------------------------|----------|----------|----------|
| Average Concentration (g/dL) | 5,50 | 6,71 | 8,40 |
| Standard Deviation (g/dL) | 0,01 | 0,02 | 0,01 |
| Coefficient of Variation (%) | 0,18 | 0,23 | 0,12 |

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of total protein. The average found was 0,022 g/dL with a standard deviation of 0,007 g/dL. The sensitivity, that indicates the method detection limit, corresponds to the average plus 3 times the standard deviation, and is equal to 0,043 g/dL.

Linearity

The reaction is linear up to 12 g/dL. For higher values dilute the serum with saline 0,85%, repeat the dosage and multiply the result obtained per the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The plasmatic levels of the Total Proteins suffer variation according to the changes of its various fractions. Its concentration is influenced by the nutritional state, hepatic, kidney and metabolic errors.

Increased values of the Total Protein can occur in dehydration, multiple myeloma, collagen diseases (erythematous lupus, rheumatoid arthritis), bacterial sub-acute endocarditis.

The hypoproteinemia can happen in the following clinical conditions: nephrotic syndrome, hepatic insufficiency (cirrhosis, chronic hepatitis, neoplasia), severe innutrition, severe anemia, extended fever state, fever infections, massive bleeding.

The Proteins in serum liquids (pleural and ascitic) is important to difference between transudates (Protein values below 50% of serum concentration) and exudates (Protein values superior to 50% of serum concentration).

The Protein levels in sensorial liquids are inferior to the blood varying from 1,2 a 2,5 g/dL.

BIBLIOGRAPHIC REFERENCES

- 1 - GORNALL, A. G.; BARDAWILL, C. J.; DAVID, M. M. J., Biol. Chem., 1977, 751.
- 2 - WEICHELBAUM, T. E.; AMER, J., Clin. Pathol., 1946, 16,40.
- 3 - SLATER, L.; CARTER, P. M.; HOBBS, J. R., Ann. Clin. Biochem., 1975, 12,333.
- 4 - BATSAKIS, J. G.; AROUSOHN, R. S.; WALKER, W. A.; BARNES, B.; AMER, J., Clin. Pathol., 1.976, 66,238.
- 5 - HOEL, P. G., em Estatística Elementar, Ed. Fundo de Cultura S/A, 1969.
- 6 - TONKS, D. B., Clin. Chem., 1983, 9,217.
- 7 - CARL, A. B. and EDWARD, R. A., Tietz Textbook of Clinical Chem. 2nd ed, 1994, 695-697.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUIMICA BASICA Ltda.

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Total Protein Monoreagent kit:
10269360116

Review: September/2013

UNIVERSAL SYMBOLOGY

| | | | |
|--|---------------------------------------|--|--|
| | CATALOG NUMBER | | MANUFACTURED BY |
| | BATCH CODE | | CONTROL |
| | DATE OF MANUFACTURE | | POSITIVE CONTROL |
| | USED BY (last day of month) | | NEGATIVE CONTROL |
| | TEMPERATURE LIMITATION (store at) | | BIOLOGICAL RISK |
| | CONTAINS SUFFICIENT FOR <N> TESTS | | INFLAMMABLE |
| | CONSULT INSTRUCTIONS FOR USE | | CORROSIVE |
| | IN VITRO DIAGNOSTIC DEVICE | | POISON |
| | EUROPEAN AUTHORIZED REPRESENTATIVE | | CE MARK |
| | KEEP AWAY FROM SUNLIGHT | | DO NOT USE IF PACKAGE IS DAMAGED |