## MAGNÉSIO MONOREAGENTE

REF K116

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## FINALIDADE

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for determination of Magnesium. Colorimetric test, for *in vitro* diagnostics use only.

#### PRINCIPLE OF ACTION

##### Methodology: Mann Yoe

Mann and Yoe dye in alkaline pH and in the presence of Magnesium develop a red color. The intensity of the red color of the complex is proportional to the concentration of Magnesium. The method does not require deproteinization, being more sensitive than the Titan Yellow method, allowing therefore the use of only 20 µL of the sample to be analyzed, what turns the technique excellent for pediatric use. Besides its simplicity, the method does not suffer interference from Calcium Gluconate. The presence of surfactants rules out the interference of lipemic serum.

#### Reagents

**Number 1 – Color Reagent** - Store between 2 and 8°C. Contains: Magon Sulfonate < 1 mmol/L, tampon, solubilizer, stabilizer, surfactant and preservative.

**Number 2 - Standard** - Store between 2 and 8°C. Contains: Magnesium 2 mg/dL e estabilizante.

#### PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 100 mL	1 x 3 mL
2	2 x 100 mL	1 x 3 mL
3	4 x 100 mL	1 x 3 mL
4	5 x 20 mL	1 x 3 mL
5	5 x 40 mL	1 x 3 mL
6	10 x 40 mL	1 x 3 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. Can be found in markets specialized in Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in the cleansing of the material should be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

6- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

7- Do not use the product in case of damaged packaging.

8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum obtained free of hemolysis, plasma collected with heparin and 24h urine. Magnesium is stable in sample for 24 hours if kept in temperatures between 15°C and 30°C and 15 days if kept between 2 and 8°C.

#### PROCESS DESCRIPTION

#### TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes: B (Blank), A (Sample), S (Standard) and proceed as follows:

	Blank	Standard	Sample
Reagent N° 1	1,0 mL	1,0 mL	1,0 mL
Sample	----	----	10 µL
Reagent N° 2	----	10 µL	----

Homogenize well and read the absorbance from Standard and from the Sample in 500 nm, correcting the zero with the Blank. The immediate color is stable for 30 minutes.

#### 24 hour urine – Sample preparation:

Homogenize urine and measure its volume. Add 1 drop of concentrated HCl in 20 mL of this urine. Dilute 1,0 mL of acidified urine with 4,0 mL of distilled water (1:5 dilution). Homogenize it well by agitating and proceed the determination of Magnesium the same way that was proposed to the serum. Multiply the results obtained by 5 (dilution factor).

#### CALCULATIONS

##### Serum

$$\text{Magnesium (mg/dL)} = \frac{\text{Sample Absorbance} \times 2}{\text{Standard Absorbance}}$$

Since the reaction follows the Lambert-Beer law, the Calibration Factor may be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (2 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{mg/dL} = \text{Sample absorbance} \times \text{Calibration Factor}$$

##### Urine

$$\text{Magnesium (mg/24h)} = \frac{\text{mg/dL} \times \text{Urinary Volume (mL)}}{100}$$

Results are expressed in mg/dL for serum and mg/24h for urine.

#### PROCESS LIMITATIONS

Usage of hemolyzed samples may result in falsely elevated results, since the erythrocytes contain 3 times more Magnesium than serum or plasma.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 929A.

#### REFERENCE VALUES

Reference values expressed in mg/dL, for the present method, were obtained through the determination of Magnesium in healthy populations of male and female.

Serum or plasma..... 1,6 to 2,4 mg/dL (All ages)  
Liquor..... 2,4 to 3,4 mg/dL  
Urine..... 32 to 150 mg/24h  
(Varies with according to eating habits)

In order to convert values of mg/dL into mmol/L (SI) multiply it by 0,41.

These values should be used as an orientation, being that each laboratory should create its own range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

## PRODUCT PERFORMANCE QUALITY CONTROL

### Accuracy

#### COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Magnesium Monoreagent was compared with other commercially available methods for measurement of Magnesium. 42 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 0,974X - 0,008$  and the correlation coefficient 0,995. With these results, we can conclude that the kit shows good methodological specificity.

### Precision

#### REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	1,91	3,77	2,87
<b>Standard Deviation (mg/dL)</b>	0,02	0,03	0,02
<b>Coefficient of Variation (%)</b>	1,07	0,80	0,65

#### REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	1,88	3,76	2,85
<b>Standard Deviation (mg/dL)</b>	0,03	0,02	0,03
<b>Coefficient of Variation (%)</b>	1,80	0,45	1,14

### Sensitivity

The sensitivity was calculated from 40 determinations of a sample free the Magnesium. The average found was 0,019 mg/dL and the standard deviation of 0,008 mg/dL. The sensitivity, that indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,043 mg/dL.

### Linearity

Reaction is linear up to the concentration of 5,0 mg/dL. For higher values, dilute the sample with NaCl 0,85%, repeat the determination and multiply the results by the dilution factor.

### DIAGNOSTIC SIGNIFICANCE

Magnesium is an electrolyte found primarily in intracellular liquid and bone. Participates as a cofactor in several enzymatic systems, in the metabolism of carbohydrates,

muscle contraction, blood clotting and is essential in preserving the molecular structure of DNA, RNA and ribosomes.

Magnesium levels decreased in plasma are associated with tetany, weakness, disorientation and drowsiness, which effect the ionized magnesium deficiency. Clinical manifestations seizures associated with hypocalcemia and hypomagnesemia, attributed to a selective defect in intestinal absorption of Magnesium, stressing the importance of dosage this ion in pediatrics.

The most frequent causes of low concentrations of Magnesium are: chronic diarrhea, acute pancreatitis, alcoholism, chronic hepatitis, diabetes mellitus, hypoparathyroidism, hyperthyroidism, hyperaldosteronism. Hypermagnesemia, although rare, can be registered in cases of renal failure, severe dehydration, intensive treatment with Magnesium salts.

### BIBLIOGRAPHIC REFERENCES

- 1 - MANN, C. K. & Yoe, J. H.: Ann. Chem. 28:202, 1956.
- 2 - MANN, C. K. & Yoe, J. H.: Anal. Chim. Acta. 16:155, 1957.
- 3 - RICE, E. W. & LAPARA, C. Z.: Clin. Chim. Acta. 10:260, 1964.
- 4 - WEISSMANN, N. & PILEGGI, V. J. (1974) in Clinical Chemistry Principles and Technics 2nd. Ed. Henry R., Cannon, D. C. e Winkelman, J. W. p.678. Haper and Row Publishers.
- 5 - TONKS, D. B., Clin. Chem. 9:217, 1963.
- 6 - Bioclin – Dados de arquivo

### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: + 55 (31) 3439.5454 - Fax + 55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

### EC REP OBELIS S.A.

Bd. Général Wahis, 53  
1030 Brussels, Belgium

### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for Magnesium Monoreagent kit: 10269360192

Review: January/2015

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY  
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION  
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT  
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS  
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED  
REPRESENTATIVE

CE MARK

KEEP AWAY  
FROM SUNLIGHT

DO NOT USE IF  
PACKAGE IS  
DAMAGED