**LACTATO**

REF **K084**

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**FINALIDADE**

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for determination of Lactate in plasma or human cephalorachidian liquid. UV Enzymatic test. For *in vitro* diagnostic only.

#### PRINCIPLE OF ACTION

**Methodology:** UV Enzymatic using the Dehydrogenase Lactate.

Lactic Dehydrogenase (LDH) catalyzes the oxidation of L-Lactate to Pyruvate, with consequent reduction of NAD<sup>+</sup> to NADH. The L-Lactate concentration is measured by the increase in absorbance at 340 nm, from the NADH formed.

#### REAGENTS

**Reagent N° 1 - Enzymatic Reagent** - Store between 2 and 8°C. Contains: Borate buffer (pH 10,0), LDH > 24 KU/L and Sodium Azide 15,38 mmol/L.

**Reagent N° 2 - Coenzyme** - Store between 2 and 8°C. Contains: NAD > 4 mmol/L, Citrate buffer 200 mmol/L (pH 3,0) and Sodium Azide 15,38 mmol/L.

**Reagent N° 3 - Standard** - Store between 2 and 8°C. Contains: Lactate 10 mg/dL and Sodium Azide 15,38 mmol/L.

#### PRESENTATION

Presentation	Reagent n°1	Reagent n°2	Reagent n°3
1	1 x 20 mL	1 x 5 mL	1 x 3 mL
2	2 x 20 mL	2 x 5 mL	1 x 3 mL
3	3 x 20 mL	3 x 5 mL	1 x 3 mL
4	4 x 20 mL	4 x 5 mL	1 x 3 mL
5	1 x 40 mL	1 x 10 mL	1 x 3 mL
6	2 x 40 mL	2 x 10 mL	1 x 3 mL
7	4 x 40 mL	4 x 10 mL	1 x 3 mL
8	2 x 40 mL	1 x 20 mL	1 x 3 mL
9	4 x 40 mL	2 x 20 mL	1 x 3 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. Can be found at the market specialized on Clinical Analysis Laboratories.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature is between 2 and 8°C and the transport temperature shall be between 15 and 30°C. not exceeding 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must to be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- The water level in water-bath must be superior to the reagents level in the test tubes.
- 6- Do not use jaundice or hemolyzed samples.
- 7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Customer Advisory Service) of Quibasa.
- 9- Do not use the product in case of damaged packaging.
- 10- It is essential that the instruments and equipments used ar properly calibrated and subjected to periodic maintenance.

#### SAMPLE

Plasma or cephalorachidian liquid. Do not use serum. As an anticoagulant use glycolytic inhibitors, example, fluoride/oxalate or fluoride/heparin. Do not use jaundice or hemolyzed samples. Plasma stability: 8 hours at 20 - 25°C or 14 days at 2 - 8°C.

#### PROCESS DESCRIPTION

#### REAGENT PREPARATION

Mixture 4 parts of Reagent N°1 with 1 part of Reagent N°2 (Ex.: 20 mL R1 + 5 mL R2) = Monoreagent  
Stability: 2 weeks at 2 - 8°C.

#### TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes with letters B (Blank), P (Standard), A (Sample) and proceed as follows:

	Blank	Standard	Sample
Standard	---	10 µL	---
Water	10 µL	---	---
Sample	---	---	10 µL
Monoreagent	1,0 mL	1,0 mL	1,0 mL

Smoothly homogenize, incubate for 5 minutes at 37°C and determine the Sample absorbance at 340 nm,

matching the zero with the Blank. Read the absorbance in, no longer, than 30 minutes.

#### CALCULATIONS

$$\text{Lactate (mg/dL)} = \frac{\text{Sample Abs.} \times \text{Standard Concentration}}{\text{Standard Absorbance}}$$

Following the Lambert-Beer law, we can use the Calibration Factor:

$$\text{Calibration Factor} = \frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$$

$$\text{mg/dL} = \text{Sample absorbance} \times \text{Calibration factor}$$

The values are expressed in mg/dL. To convert the values to mmol/L, multiply the obtained value by 0,1109:  
Lactate (mmol/L) = Lactate (mg/dL) x 0,1109

#### PROCEDURE LIMITATIONS

Ascorbic Acid above 30 mg/dL, direct and indirect Bilirubin above 60 mg/dL, Lipemia above 2000 mg/dL of Triglycerides, Hemoglobin above 1000 mg/dL, Dopamine above 10 mg/dL, L-Dopamine above 20 mg/L, Methyl Dopamine above 10 mg/L and Glycolic Acid above 1200 mg/L can interfere in Lactate dosage.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of Lactate in healthy populations of male and female.

PLASMA		
	mg/dL	mmol/L
Vein	4,5 - 19,8	0,5 - 2,2
Arterial	4,5 - 14,4	0,5 - 1,6
CEPHALORACHIDIAN LIQUID		
	mg/dL	mmol/L
Adults	10 - 22	1,1 - 2,4
Newborns	10 - 60	1,1 - 6,7
3 - 10 days	10 - 40	1,1 - 4,4
> 10 days	10 - 25	1,1 - 2,8

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

## PRODUCT PERFORMANCE QUALITY CONTROL

### Accuracy RECOVERY

The recovery analysis was made with 05 sample determinations. The accuracy was calculated and was found in a good concordance with the reference values, obtaining a recovery between 94 and 104%.

### COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Lactate kit was compared with other method commercially available for dosage of Lactate. 07 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 1,000X - 0,011$  and the correlation coefficient 0,999. With these results, we can conclude that the kit shows good methodological specificity.

### Precision REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	8,09	15,93	21,51
<b>Standard Deviation (mg/dL)</b>	0,17	0,19	0,08
<b>Coefficient of Variation (%)</b>	2,04	1,21	0,37

### REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	8,08	15,99	21,50
<b>Standard Deviation (mg/dL)</b>	0,02	0,05	0,02
<b>Coefficient of Variation (%)</b>	0,19	0,33	0,10

### Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Lactate. The average found was 0,9855 mg/dL with standard deviation of 0,0060 mg/dL. The sensitivity, that indicates the method detection limit,

corresponds the average plus 3 times the standard deviation and is equal to 1,0036 mg/dL.

### Linearity

The reaction is linear up to 120 mg/dL. For values above 120 mg/dL, dilute the sample in Sodium Chloride solution 9 g/L, repeat the dosage and multiply the obtained result by the dilution factor.

### DIAGNOSTIC SIGNIFICANCE

Lactate is the final product of the anaerobic glycolysis and it is used as an indicative of the oxygen state in cellular tissues.

The elevated Lactate levels in blood occur in anoxia due to shock, congestive cardiac stroke, intoxication and thiamine's deficiency. That is why the Lactate is dosed with extreme caution.

The determination of muscular Lactate is used to evaluate the metabolism variation in athletes training.

### BIBLIOGRAPHIC REFERENCES

1. Thomas L. Clinical Laboratory Diagnostics 1st ed. Frankfurt: THBooks Verlagsgesellschaft; 1998. p. 160-166.
2. David B Sacks, M. B., Ch.B. F.A.C.P. Carbohydrates In: Burtis CA, Ashwood ER, editors. Tietz Book of Clinical Chemistry. 3 ed. Philadelphia: W.B Saunders Company; 1999. p. 787-790.
3. Westgard J.O., Lahmeyer B.L., Bimbaum M.L. Use of the Du Pont "Automatic Clinical Analyser" in Direct Determination of Lactic Acid in Plasma Stabilized with sodium Fluoride. Clin chem 1972; 18:1334-8.

### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

### **OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium

### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for the Lactate kit: 10269360143

**Review:** August/2014

## UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED