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REF K132

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FINALIDADE

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USAGE INSTRUCTIONS


FUNCTION

Method for the quantitative determination of Sodium in human serum. Enzymatic test, for *In vitro* Diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic

The Sodium is enzymatic determinate through the activity of the β -Galactosidase, depending of the Sodium with the ONPG as substrate. The absorbance at 405 nm of the product o-Nitrophenyl is proportional to the Sodium concentration.

ONPG = o-Nitrophenyl - β -D-Galactopyranose

REAGENTS

Reagent N° 1: Buffer - Store between 2 and 8°C. Contains: Tris buffer, β -D-Galactosidase and preservative.

Reagent N° 2: Enzymatic Reagent - Store between 2 and 8°C. Contains: Tris buffer, o-Nitrophenyl- β -D-Galactopyranose and preservative.

Reagent N° 3: Standard Level 1 - Store between 2 and 8°C. Contains: Sodium Chloride (concentration varies to each lot).

Reagent N° 4: Standard Level 2 - Store between 2 and 8°C. Contains: Sodium Chloride (concentration varies to each lot).

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2	Reagent N° 3	Reagent N° 4
1	40 mL	20 mL	3 mL	3 mL
2	2 x 40 mL	2 x 20 mL	3 mL	3 mL
3	3 x 40 mL	3 x 20 mL	3 mL	3 mL
4	4 x 40 mL	4 x 20 mL	3 mL	3 mL
5	30 mL	10 mL	3 mL	3 mL
6	2 x 30 mL	2 x 10 mL	3 mL	3 mL
7	3 x 30 mL	3 x 10 mL	3 mL	3 mL
8	4 x 30 mL	4 x 10 mL	3 mL	3 mL
9	60 mL	20 mL	3 mL	3 mL
10	2 x 60 mL	2 x 20 mL	3 mL	3 mL
11	3 x 60 mL	3 x 20 mL	3 mL	3 mL
12	4 x 60 mL	4 x 20 mL	3 mL	3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Although the substrate contains preservative, must avoid bacterial contamination.

5- Do not ingest. Avoid contact with eyes and skin.

6- The samples must be treated as potentially infectious.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum. The analyte is stable for 4 days between 2 and 8°C and 14 days at -20°C.

PROCESS DESCRIPTION

REAGENT PREPARATION

The reagents are ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

The kit is indicated for biochemistry automated analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

RESULTS

The results are expressed in mmol/L.

PROCEDURE LIMITATIONS

Test periodically and calibrate to detect possible changes in the response of the equipment.

INTERFERENCE

The test does not suffer interference from the substances with the top concentrations of:

NH₄Cl 1,5 mmol/L

KP_i 2,0 mmol/L

CaCl₂ 7,5 mmol/L

KCl 10 mmol/L

CuCl₂ 0,5 mmol/L

ZnCl₂ 0,5 mmol/L

FeCl₃ 0,5 mmol/L

Glucose 5 mmol/L

Ascorbate 10 mmol/L

Bilirubin 40mg/dL

Bilirubin conjugated 40 mg/dL

Hemoglobin 500 mg/dL

Triglycerides 1000 mg/dL

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standards are traceable to the reference material NIST (National Institute of Standards and Technology) SRM 919B.

REFERENCE VALUES

The reference values in mmol/L, for this method, were obtained through the determination of sodium in healthy populations of male and female.

Serum	136 - 146 mmol/L
	313 - 336 mg/dL

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE**QUALITY CONTROL****Accuracy****COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY**

The Sodium Ions Automated was compared with another method commercially available to measure sodium. 07 analysis were conducted and the results were evaluated. The linear equation obtained was $Y = 1,051X - 7,383$, with a correlation coefficient 0,998. With these results we can conclude that the kit has good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mmol/L)	122,15	141,75	150,95
Standard Deviation (mmol/L)	0,67	0,72	0,76
Coefficient of Variation (%)	0,55	0,51	0,50

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mmol/L)	121,92	141,82	151,08
Standard Deviation (mmol/L)	0,21	0,16	0,13
Coefficient of Variation (%)	0,17	0,11	0,08

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Sodium. The average found was 78,50 mmol/L with standard deviation of 0,513 mmol/L. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 80,04 mmol/L.

Linearity

The Sodium Ions Automated kit is linear between 80 and 180 mmol/L.

DIAGNOSTIC SIGNIFICANCE

In a healthy subject, the level of the Sodium extra cellular is regulated between 136-146 mmol/L (313-336 mg/dL). Small deviations of the regular level may cause great consequences in health. Sodium dosages are used in the diagnostic and treatment of aldosteronism, diabetes insipidus, adrenal hypertension, Addison disease, dehydration, inadequate secretion of the antidiuretic hormone, or other disease related to the electrolytic disequilibrium, metabolic and cardiovascular disorders. Therefore, the monitoring of the Sodium concentration in serum is important in both routine and emergency checking.

BIBLIOGRAPHIC REFERENCES

- Berry, M. N. et al., (1988) Clin. Chem. 34,2295.
- Tietz, N. W. (1983) Clinical guide to Laboratory Tests, p. 384 W.B. Saunders Co., Philadelphia.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

**OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Sodium Ions Automated kit:
10269360211

Review: September/2015

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED