DESIDROGENASE LÁTICA LDH UV

REF K014

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Lactate Dehydrogenase (LDH). Kinetic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Kinetic

The Lactate Dehydrogenase (LDH) catalyzes the reduction of Pyruvate with NADH, and NAD⁺ to Lactate. The catalytic concentration is determined by the speed decomposition of NADH, measured by the drop in absorptivity at 340 nm.

REAGENTS

Number 1 – Buffered Substrate - Store between 2 and 8°C. Contains: Tris Buffer < 200 mmol/L, Pyruvate < 6 mmol/L, stabilizer and preservative.

Number 2 - Coenzyme - Store between 2 and 8°C. Contains: Buffer < 20 mmol/L, NADH < 5,0 mmol/L and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	54 mL	6 mL
2	1 x 40 mL	1 x 10 mL
3	2 x 40 mL	2 x 10 mL
4	4 x 40 mL	4 x 10 mL
5	2 x 40 mL	1 x 20 mL
6	4 x 40 mL	2 x 20 mL
7	4 x 27 mL	1 x 12 mL
8	6 x 36 mL	1 x 24 mL
9	6 x 36 mL	2 x 12 mL
10	3 x 36 mL	1 x 12 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatch, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- It is important, for the good development of the test, a rigorous control of time and temperature.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis or plasma collected with EDTA or heparin. The LDH in serum or plasma for 2 days if kept in temperature between 2 and 8°C.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENT

Mix nine (9) parts of Reagent N° 1 with one (1) part of Reagent N° 2. Working Reagent is stable during 14 days, if kept in temperatures between 2 and 8°C.

REACTION CONDITIONS

It is indispensable to use thermostated cuvette at 37°C, and 1 cm optical path reading at 340 nm.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Add 20 μ L of Sample to 1,0 mL of Working Reagent, mixing and transferring to a thermostated cuvette at 37°C and wait 1 minute. Make the initial reading simultaneously starting the stopwatch. Repeat readings after 1, 2 and 3 minutes. Calculate the mean differences in absorbance per minute ($\Delta A/\text{min.}$) and use it to calculate the result.

CALCULATIONS

LDH (U/L) = $\Delta A/\text{min.} \times 8016$
Results are expressed as U/L.

PROCEDURE LIMITATIONS

The specifications below refers to semi-automated equipments:

The kinetic method is based on the absorptivity and by this reason, the readings must be conducted in a spectrophotometer that satisfies the following conditions: Wavelength 340 nm.

Semi trajectory of the pass band 10 nm

Stray light less than 0,5%

1cm thermostated cuvette

INTERFERENCES

Hemolysis or late separation of serum could cause elevated results due to a high concentration of LDH in erythrocytes. Lipemia (Tryglicerides > 1000 mg/dL) and Bilirubin (> 20 mg/dL) could result in falsely elevated results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The calibration kit can be made using the theoretical calibration factor based on the molar absorptivity of NADH, or by calibrator BIOCAL. The Bioclin recommends using BIOCAL calibrator which is traceable to the reference material ERM-AD453/IFCC.

REFERENCE VALUES

The reference values in U/L, for this method were obtained through the determination of LDH in healthy populations of male and female.

Children and Adolescents	
1 - 3 years	490 - 730 U/L
4 - 9 years	320 - 520 U/L
10 - 13 years	250 - 500 U/L
Adults	
200 - 480 U/L	

These values should be use as reference, each laboratory must create their own range of values for reference, according to the served population.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE**QUALITY CONTROL****Accuracy****COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY.**

The kit was compared with other method commercially available for dosage of LDH. 42 analysis were performed and the results were evaluated. The linear equation obtained was $Y = 0,939X + 5,054$ and the coefficient of correlation 0,990. With these results we can conclude that the kit shows good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	283,13	537,10	418,18
Standard Deviation (U/L)	1,42	3,68	2,00
Coefficient of Variation (%)	0,50	0,68	0,48

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	289,30	538,41	416,19
Standard Deviation (U/L)	5,79	8,07	6,41
Coefficient of Variation (%)	2,00	1,50	1,54

Sensitivity

The sensitivity was calculated from 40 determinations from a sample free of LDH. The average found was 6,00 U/L with a standard deviation of 0,51 U/L. The sensitivity, that indicates the method detection limit corresponds the average plus 3 times the standard deviation and is equal to 7,52 U/L.

Linearity

The reaction is linear to the concentration of 2000 U/L. For an average change in absorbance at 340 nm greater than 0,12, repeat the determination, diluting the sample with 0,85% NaCl. Multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Lactate Dehydrogenase (LDH) is present in all body cells, and the highest concentrations are in liver, heart, kidney, muscle skeletal and erythrocytes. The concentration of LDH in serum or plasma is increased in patients with liver diseases, renal disorders, myocardial infarction, many malignant diseases, progressive muscular dystrophy and in any case of hemolysis. The diagnosis should not be done taking into account only the result of a single assay of LDH, but it should be integrate clinical data with the laboratory.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until piration date stated on the presentation packaing, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium


CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Lactate Dehydrogenase LDH UV kit: 10269360074

Review: March/2015

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED