COLESTEROL MONOREAGENTE

REF K083

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Cholesterol. Enzymatic colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic Colorimetric - COD - PAP

The enzymatic determination of Cholesterol is made accordingly with the following reactions:

The intensity of formed cherry color is directly proportional to the concentration of Cholesterol in the sample.

REAGENTS

Reagent N° 1 - Enzymatic Reagent - Store between 2 and 8°C. Buffer (pH 7,0) 75 mmol/L, Phenol 4,5 mmol/L, 4-Aminoantipyrine 0,3 mmol/L, Cholesterol Oxidase > 200 U/L, Lipoprotein Lipase > 700 U/L, Peroxidase > 300 U/L, Sodium Azide 14,6 mmol/L, stabilizers and surfactants.

Reagent N° 2 - Standard - Store between 2 and 8°C. Contains: Cholesterol 200,0 mg/dL, stabilizers and solubilizers.

PRESENTATION

PRESENTATION	REAGENT N° 1	REAGENT N° 2
1	100 mL	3 mL
2	2 x 100 mL	3 mL
3	4 x 100 mL	3 mL
4	5 x 20 mL	3 mL
5	5 x 40 mL	3 mL
6	10 x 40 mL	3 mL
7	2 x 60 mL	3 mL
8	4 x 60 mL	3 mL
9	6 x 60 mL	3 mL
10	8 x 60 mL	3 mL
11	10 x 60 mL	3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water bath (37°C), watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Water level in water bath should be above the level of reagents in test tubes.

6- Plasma collected with anticoagulants such as oxalate, EDTA or citrate, generate wrongly lowered results.

7- Handle reagents with care. Reagent N° 1 contains Sodium Azide, irritant to skin and mucous.

8- Slightly pink color development in the Reagent N° 1 does not damage their quality, and stability as long as the corresponding Blank and Standard dosages are used.

9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

10- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

11- Do not use the product in case of damaged packaging.

12- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolyses or plasma collected with heparin. Cholesterol is stable in plasma or serum for up to 7 days between 2 and 8°C and up to 05 days between 15 and 30°C.

PROCESS DESCRIPTION

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark three test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Sample	--	--	10 µL
Reagent N°2	--	10 µL	--
Reagent N° 1	1,0 mL	1,0 mL	1,0 mL

Homogenize well and put it in water bath 37°C for 10 minutes. Read the absorbance of Sample and Standard at 500 nm (490 - 540 nm), matching the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

$$\text{Cholesterol (mg/dL)} = \frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times 200$$

Since reaction follows the Lambert-Beer Law, the Calibration Factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (200 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{mg/dL} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

The results shall be express as mg/dL.

PROCESS LIMITATIONS

The equilibrium of the assay is affected by some interference substances such as the Ascorbic Acid (even in small concentrations), Hemoglobin higher than 150 mg/dL and Bilirubin higher than 20 mg/dL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 911C.

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of Cholesterol in populations of male and female.

Adults	Desirable	< 200 mg/dL
	Limit	200 - 239 mg/dL
	Elevated	≥ 240 mg/dL

Childrens and adolescents (2 to 19 years)	Desirable	< 150 mg/dL
	Limit	150 - 169 mg/dL
	Elevated	≥ 170 mg/dL

To convert the values from mg/dL into mmol/L (SI) multiply by 0,026.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

RECOVERY

The recovery analysis was mad through 5 determinations of samples. Accuracies were calculated and were found in good concordance to the reference values, obtaining a recuperation between 95% and 105%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Cholesterol Monoreagent kit was compared with another method commercially available to measure the dosage of

Cholesterol. 07 analyzes were performed and results were evaluated. The linear equation obtained was $Y = 0,992X + 2,032$ with correlation coefficient 0,999. With these results we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	101,75	204,75	246,35
Standard Deviation (mg/dL)	0,72	0,79	0,67
Correlation Coefficient (%)	0,70	0,38	0,27

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	101,70	204,60	246,40
Standard Deviation (mg/dL)	0,05	0,13	0,05
Correlation Coefficient (%)	0,05	0,06	0,02

Sensitivity

Sensitivity was calculated through 20 determinations of a sample free of cholesterol. The average found was 1,461 mg/dL with a standard deviation of 0,004 mg/dL, the sensitivity, which indicates the limit of detection from the method, corresponds the average plus 3 times the standard deviation, and is equal to 1,472 mg/dL.

Linearity

Reaction is linear up to 500 mg/dL. For samples with higher values than 500 mg/dL, or optical density higher than 0,8, dilute the sample with Sodium Chloride 0,85%, repeat the dosage and multiply the results by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The concentration of plasmatic Cholesterol is influenced by hereditary characters, endocrine function, nutrition and integrity of vital organs such as liver and kidneys. Numerous investigations confirm the relationship between Total Cholesterol and the development of atherosclerotic coronary disease.

The Cholesterol is increased in diabetes, nephrotic syndrome, cirrhosis bile in hypothyroidism and in hyperlipoproteinemia type IIa, IIb and III. Research shows that high levels of LDL Cholesterol (Cholesterol Bound to Lipoproteins (LP) Low-Density) are closely related with Ischemic Coronary Heart Disease (IHD). Rather, the elevation HDL Cholesterol (Cholesterol Bound to High Density Lipoprotein) represents a protective factor against DCI. Hypertension, smoking, obesity are other causes responsible for atherosclerosis and DCI.

Lower values may occur in the presence of diseases that affect the liver parenchyma (viral hepatitis, toxic hepatitis), occasionally acute infections (pneumonia, typhoid fever), hyperthyroidism, anemia, malnutrition.

BIBLIOGRAPHIC REFERENCES

- 1 - ALLAIN, C.C., and al. Clin. Chem., 1974. 20-470.
- 2 - TONKS, D.B., Quality Control in Clinical laboratories, 1983.
- 3 - TRINDER, P., Ann. Clin. bioche.,, 1969, 6-24.
- 4 - HENRY, J.B., Diagnosis and Management by laboratory Methods, 19th ed., 1996.
- 5 - CARL, A.B. and EDWARDR .A., Tiedz Textbook of Clinical Chem. 2nd ed., 1994,1002-1081.
- 6- Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium


CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Cholesterol Monoreagent kit:
10269360141

Review: August/2014

UNIVERSAL SYMBOLOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED