Bioclin

CK NAC UV

REF K010

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FINALIDADE

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of the creatine kinase (CK). Kinetic test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Kinetic UV

Determination of creatine kinase according to the reactions:

CK - Creatine Kinase

G-6-PDH – Glucose-6-Phosphate-Dehydrogenase

HK - Hexokinase

The reduction speed of NADP⁺ to NADPH is proportional to the activity of the CK in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Imidazole Acetate (pH 6,7) 100 mmol/L, Glucose 20 mmol/L, EDTA 2 mmol/L, NADP 2 mmol/L, Hexokinase 3500 U/L, Magnesium Acetate 10 mmol/L and N-Acetylcysteine 20 mmol/L.

Number 2 - Enzyme - Substrate - Store between 2 and 8°C. Contains: Glucose-6-Phosphate-Dehydrogenase 2000 U/L, Creatine Phosphate 30 mmol/L, ADP 2 mmol/L, AMP 5 mmol/L and Diadenosine Pentaphosphate 10 mmol/L.

PRESENTATION

Reagent	Presentation				
	1	2	3	4	5
N° 1	1 x 40 mL	2 x 40 mL	4 x 40 mL	2 x 40 mL	4 x 40 mL
N° 2	1 x 10 mL	2 x 10 mL	4 x 10 mL	1 x 20 mL	2 x 20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, test tubes, watch or stopwatch, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperature between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- It is important for the proper performance of the test, a strict control of time, temperature and pH.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma collected with heparin or EDTA. The sample is stable for 24 hours between 15 and 30°C and a week if kept between 2 and 8°C.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENT

Mix 4 parts of reagent N° 1 with 1 part of reagent N° 2. Stable for 2 days at temperatures between 15 to 30°C and 15 days if kept between 2 to 8°C. Do not use hemolyzed sample.

REACTION CONDITIONS

It is indispensable the use of a thermostated cuvette at 37°C, light path of 1cm and reading at 340 nm.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Working Reagent	1,0 mL
Sample	20 µL

Mix and transfer to a thermostated cuvette at 37°C and wait 2 minutes. Make the initial reading, triggering the timer. Repeat the readings at 1, 2 and 3 minutes.

CALCULATIONS

Calculate the mean difference in absorbance per minute ($\Delta A/\text{min.}$) and use this value to calculate the result.

CK (U/L) 340 nm = $\Delta A/\text{min} \times 8095$

PROCEDURE LIMITATIONS

The readings should be taken in a spectrophotometer that satisfies the following conditions:

Wavelength 340 nm

Semi trajectory of the pass band 10 nm

Spurious light lower than 0,5%

10 mm Thermostated cuvette.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the NADP, or through the BIOCAL calibrator which is traceable to the reference material ERM-AD455 / IFCC and the IFCC reference method.

REFERENCE VALUES

The reference values for this method were obtained by determining CK in healthy populations of male and female.

Adults (37°C)	
Male	24 - 195 U/L
Female	24 - 170 U/L

Children and Teenagers (37°C)		
Age	Male	Female
1 - 3 month	29 - 303 U/L	40 - 474 U/L
4 - 12 month	25 - 172 U/L	27 - 242 U/L
13 - 24 month	28 - 162 U/L	25 - 177 U/L
2 - 10 years	31 - 152 U/L	25 - 177 U/L
11 - 14 years	31 - 152 U/L	31 - 172 U/L
15 - 18 years	34 - 147 U/L	28 - 142 U/L

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The CK NAC UV kit was compared to other method for CK NAC dosage available in the market. 07 analysis were performed and the results were evaluated. The linear equation obtained was $Y = 0,994X + 1,043$, with linear correlation coefficient equal to 0,999. With these results we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	27,65	101,90	191,95
Standard Deviation (U/L)	0,75	0,79	0,89
Variation Coefficient (%)	2,69	0,77	0,46

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	27,75	101,78	191,98
Standard Deviation (U/L)	0,13	0,13	0,03
Variation Coefficient (%)	0,48	0,12	0,02

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of the analyte (CK NAC). The mean concentration of CK NAC UV was 5,45 U/L with standard deviation of 1,32 U/L. The sensitivity, which indicates the limit of detection, which corresponds to the average plus three times the standard deviation is equal to 9,40 U/L.

Linearity

The method is linear to 2000 U/L. If the variation of molar extinction per minute is larger than 0,250, repeat the test using a dilution of 1:10 with physiologic serum. The result obtained is multiplied by 10.

DIAGNOSTIC SIGNIFICANCE

The enzyme Creatine Phosphokinase (CPK) and Creatine Kinase (CK) is found in high concentrations in cardiac and skeletal muscles and brain.

The values of CK are elevated in the presence of muscle damage as in myocardial infarction. Following the attack, the serum concentration of this enzyme increases rapidly (3 to 6 hours), reaching a peak after 12 to 24 hours and remains elevated for a short period of 2 to 3 days. Other potential sources of elevated total CK are: hypothyroidism, muscle disease (myopathy, polymyositis), stroke, seizures, cardiac catheterization, multiple trauma, strenuous and prolonged immobilization.

BIBLIOGRAPHIC REFERENCES

- 1 -SCANDINAVIAN Society For Clinical Chemistry Scand. J. Clin. Lab. INVEST., 1974, 33, 291.
- 2 - Ann. Biol. Clin., 1982, 40,99.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for CK NAC UV kit: 10269360092

Review: February/2013

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED