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REF K107

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Direct Bilirubin in samples serum or plasma. Colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Colorimetric test

Direct Bilirubin, through the coupling reaction with 2,4 Dichloroaniline diazotised, forms an azocompound, a red colored complex with maximum absorption at 546 nm. The intensity of color formed is directly proportional to the concentration of Direct Bilirubin in the sample.

REAGENTS

Reagent Nº 1 Buffer - Store between 2 and 8°C. Contains: Sulfamic Acid 70 mmol/L, stabilizer and preservative.

Reagent Nº 2 Color Reagent - Store between 2 and 8°C. Contains: 2,4 Dichlorophenyl-Diazonium Salt 0,09 mmol/L, HCl 130 mmol/L, stabilizer.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2
1	1 x 40 mL	1 x 10 mL
2	2 X 40 mL	2 x 10 mL
3	3 X 40 mL	3 x 10 mL
4	4 X 40 mL	4 x 10 mL
5	5 X 40 mL	5 x 10 mL
6	2 x 40 mL	1 x 20 mL
7	4 x 40 mL	2 x 20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the biochemistry automated equipment, Biocal kits and Biocontrol Bioclin.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- The reagent Nº 2 must be kept protected from light.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Calibration must be repeated periodically to verify any alterations in colorimetric responses or spectrophotometer.

6- Hemolysis, even if discreet, interfere in dosage.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma free of hemolysis.

Sample must be protected from light.

In samples stored at temperatures between 2 to 8°C and protected from light, Bilirubin is stable for three (03) days. Sample for therapeutic control should always be collected at the same hour.

PROCESS DESCRIPTION

PREPARATION OF REAGENTS

Reagents are ready for use.

TECHNIQUE

For the use of the kit, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol Bioclin kit.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

INTERFERENCES

Hemolysis, even if slight, interfere with measurement of Bilirubin. No interference was observed for Ascorbic Acid up to 30 mg/dL and lipemia up to 1000 mg/L from Triglycerides. According to bibliographic data, some drugs, as Anabolic Steroids, Ascorbic Acid, Salicylates and Vitamin A, can raise bilirubin values. Caffeine and Penicillins may lead to falsely decreased.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention

that all measurement systems present an analytical variety, and it must be monitored by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 916A.

REFERENCE VALUES

Reference values were obtained through determinations of Direct Bilirubin in healthy populations of male and female.

For all ages:

Direct Bilirubin: up to 0,4 mg/dL

To convert the values to mg/dL to mmol/L (SI), multiply the results obtained by 0,0171.

These values should be used only as a guide. Each laboratory should establish its own reference range, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS

Bilirubin Direct Automated kit was compared with other methods for determination of Direct Bilirubin dosage commercially available measurement. 42 analyzes were performed and the results were evaluated. The equation linear obtained was $Y = 0,999X - 0,001$ and the correlation coefficient 0,99996. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0,49	1,45	0,86
Standard Deviation (mg/dL)	0,01	0,01	0,01
Coefficient of Variation (%)	1,63	0,48	0,72

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0,50	1,46	0,86
Standard Deviation (mg/dL)	0,004	0,01	0,01
Coefficient of Variation (%)	0,78	0,52	0,58

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of Direct Bilirubin. The average found was 0,016 mg/dL, with deviation standard of 0,003 mg/dL. The sensitivity, which indicates the method limit detection, corresponds the average plus 3 times the standard deviation, and is equal to 0,025 mg/dL.

Linearity

The reaction is linear to the concentration of 12 mg/dL. For higher values dilute the sample with NaCl 0,85% and repeat the dose. Multiply the resulting number by the dilution factor used.

DIAGNOSTIC SIGNIFICANCE

Bilirubin is a breakdown product of hemoglobin in reticuloendothelial system. It is conjugated in the liver to the then be excreted in bile.

The level of Bilirubin is useful for differential diagnosis of hepatobiliary disease and other causes of jaundice which becomes clinically evident when the concentration of Total Bilirubin greater than 2,5 mg/dL.

The Direct Bilirubin rises in plasma in the presence of hereditary liver diseases such as disease Dubin-Johnson and Rotor, hepatocyte injury (viral, toxic or alcoholic), biliary obstruction (stones, neoplasm), acute or chronic hepatitis and toxic reactions to various drugs (such as Chlorpromazine, and organic arsenic compounds, Methyltestosterone, among others). Levels of Direct Bilirubin over 50% of the Total Bilirubin concentration indicate the existence of post-liver.

The use of drugs that activate the hepatic microsomal system can reduce the concentration of Total and Direct Bilirubin in plasma.

BIBLIOGRAPHIC REFERENCES

1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:192-202.
2. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors, Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
3. Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962; 6:570-8.
4. Bioclin – Dados de arquivo

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454.
E-mail: sac@bioclin.com.br

ANVISA registration for Bilirubin Direct Automated kit:
10269360176

Review: August/2014

UNIVERSAL SYMBOLOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED