## DENGUE BIO

REF K105

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#### FINALIDADE

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- Halstead SB. Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB. Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.]
- Bioclin - Dados de Arquivo.

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- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.]
- Bioclin - Datos de Archivo.

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### USAGE INSTRUCTIONS

#### FUNCTION

Quick Immunochromatographic. Test for determination of antibodies IgM and IgG for the Dengue virus in whole blood, serum or plasma samples.

#### PRINCIPLE OF ACTION

**Methodology:** immunochromatography

The Dengue Bio kit is an immunochromatographic test who detects separately classes antibodies IgM and IgG for the Dengue for virus. The kit contains antibodies of goat anti-human IgG in the test line for IgG, and goat anti-human IgM. In test line for IgM. After the sample adding, antigens of the Dengue virus conjugated to gold colloidal particles, present in the kit, bind to the sample antibodies and migrate by capillarity through the membrane. In case of positive samples for IgG, it will form a red line in the IgM test area. Positive samples for IgM will form a red line in the IgM test area. The mixture of the antigen-sample continues to run through the membrane until the control line, where occurs a formation of a red line, for the positive and negative samples. The blue line formation demonstrates that the kit is adequate for use.

#### REAGENTS

**Reagent N° 1: Cassette:** Store between 15 and 30°C. Do not freeze. Contains:

- C control region: rabbit IgG anti goat antibody.
- G test region (T1): goat anti-human IgG antibody.
- M test region (T2): goat antihuman IgM antibody.
- Conjugate: Dengue virus antigen (I, II, III and IV) conjugated with the gold colloidal particles.

**Reagent N° 2: Diluent** - Store between 15 and 30°C. Do not freeze. Contains: Tris buffer 100 mmol/L pH 7,4, preservative.

#### PRESENTATION

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	1	10	1,00 mL
	2	25	2,50 mL
	3	50	2 x 2,50 mL
	4	100	4 x 2,50 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Watches and stopwatches. They can be found at markets specialized in laboratories of clinical analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 15 and 30°C. Protect from light and avoid moisture.

#### SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- The cassettes do not present contaminant substances, but once the samples are biological, the safety routine must be followed strictly and the use of disposable gloves is indispensable.
- 4- Use clean and dry recipients to settle the samples.
- 5- The reagent N° 2 contain sodium azide, irritating to skin and mucous membranes. Handle the reagent with care.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLE

Whole Blood: shall be stored between 2 and 8°C and tested In 24 hours after collection.

Plasma or serum: shall be separate from blood interface immediately after centrifugation, and must occur instantly after the collection. The plasma or serum can be stored for 3 days between 2 and 8°C and 6 months at -20°C.

#### PROCESS DESCRIPTION

- 1- The sample must be in a temperature between 15 and 30°C before begin the test;
- 2- Take off the reactive strip of the protection case and identify it in an adequate form;
- 3- Transfer 10 µL of the total blood, plasma or serum inside the sample window;
- 4- Hold the flask vertically and apply 2 drops (100 µL) of diluent in the sample pool;
- 5- Wait the formation of the lines after 15-20 minutes. Do not interpret it after 20 minutes.

#### RESULTS INTERPRETATION



**IgG positive test:** a red line forms in the test area G (T1) and another line in the control area (C) in the first 15 to 20 minutes.

**IgM positive test:** a red line forms in the test area M (T2) and another line in the control area (C) in the first 15 to 20 minutes.

**IgM and IgG positive test:** two red line forms in the test area (G = T1 and M = T2) and a red line in the control area (C) in the first 15 to 20 minutes.

**Negative test:** a red line forms in the control area (C) and complete absence of red lines in the test areas (G = T1 and M = T2).

**Invalid Test:** complete absence of red lines in the control area (C) and test regions (G = T1 and M = T2).

#### NOTES:

- 1- The negative results must only be read in 20 minutes, and not before or after this time.
- 2- The line in the control area must show before 10 minutes of incubation test, and that does not mean the negative results can be interpreted before the time.
- 3- The results must not be interpreted after 20 minutes in order to avoid incorrect results.
- 4- The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

#### INTERFERENCES

Cross reactivity with other virus of the flavivirus can be common, according to bibliographic references, however there was no evidence in the product developing tests.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### LIMITATIONS

- 1- The Bioclin's Bio Dengue kit can only be used for the *in vitro* diagnosis.
- 2- The Bioclin's Bio Dengue Kit determines only the the presence of antibodies against the dengue virus (I, II, III and IV) and shall not be used isolated as a diagnosis criteria. Positive results should be evaluated with other clinical information available. The isolated results of the test must not be used for a definitive diagnosis.
- 3- A negative result does not extinct the possibility of infection by the Dengue virus (I, II, III and IV). Scientific studies reveal that infected people can present an immunological window of 5 to 10 days for IgM and IgG detection, respectively. The patient shall be evaluated again during this period.
- 4- Immunosuppressed patients' results should be carried out with caution.
- 5- In every diagnosis test, the results must be interpreted within other patients clinical data.
- 6- If the negative result continues and the clinical symptoms persists, an additional test can be carried out by another method. A negative result does not extinct, in any way, the possibility of infection by the dengue virus.

## PRODUCT PERFORMANCE

### QUALITY CONTROL

#### Accuracy

#### METHODS COMPARISON, METHODOLOGICAL SPECIFICITY AND SENSITIVITY

To establish the sensitivity and specificity of the Dengue Bio kit in relation to other Dengue test, 450 known clinical samples were studied. Other qualitative kit available in the market was used to compare with the kit and for determination of relative sensitivity and specificity. The found clinical sensitivity was 96% and the specificity found was of 95,8%.

#### Precision

#### REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 5 samples with different concentrations, obtaining the following results:

Samples	Repetition Nº	Expected Result	Found result
Blood 01	20	100% Positive IgG	100% Positive IgG
Blood 02	20	100% Negative	100% Negative
Serum 01	20	100% Negative	100% Negative
Serum 02	20	100% Positive IgM	100% Positive IgM
Plasma 01	20	100% Positive IgG/IgM	100% Positive IgG/IgM

#### REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Samples	Nº of Days	Repetition Nº	Expected Result	Found Result
Blood 01	3	20	100% Negative	100% Negative
Blood 02	3	20	100% Positive IgM	100% Positive IgM
Serum 01	3	20	100% Positive IgG/IgM	100% Positive IgG/IgM
Serum 02	3	20	100% Negative	100% Negative
Plasma 01	3	20	100% Positive IgG/IgM	100% Positive IgG/IgM

#### DIAGNOSTIC SIGNIFICANCE

The Dengue virus is a flavivirus transmitted by the Aedes aegypti mosquito. It is distributed in the tropical and subtropical areas worldwide, causing over 100 million infections yearly. The classic Dengue is characterized by high fever (normally between 38 and 40°C of suddenly beginning, low appetite, minor illnesses, caphaleia, muscle and eyes sore. Hemorrhagic Dengue can cause gingivorrhagy and epistaxis, internal bleeding and disseminated intravascular coagulation, with many vital organs damage, occasioning in death. The primary infection of dengue causes a increase in the IgM antibodies after 3 to 5 days in the beginning of the fever. IgM antibodies normally remain in the circulation for 30 to 90 days. Patients of endemic regions can present the secondary infections, that result in high levels of IgG antibody, isolated or simultaneously as a response of IgM. Therefore, an specific detection of antibodies igM or IgG can help to difference the primary or secondary.

## BIBLIOGRAPHIC REFERENCES

- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.
- Bioclin - Dados de Arquivo.

#### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

#### OBELIS S.A.

Bd. Général Wahis, 53  
1030 Brussels, Belgium

#### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for Dengue Bio kit: 10269360174

**Review:** September/2013

## UNIVERSAL SYMBOLOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY  
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION  
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT  
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS  
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED  
REPRESENTATIVE

CE MARK

KEEP AWAY  
FROM SUNLIGHT

DO NOT USE IF  
PACKAGE IS  
DAMAGED