## FERRITINA

REF K081

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#### FINALIDADE

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## REFERÊNCIAS

1. Kaltwasser JP, Werner E. Diagnosis and clinical evaluation of iron overload. Baillieres Clin Haematol 1989; 2; 363-89.
2. Wick M, Pingerra W, Lehmann P. Iron metabolism: diagnosis and therapy of anemias. 3<sup>rd</sup> ed. Vienna, New York: Springer Verlag, 1996.
3. Lee MH, Means RT Jr. Extremely elevated serum ferritin levels in a university hospital: associated diseases and clinical significance. Am J Med 1996; 98:566-71.
4. Worwood M. The laboratory assessment of iron status – an update. Clin Chim Acta 1997; 259:3-23.
5. Baynes RD, Cook JD. Current issues in iron deficiency. Curr Opin Hematol 1996; 3:145-9.

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for quantitative determination of Ferritin in serum, in photometric systems. Immunoturbidimetric test, for *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Immunoturbidimetry

Determination of ferritin concentration's ending point through photometric measuring of antigen-antibody reaction, between latex particles marked with antibodies anti-ferritin and ferritin present in sample.

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: Glycine 170 mmol/L (pH 8,3), NaCl 100 mmol/L and Albumin 3 g/L.

**Number 2 - Ferritin Latex** - Store between 2 and 8°C. Contains: Latex particles coated with bunny antibodies anti-ferritin 0,23 g/L, Glycine 170 mmol/L and NaCl 100 mmol/L.

**Number 3 - Calibrator** - Store between 2 and 8°C. Contains: Ferritin and Sodium Azide 15,38 mmol/L. **Potentially infectious.**

**Number 4 - Control** - Store between 2 and 8°C. Contains: Ferritin and Sodium Azide 15,38 mmol/L. **Potentially infectious.**

**Warning:** Concentration of Ferritin, control and calibrator may vary according to the lot - See bottle label.

#### PRESENTATION

Reagent	Presentation				
	1	2	3	4	5
Nº1	18 mL	45 mL	2 x 45 mL	40 mL	2 x 45 mL
Nº2	2 mL	5 mL	2 x 5 mL	10 mL	2 x 10 mL
Nº3	0,5 mL	1,0 mL	2 x 1,0 mL	1,0 mL	2 x 1,0 mL
Nº4	0,5 mL	1,0 mL	2 x 1,0 mL	1,0 mL	2 x 1,0 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

This kit must be used in automated equipments only. They can be found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For professional *in vitro* diagnostic use only.**

**2- Strictly follow the methodology proposed to obtain exact results.**

**3- Water used in material cleaning must be recent and free of contaminants.**

**4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.**

**5- Water level in water bath must be superior to the levels of reagents in test tubes.**

**6- Reagents Nº 3 and 4 should be handled cautiously as they are subject to biological contamination.**

**7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.**

**8-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Advisory Service Customer) of Quibasa.

**9-** Do not use the product in case of damaged packaging.

**10-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum. Ferritin is stable for 07 days if kept between 2 and 8°C or up to 3 weeks if kept at - 20°C. Sample must be frozen only once. Do not use highly lipemic or hemolysed samples.

#### TECHNIQUE

For reaction control and calibration, use only calibrators and controls from this kit.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC 0800 031 5454.**

#### PROCEDURE LIMITATIONS

No interference was observed by Ascorbic Acid up to 30 mg/dL, Bilirubin up to 20 mg/dL and Hemoglobin up to 1000 mg/dL. The photometric scale of some instruments can be exceeded by highly lipemic samples, containing levels elevated ferritin.

The performance of the product depends directly on the performance of the automated equipment used. Therefore, to ensure best results of the test is very important to check periodically the optical and pipetting system of the equipment.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

**Do not compare the obtained results by different methodologies.**

#### TRACEABILITY

The kit's calibrator is traceable to the reference material NIBSC 94/572 (Ferritin Human Recombinant - WHO International Standard).

#### REFERENCE VALUES

The reference values in  $\mu\text{g/L}$  or  $\text{ng/mL}$ , for this method were obtained through the determination of Ferritin in healthy populations of male and female.

Men: 30 - 200  $\mu\text{g/L}$  or  $\text{ng/mL}$

Women: 20 - 110  $\mu\text{g/L}$  or  $\text{ng/mL}$

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

#### PRODUCT PERFORMANCE

##### QUALITY CONTROL

##### Accuracy

##### RECOVERY

The recovery analysis was performed with samples of 05 determinations. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 94% and 107%.

##### COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Ferritin kit was compared with other commercially available methods for measurement of ferritin. 07 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 0,995X + 0,676$ , and correlation coefficient 0,999. With these results we can conclude the kit shows good methodological specificity.

##### Precision

##### REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration ( $\mu\text{g/L}$ )	60,65	110,85	220,00
Standard Deviation ( $\mu\text{g/L}$ )	0,75	0,88	0,86
Coefficient of Variation (%)	1,23	0,79	0,39

##### REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration ( $\mu\text{g/L}$ )	60,77	110,90	220,02
Standard Deviation ( $\mu\text{g/L}$ )	0,16	0,05	0,03
Coefficient of Variation (%)	0,26	0,05	0,01

##### Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Ferritin. The average found was of 5,090  $\mu\text{g/dL}$  with a standard deviation of 0,0013  $\mu\text{g/dL}$ . The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 5,094  $\mu\text{g/L}$ .

##### Linearity

This kit is linear between 15  $\mu\text{g/L}$  and the highest point in the calibration curve. For higher values, dilute the sample with a solution of NaCl 9 g/L. Multiply the results by the dilution factor.

##### SIGNIFICANCE DIAGNOSTIC

Ferritin is the most important iron storage protein and is found in all cells, especially those involved in the synthesis of ferric compounds, metabolism and it's reserve.

The level of ferritin is the most accurate indicator of the amount of iron stored in the body.

His clinic is very useful in differential diagnosis of microcytic and hypochromic anemia iron deficiency, anemia of other etiologies. In these cases, the ferritin decreases before changes of serum iron levels and red blood cell morphological changes.

However, to be part of the group of acute phase proteins, ferritin is increased in response to infection, trauma and acute inflammation. The elevation occurs in 24 to 48 hours, with a peak in the third day, and it stays for a few weeks, which may hinder its interpretation. Its levels can rise in excess of iron in transfused patients and cancer, especially in leukemias and lymphomas and breast cancer, liver, lung, colon and prostate.

Elevates also in megaloblastic and hemolytic anemias and in liver injury, especially alcoholic beverages. About 25% of patients with chronic hepatitis have increased ferritin.

#### BIBLIOGRAPHIC REFERENCES

1. Kaltwasser JP, Werner E. Diagnosis and clinical evaluation of iron overload. *Baillieres Clin Haematol* 1989; 2; 363-89.
2. Wick M, Pingerra W, Lehmann P. Iron metabolism: diagnosis and therapy of anemias. 3rd ed. Vienna, New York: Springer Verlag, 1996.
3. Lee MH, Means RT Jr. Extremely elevated serum ferritin levels in a university hospital: associated diseases and clinical significance. *Am J Med* 1996; 98:566-71.
4. Worwood M. The laboratory assessment of iron status – an update. *Clin Chim Acta* 1997; 259:3-23.
5. Baynes RD, Cook JD. Current issues in iron deficiency. *Curr Opin Hematol* 1996; 3:145-9.

#### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: [bioclin@bioclin.com.br](mailto:bioclin@bioclin.com.br)  
CNPJ: 19.400.787/0001-07 - Made in Brazil

#### **OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium



#### CUSTOMER SERVICE

Customer Advisory Service  
Phone: 0800 0315454  
E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

ANVISA registration for Ferritin kit: 10269360135

**Review:** August/2014

## UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED