ALFA 1 GLICOPROTEÍNA ÁCIDA

REF: K060

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for the determination of Alpha 1 Acid Glycoprotein. Immunoturbidimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry

The reagent allows the quantitative determination of Alpha 1 Acid Glycoprotein in human serum by immunoprecipitation, in the presence of an activator polymer that increases the sensitivity and immunoturbidimetric test speed. The Alpha1 Glycoprotein Acid to forms along with the specific antiserum an insoluble complex, producing turbidity, whose absorbance is proportional to the concentration of Alpha 1 Acid Glycoprotein in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Sodium Chloride 0,15 mol/L, Tris 50 mmol/L, PEG 6.000 50 g/L, Sodium Azide 15,38 mmol/L and surfactant.

Number 2 - Antiserum - Store between 2 and 8°C. Contains: Antiserum Anti Alpha 1 Acid Glycoprotein and Sodium Azide 15,38 mmol/L.

PRESENTATION

Reagent	Presentation			
	1	2	3	4
Nº 1	25 mL	2 x 25 mL	3 x 25 mL	4 x 25 mL
Nº 2	2 mL	2 x 2 mL	3 x 2 mL	4 x 2 mL

Calibrator not included in kit.

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the biochemistry automated equipment, Multical kit and Multicontrol Bioclin.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- For professional *in vitro* diagnostic use only.
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Water level in water bath should be above the level of reagents in test tubes.

6- Reagent 1 and 2 should be handled cautiously because they are subject to biological contamination.

7- Handle with care as the reagent contains Sodium Azide, which is irritant to skin and mucous.

8- The biological originated material has been tested for HIV and HBs Antigen using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. Potentially Infectious.

9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

10- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

11- Do not use the product in case of damaged packaging.

12- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis and intense lipemia. Samples obtained are stable in temperatures between 2 and 8°C for up to 7 days. Freezing the samples can result in deterioration and should be avoided.

PROCESS DESCRIPTION

WORKING REAGENT PREPARATION

Dilute a portion of the Antiserum (Reagent Nº 2) with 12,5 portions of Buffer (Reagent Nº1).

Example: 1 mL of Antiserum + 12,5 mL of Buffer.

The diluted Antiserum, stored in an amber flask, it is stable for 07 days in temperatures between 2 and 8°C and 01 day in temperature between 15 and 30°C.

TECHNIQUE

For calibration and control of the reaction, utilize the Multical kits and Multicontrol Bioclin.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

PROCEDURE LIMITATIONS

Hemolysis, icterus and lipemia (low or moderate) do not interfere in performance of the tests, using the recommended dilution for sample. Gross lipemia and turbid samples should be processed with a higher dilution.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

Serum Adults 30 - 120 mg / dL

These values should be used as a guideline, and each laboratory should establish its range of reference values, according with the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy RECOVERY

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 98 and 106%.

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

Alpha 1 Acid Glycoprotein kit was compared with other method for determination of Alpha 1 Acid Glycoprotein dosage commercially available. 07 tests were performed and the results were evaluated. The linear equation obtained was $Y = 1,0498 + 0,9995X$ and correlation coefficient 0,99946. With these results we can conclude that the kit has good methodological specificity

Precision REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	153,4	75,25	116,20
Standard Deviation (mg/dL)	3,050	1,713	1,322
Coefficient of Variation (%)	1,989	2,276	1,138

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	154,37	75,38	116,0
Standard Deviation (mg/dL)	1,19	0,37	0,35
Coefficient of Variation (%)	0,77	0,49	0,30

Sensitivity

Sensitivity was calculated from 20 determinations of a sample free of the Alpha 1 Acid Glycoprotein. The average found was 1,55 mg/dL, with standard deviation of 0,60 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 3,35 mg/dL.

Linearity

The reaction is linear up to the level of the highest calibrator concentration.

DIAGNOSTIC SIGNIFICANCE

Main component of the Winzler Mucoproteins, rises considerably in acute inflammatory processes. It is an important index of rheumatic activity, because it remains high when other tests normalize. The increase in values also happens in erythematosus lupus, neoplasm metastasis and administration of estrogen. Decrease of values occurs in liver failure, adrenal insufficiency, malnutrition and pituitary insufficiency .

NUMBER OF TESTS**Presentation Nº1**

25 Tests / 50 µL of Sample / 1000 µL of Reagent
50 Tests / 25 µL of Sample / 500 µL of Reagent

Presentation Nº2

50 Tests / 50 µL of Sample / 1000 µL of Reagent
100 Tests / 25 µL of Sample / 500 µL of Reagent

Presentation Nº3

75 Tests / 50 µL of Sample / 1000 µL of Reagent
150 Tests / 25 µL of Sample / 500 µL of Reagent

Presentation Nº4

100 Tests / 50 µL of Sample / 1000 µL of Reagent
200 Tests / 25 µL of Sample / 500 µL of Reagent

BIBLIOGRAPHIC REFERENCES

1 - KILLINGSWORTH, L. M.; SAVORY, J.: Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem., 1971; 17, 936.
2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry, C. V. Mosby Company, 1987.

3 - HELLSING, K., Profides in the Biological Fluids, 1973; 23, 579.

4 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22, 657.

5 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980; 26, 1459.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium


CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
Email: sac@bioclin.com.br

ANVISA registration for Alpha 1 Acid Glycoprotein kit:
10269360061

Review: April/2014

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED