HBsAg

REF K145

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FINALIDADE

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1. Blumberg, B.S. The Discovery of Australian Antigen and its Relation to Viral Hepatitis. Vitro. 1971;7:223.
2. Krugman, S. Glies J.P. Viral Hepatitis, Type B (MS-2-Strain). Further Observations on Natural History and Prevention. New England Journal of Medicine. 288, 755.
3. Krugman, S. Overby L.R, et al. Viral Hepatitis Type B Studies On Natural History and Prevention Reexamined. New England Journal of Medicine. 300, 101.
4. Bioclin – Dados de arquivos

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1. Blumberg, B.S. The Discovery of Australian Antigen and its Relation to Viral Hepatitis. *Vitro*. 1971;7:223.
2. Krugman, S. Glikes J.P. Viral Hepatitis, Type B (MS-2-Strain). Further Observations on Natural History and Prevention. *New England Journal of Medicine*. 288, 755.
3. Krugman, S. Overby L.R, et al. Viral Hepatitis Type B Studies On Natural History and Prevention Reexamined. *New England Journal of Medicine*. 300, 101.
4. Bioclin – Datos de archivos

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USAGE INSTRUCTIONS

FUNCTION

Immunochromatographic test for rapid determination of qualitative Hepatitis B surface antigen (HBsAg) in serum samples.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The HBsAg kit is an immunochromatographic test that contains IgG conjugate (anti HBs) - colloidal gold conjugate that binds to the HBsAg surface antigen presented in the sample and forms a complex antigen-antibody. The complex will migrate through the membrane by capillarity action in the direction of the anti-HBsAg antibody in the (T) Test region, when it is linked occurs the appearance of the strip that determinates the positive reaction. In the absence of the HBsAg surface antigen, the strip will not appear. The complex will continue to migrate until it reaches the (C) Control region. The conjugate not linked to the antigen binds to the anti-IgG antibody of the C region forming a colored strip indicating that the test is working as planned and that the result is valid.

REAGENTS:

Reagent Nº 1: Cassette: Store between 15 and 30°C. Do not freeze. Contains:

- Control Region (C): Antibody anti-IgG
- Region Test (T): Anti-HBsAg goat antibody
- Conjugate: Conjugate IgG (Anti-HBs) - colloidal gold

PRESENTATION

Reagent	Presentation	Number of Cassettes for Package
Cassette	1	1
	2	5
	3	10
	4	15
	5	20
	6	25
	7	30
	8	40
	9	50
	10	100

EQUIPMENTS AND OPERATIONAL INPUTS

Watch or stopwatch. Can be found in specialized markets for articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The temperature of storage and transportation must be 15 to 30°C. Avoid moisture and protect from light.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The cassettes do not show contaminants, but how the samples are biological safety routines should be followed accurately and the use of disposable gloves is essential.
- 4- Use clean and dry containers for storage of samples.

5- Treat all samples as potentially infecting.

6- The kit should not be frozen or exposed to high temperatures, as this will cause deterioration.

7- The reagents and the sample must be equilibrated environment before starting the test.

8- Open the envelope of aluminum, to remove the cassette, only in time of use.

9- It is interesting to conduct the test immediately after sample collection.

10- The kit must not be used after expiration date.

11- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

12- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

13- Do not use the product in case of damaged packaging.

14- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum: The samples can be stored under refrigeration, between 2 and 8°C for a maximum period of 3 days. If the samples can not be analyzed within those three days, they can be stored at -20°C for up to 30 days.

DESCRIPTION OF PROCESS

TECHNIQUE

1- The sample should be at room temperature between 15 and 30°C before starting the test.


2- Remove the cassette from its protective packaging, place it on a clean and leveled surface and appropriately identify it so.

3- Transfer 100 µL serum in the sample well.


4- Wait for the formation of lines after resting from 10 to 20 minutes. Do not take a reading 30 minutes past the reaction.

INTERPRETATION OF RESULTS


Testing Positive: Formation of two lines in the (C) Control region and (T) Test region in the first 10 to 20 minutes. Do not take a reading 30 minutes past the reaction

Test Negative: formation of one line in the (C) Control region and complete absence of a line in the (T) Test region.

Inadequate Testing: Complete absence of line on the (C) Control region with or without a line in the (T) Test region. Repeat the sample test.

Notes:

- 1- Negative results can only be read in 20 minutes, nor before or later than that.
- 2- The line in the Control region will appear before the 20 minutes of the test incubation, and that does not mean that the negative results can be interpreted before the time.
- 3- The positive result will not change after 30 minutes, but the result must not be interpreted after this time to avoid incorrect results.
- 4- The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

INTERFERENCES

Hemolyses with hemoglobin values up to 10 g/L do not interfere in the HBsAg dosage. No significant interference were found in up to 30 g/L of lipids, 0,2 g/L of bilirubin and 60 g/L of proteins.

INTERNAL QUALITY CONTROL

The Clinical laboratory must have an internal quality control program where procedures, regulations, limits, and tolerance to variations are clearly established. It is important that every measurement system present a analytical variability, that must be monitored by the laboratories. Therefore, it is recommended the use of controls that allow the dosages accuracy and precision evaluation.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS, METHODOLOGY AND SPECIFICITY SENSITIVITY

The HBsAg kit was compared to other commercially available kits for the HBsAg testing. For this test two hundred and twenty eight (228) clinical samples were used. The following results were found:

HBsAg Bioclin	Kit Results Reference		Subtotal
	+	-	
Tested Kit Results	+	129	130
	-	0	98
Subtotal	129	99	228

The found results suggest the following values:

Clinical Sensitivity: 100% (129/129)

Specificity: 98,99% (98/99)

Analytical Sensitivity: 10 UI/mL

Accuracy**REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 5 different samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Results	Found Results
01	20	100% Negative	100% Negative
02	20	100% Positive	100% Positive
03	20	100% Positive	100% Positive
04	20	100% Negative	100% Negative
05	20	100% Positive	100% Positive

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

Samples	Nº of Days	Nº of Repetitions	Expected Results	Found Results
01	03	20	100% Negative	100% Negative
02	03	20	100% Positive	100% Positive
03	03	20	100% Positive	100% Positive
04	03	20	100% Negative	100% Negative
05	03	20	100% Positive	100% Positive

DIAGNOSTIC SIGNIFICANCE

HBsAg is one of the first markers to appear on blood after the infection by the Hepatitis B virus (HBV). This liver infection is transmitted by sexual contact, blood exposure, from mother to fetus during the birth or by sharing cutting/needled objects. The main subtypes of HBsAg include ad and Ay, sharing the determinant 'a'. The infection by HBV causes a great variety of liver damage, such as acute self-limited, fulminating hepatitis, chronic hepatitis, and asymptomatic chronic carrier state. The HBV virus in infected people, persist for the rest of their lives and can be transmitted to other people. Therefore, Hepatitis B had become a public health matter. The HBV infection results in a number of serum markers and one of the first markers is the Surface Antigen of Hepatitis B (HBsAg). The HBsAg appears in 1 - 10 weeks after the exposure and before the biochemistry evidences of the hepatic disease or icterus. Three weeks after the beginning of the acute hepatitis, almost half of the patients will be positive for HBsAg. In the chronic carrier state, the HBsAg persist for 6 - 12 months, in serum, without conversion-serum for the for the correspondent antibody.

For that reason, the classification for HBsAg is highly recommended for all donors, pregnant women and people in the risk groups.

NUMBER OF TESTS

Reagents	Presentation	Number of Tests
Cassettes	1	1
	2	5
	3	10
	4	15
	5	20
	6	25
	7	30
	8	40
	9	50
	10	100

BIBLIOGRAPHIC REFERENCES

1. Blumberg, B.S. The Discovery of Australian Antigen and its Relation to Viral Hepatitis. *Vitro*. 1971;7:223.
2. Krugman, S. Glies J.P. Viral Hepatitis, Type B (MS-2-Strain). Further Observations on Natural History and Prevention. *New England Journal of Medicine*. 288, 755.
3. Krugman, S. Overby L.R, et al. Viral Hepatitis Type B Studies On Natural History and Prevention Reexamined. *New England Journal of Medicine*. 300, 101.
4. Bioclin – Dados de arquivos

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for HBsAg kit: 10269360234

Review: March/2015

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED